

Anesthesia Labels & Tape

PDC Healthcare's Anesthesia Labels & Tape help address errors caused by the use of unlabeled medications and other procedural breakdowns that can put patients at risk. They adhere to specific guidelines set forth by The Joint Commission, United States Pharmacopoeia (USP), and the American Society for Testing and Materials (ASTM), including:

1. Specific color and format standards for each drug class
2. Appropriate information on each tape/label (drug name, strength, amount, date, time, initials)
3. TALLman lettering to clearly differentiate common look-alike, sound-alike drugs
4. Tape and labels with standard pre-printed concentrations

Our stock, pre-printed Anesthesia Tape and Labels help increase patient safety by delivering higher accuracy in medication administration. If you do not see the anesthesia products needed for your facility, PDC Healthcare can customize any tape or label with the desired information.



Table of Contents

Anesthesia Labels & Tapes	80
Anesthesia Labels & Tapes Custom.....	86
Dispensers.....	86

Anesthesia Tape & Labels help healthcare organizations meet applicable patient safety regulatory guidelines, including:



- **The Joint Commission National Patient Safety Goals:** Identify Patients Correctly, Improve Staff Communication, Use Medicines Safely, Prevent Infection, Identify Patient Safety Risks, and Prevent Mistakes in Surgery.
- **Five Rights of Patient Safety:** to ensure right patient, medication, dose, time, and method.

ANESTHESIA

ANESTHESIA LABELS & TAPE ALPHABETICALLY LISTED

The following stock drug tapes and labels are available for accurate and reliable labeling of all anesthesia drug syringes. Each tape/label conforms to the nationally recognized ASTM Standard D4774-11(2017) Specification for User Applied Drug Labels in Anesthesiology. This specification covers the size, color, pattern, and type used on labels applied to unlabeled syringes, filled by the users or their agents, to identify the drug content. For custom imprinting options, please see page 86. Copyright ASTM International. All rights reserved.

In addition, tapes and labels meet the published 2013 Joint Commission National Patient Safety Goal 3 to improve the safety of using medications. Common pre-printed medication strengths are available. Tall man lettering helps reduce error due to look-alike/sound-alike drugs.

<u>Tape Products</u>	<u>Label Products</u>
Product Code: AN-___ or ANS-___ (Use "AN" or "ANS" when ordering tape products)	Product Code: LAN-___ or LANS-___ (Use "LAN" or "LANS" when ordering label products)
Materials: Coated paper tape with removable adhesive	Materials: Litho Permanent
Quantity: 1/2" x 500" roll - 333 imprints per roll with perforations	Quantity: 600 labels per roll
Label Size: 1-1/2" x 1/2"	Label Size: 1-1/2" x 1/2"
Core size: 1"	Core size: 1"
	

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description	Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description																																																																																																																
2-Chloroprocaine (alk) ___ % Date ___ Time ___ Int. ___	AN-214	Tape	2-Chloroprocaine (alk) Local Anesthetics (Color: Gray)	Benadryl ___ mg/ml Date ___ Time ___ Int. ___	AN-109	Tape	Benadryl Miscellaneous (Color: White)																																																																																																																
	LAN-214	Label			LAN-109	Label		Adrenalin ___ mg/ml Date ___ Time ___ Int. ___	AN-31	Tape	Adrenalin Vasopressors (Color: Violet/Black)	___ mg/ml Date ___ Time ___ Int. ___	AN-84	Tape	Blank (Write drug name) Miscellaneous (Color: White)	LAN-31	Label	LAN-84	Label	Alfentanil ___ mcg/ml Date ___ Time ___ Int. ___	AN-106	Tape	Alfentanil Narcotics (Color: Lt. Blue)	Brevital 10 mg/ml Date ___ Time ___ Int. ___	AN-41	Tape	Brevital 10 mg/ml Induction Agents (Color: Yellow)	LAN-106	Label	LAN-41	Label	Alfentanil 500 mcg/ml Date ___ Time ___ Int. ___	AN-106D500	Tape	Alfentanil 500 mcg/ml Narcotics (Color: Lt. Blue)	Bupivacaine ___ mg/ml Date ___ Time ___ Int. ___	AN-42	Tape	Bupivacaine Local Anesthetics (Color: Gray)	LAN-106D500	Label	LAN-42	Label	Anectine ___ mg/ml Date ___ Time ___ Int. ___	AN-33	Tape	Anectine Muscle Relaxants (Color: Fl. Red/Black)	Bupivacaine ___ % Date ___ Time ___ Int. ___	AN-42PC	Tape	Bupivacaine ___ % Local Anesthetics (Color: Gray)	LAN-33	Label	LAN-42PC	Label	Antibiotic ___ mg/ml Date ___ Time ___ Int. ___	AN-35	Tape	Antibiotic Miscellaneous (Color: White)	Bupivacaine 0.5% Date ___ Time ___ Int. ___	AN-42D05	Tape	Bupivacaine 0.5% Local Anesthetics (Color: Gray)	LAN-35	Label	LAN-42D05	Label	ATRAcurium ___ mg/ml Date ___ Time ___ Int. ___	AN-39	Tape	ATRAcurium Muscle Relaxants (Color: Fl. Red)	Bupivacaine 0.25% 2.5 mg/ml Date ___ Time ___ Int. ___	AN-42D25	Tape	Bupivacaine 0.25% 2.5 mg/ml Local Anesthetics (Color: Gray)	LAN-39	Label	LAN-42D25	Label	ATRAcurium 10 mg/ml Date ___ Time ___ Int. ___	AN-39D10	Tape	ATRAcurium 10 mg/ml Muscle Relaxants (Color: Fl. Red)	Calcium Chloride ___ mg/ml Date ___ Time ___ Int. ___	AN-110	Tape	Calcium Chloride Miscellaneous (Color: White)	LAN-39D10	Label	LAN-110	Label	Atropine ___ mg/ml Date ___ Time ___ Int. ___	AN-1	Tape	Atropine Anticholinergic Agents (Color: Green)	Cisatracurium ___ mg/ml Date ___ Time ___ Int. ___	AN-162	Tape	Cisatracurium Muscle Relaxants (Color: Fl. Red)	LAN-1	Label	LAN-162	Label	Atropine 0.4 mg/ml Date ___ Time ___ Int. ___	AN-1D04	Tape	Atropine 0.4 mg/ml Anticholinergic Agents (Color: Green)	Cisatracurium 2 mg/ml Date ___ Time ___ Int. ___	AN-162D2	Tape	Cisatracurium 2 mg/ml Muscle Relaxants (Color: Fl. Red)	LAN-1D04	Label	LAN-162D2	Label	Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)
Adrenalin ___ mg/ml Date ___ Time ___ Int. ___	AN-31	Tape	Adrenalin Vasopressors (Color: Violet/Black)	___ mg/ml Date ___ Time ___ Int. ___	AN-84	Tape	Blank (Write drug name) Miscellaneous (Color: White)																																																																																																																
	LAN-31	Label			LAN-84	Label		Alfentanil ___ mcg/ml Date ___ Time ___ Int. ___	AN-106	Tape	Alfentanil Narcotics (Color: Lt. Blue)	Brevital 10 mg/ml Date ___ Time ___ Int. ___	AN-41	Tape	Brevital 10 mg/ml Induction Agents (Color: Yellow)	LAN-106	Label	LAN-41	Label	Alfentanil 500 mcg/ml Date ___ Time ___ Int. ___	AN-106D500	Tape	Alfentanil 500 mcg/ml Narcotics (Color: Lt. Blue)	Bupivacaine ___ mg/ml Date ___ Time ___ Int. ___	AN-42	Tape	Bupivacaine Local Anesthetics (Color: Gray)	LAN-106D500	Label	LAN-42	Label	Anectine ___ mg/ml Date ___ Time ___ Int. ___	AN-33	Tape	Anectine Muscle Relaxants (Color: Fl. Red/Black)	Bupivacaine ___ % Date ___ Time ___ Int. ___	AN-42PC	Tape	Bupivacaine ___ % Local Anesthetics (Color: Gray)	LAN-33	Label	LAN-42PC	Label	Antibiotic ___ mg/ml Date ___ Time ___ Int. ___	AN-35	Tape	Antibiotic Miscellaneous (Color: White)	Bupivacaine 0.5% Date ___ Time ___ Int. ___	AN-42D05	Tape	Bupivacaine 0.5% Local Anesthetics (Color: Gray)	LAN-35	Label	LAN-42D05	Label	ATRAcurium ___ mg/ml Date ___ Time ___ Int. ___	AN-39	Tape	ATRAcurium Muscle Relaxants (Color: Fl. Red)	Bupivacaine 0.25% 2.5 mg/ml Date ___ Time ___ Int. ___	AN-42D25	Tape	Bupivacaine 0.25% 2.5 mg/ml Local Anesthetics (Color: Gray)	LAN-39	Label	LAN-42D25	Label	ATRAcurium 10 mg/ml Date ___ Time ___ Int. ___	AN-39D10	Tape	ATRAcurium 10 mg/ml Muscle Relaxants (Color: Fl. Red)	Calcium Chloride ___ mg/ml Date ___ Time ___ Int. ___	AN-110	Tape	Calcium Chloride Miscellaneous (Color: White)	LAN-39D10	Label	LAN-110	Label	Atropine ___ mg/ml Date ___ Time ___ Int. ___	AN-1	Tape	Atropine Anticholinergic Agents (Color: Green)	Cisatracurium ___ mg/ml Date ___ Time ___ Int. ___	AN-162	Tape	Cisatracurium Muscle Relaxants (Color: Fl. Red)	LAN-1	Label	LAN-162	Label	Atropine 0.4 mg/ml Date ___ Time ___ Int. ___	AN-1D04	Tape	Atropine 0.4 mg/ml Anticholinergic Agents (Color: Green)	Cisatracurium 2 mg/ml Date ___ Time ___ Int. ___	AN-162D2	Tape	Cisatracurium 2 mg/ml Muscle Relaxants (Color: Fl. Red)	LAN-1D04	Label	LAN-162D2	Label	Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)	Curare ___ mg/ml Date ___ Time ___ Int. ___	AN-2	Tape	Curare Muscle Relaxants (Color: Fl. Red)								
Alfentanil ___ mcg/ml Date ___ Time ___ Int. ___	AN-106	Tape	Alfentanil Narcotics (Color: Lt. Blue)	Brevital 10 mg/ml Date ___ Time ___ Int. ___	AN-41	Tape	Brevital 10 mg/ml Induction Agents (Color: Yellow)																																																																																																																
	LAN-106	Label			LAN-41	Label		Alfentanil 500 mcg/ml Date ___ Time ___ Int. ___	AN-106D500	Tape	Alfentanil 500 mcg/ml Narcotics (Color: Lt. Blue)	Bupivacaine ___ mg/ml Date ___ Time ___ Int. ___	AN-42	Tape	Bupivacaine Local Anesthetics (Color: Gray)	LAN-106D500	Label	LAN-42	Label	Anectine ___ mg/ml Date ___ Time ___ Int. ___	AN-33	Tape	Anectine Muscle Relaxants (Color: Fl. Red/Black)	Bupivacaine ___ % Date ___ Time ___ Int. ___	AN-42PC	Tape	Bupivacaine ___ % Local Anesthetics (Color: Gray)	LAN-33	Label	LAN-42PC	Label	Antibiotic ___ mg/ml Date ___ Time ___ Int. ___	AN-35	Tape	Antibiotic Miscellaneous (Color: White)	Bupivacaine 0.5% Date ___ Time ___ Int. ___	AN-42D05	Tape	Bupivacaine 0.5% Local Anesthetics (Color: Gray)	LAN-35	Label	LAN-42D05	Label	ATRAcurium ___ mg/ml Date ___ Time ___ Int. ___	AN-39	Tape	ATRAcurium Muscle Relaxants (Color: Fl. Red)	Bupivacaine 0.25% 2.5 mg/ml Date ___ Time ___ Int. ___	AN-42D25	Tape	Bupivacaine 0.25% 2.5 mg/ml Local Anesthetics (Color: Gray)	LAN-39	Label	LAN-42D25	Label	ATRAcurium 10 mg/ml Date ___ Time ___ Int. ___	AN-39D10	Tape	ATRAcurium 10 mg/ml Muscle Relaxants (Color: Fl. Red)	Calcium Chloride ___ mg/ml Date ___ Time ___ Int. ___	AN-110	Tape	Calcium Chloride Miscellaneous (Color: White)	LAN-39D10	Label	LAN-110	Label	Atropine ___ mg/ml Date ___ Time ___ Int. ___	AN-1	Tape	Atropine Anticholinergic Agents (Color: Green)	Cisatracurium ___ mg/ml Date ___ Time ___ Int. ___	AN-162	Tape	Cisatracurium Muscle Relaxants (Color: Fl. Red)	LAN-1	Label	LAN-162	Label	Atropine 0.4 mg/ml Date ___ Time ___ Int. ___	AN-1D04	Tape	Atropine 0.4 mg/ml Anticholinergic Agents (Color: Green)	Cisatracurium 2 mg/ml Date ___ Time ___ Int. ___	AN-162D2	Tape	Cisatracurium 2 mg/ml Muscle Relaxants (Color: Fl. Red)	LAN-1D04	Label	LAN-162D2	Label	Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)	Curare ___ mg/ml Date ___ Time ___ Int. ___	AN-2	Tape	Curare Muscle Relaxants (Color: Fl. Red)																				
Alfentanil 500 mcg/ml Date ___ Time ___ Int. ___	AN-106D500	Tape	Alfentanil 500 mcg/ml Narcotics (Color: Lt. Blue)	Bupivacaine ___ mg/ml Date ___ Time ___ Int. ___	AN-42	Tape	Bupivacaine Local Anesthetics (Color: Gray)																																																																																																																
	LAN-106D500	Label			LAN-42	Label		Anectine ___ mg/ml Date ___ Time ___ Int. ___	AN-33	Tape	Anectine Muscle Relaxants (Color: Fl. Red/Black)	Bupivacaine ___ % Date ___ Time ___ Int. ___	AN-42PC	Tape	Bupivacaine ___ % Local Anesthetics (Color: Gray)	LAN-33	Label	LAN-42PC	Label	Antibiotic ___ mg/ml Date ___ Time ___ Int. ___	AN-35	Tape	Antibiotic Miscellaneous (Color: White)	Bupivacaine 0.5% Date ___ Time ___ Int. ___	AN-42D05	Tape	Bupivacaine 0.5% Local Anesthetics (Color: Gray)	LAN-35	Label	LAN-42D05	Label	ATRAcurium ___ mg/ml Date ___ Time ___ Int. ___	AN-39	Tape	ATRAcurium Muscle Relaxants (Color: Fl. Red)	Bupivacaine 0.25% 2.5 mg/ml Date ___ Time ___ Int. ___	AN-42D25	Tape	Bupivacaine 0.25% 2.5 mg/ml Local Anesthetics (Color: Gray)	LAN-39	Label	LAN-42D25	Label	ATRAcurium 10 mg/ml Date ___ Time ___ Int. ___	AN-39D10	Tape	ATRAcurium 10 mg/ml Muscle Relaxants (Color: Fl. Red)	Calcium Chloride ___ mg/ml Date ___ Time ___ Int. ___	AN-110	Tape	Calcium Chloride Miscellaneous (Color: White)	LAN-39D10	Label	LAN-110	Label	Atropine ___ mg/ml Date ___ Time ___ Int. ___	AN-1	Tape	Atropine Anticholinergic Agents (Color: Green)	Cisatracurium ___ mg/ml Date ___ Time ___ Int. ___	AN-162	Tape	Cisatracurium Muscle Relaxants (Color: Fl. Red)	LAN-1	Label	LAN-162	Label	Atropine 0.4 mg/ml Date ___ Time ___ Int. ___	AN-1D04	Tape	Atropine 0.4 mg/ml Anticholinergic Agents (Color: Green)	Cisatracurium 2 mg/ml Date ___ Time ___ Int. ___	AN-162D2	Tape	Cisatracurium 2 mg/ml Muscle Relaxants (Color: Fl. Red)	LAN-1D04	Label	LAN-162D2	Label	Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)	Curare ___ mg/ml Date ___ Time ___ Int. ___	AN-2	Tape	Curare Muscle Relaxants (Color: Fl. Red)																																
Anectine ___ mg/ml Date ___ Time ___ Int. ___	AN-33	Tape	Anectine Muscle Relaxants (Color: Fl. Red/Black)	Bupivacaine ___ % Date ___ Time ___ Int. ___	AN-42PC	Tape	Bupivacaine ___ % Local Anesthetics (Color: Gray)																																																																																																																
	LAN-33	Label			LAN-42PC	Label		Antibiotic ___ mg/ml Date ___ Time ___ Int. ___	AN-35	Tape	Antibiotic Miscellaneous (Color: White)	Bupivacaine 0.5% Date ___ Time ___ Int. ___	AN-42D05	Tape	Bupivacaine 0.5% Local Anesthetics (Color: Gray)	LAN-35	Label	LAN-42D05	Label	ATRAcurium ___ mg/ml Date ___ Time ___ Int. ___	AN-39	Tape	ATRAcurium Muscle Relaxants (Color: Fl. Red)	Bupivacaine 0.25% 2.5 mg/ml Date ___ Time ___ Int. ___	AN-42D25	Tape	Bupivacaine 0.25% 2.5 mg/ml Local Anesthetics (Color: Gray)	LAN-39	Label	LAN-42D25	Label	ATRAcurium 10 mg/ml Date ___ Time ___ Int. ___	AN-39D10	Tape	ATRAcurium 10 mg/ml Muscle Relaxants (Color: Fl. Red)	Calcium Chloride ___ mg/ml Date ___ Time ___ Int. ___	AN-110	Tape	Calcium Chloride Miscellaneous (Color: White)	LAN-39D10	Label	LAN-110	Label	Atropine ___ mg/ml Date ___ Time ___ Int. ___	AN-1	Tape	Atropine Anticholinergic Agents (Color: Green)	Cisatracurium ___ mg/ml Date ___ Time ___ Int. ___	AN-162	Tape	Cisatracurium Muscle Relaxants (Color: Fl. Red)	LAN-1	Label	LAN-162	Label	Atropine 0.4 mg/ml Date ___ Time ___ Int. ___	AN-1D04	Tape	Atropine 0.4 mg/ml Anticholinergic Agents (Color: Green)	Cisatracurium 2 mg/ml Date ___ Time ___ Int. ___	AN-162D2	Tape	Cisatracurium 2 mg/ml Muscle Relaxants (Color: Fl. Red)	LAN-1D04	Label	LAN-162D2	Label	Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)	Curare ___ mg/ml Date ___ Time ___ Int. ___	AN-2	Tape	Curare Muscle Relaxants (Color: Fl. Red)																																												
Antibiotic ___ mg/ml Date ___ Time ___ Int. ___	AN-35	Tape	Antibiotic Miscellaneous (Color: White)	Bupivacaine 0.5% Date ___ Time ___ Int. ___	AN-42D05	Tape	Bupivacaine 0.5% Local Anesthetics (Color: Gray)																																																																																																																
	LAN-35	Label			LAN-42D05	Label		ATRAcurium ___ mg/ml Date ___ Time ___ Int. ___	AN-39	Tape	ATRAcurium Muscle Relaxants (Color: Fl. Red)	Bupivacaine 0.25% 2.5 mg/ml Date ___ Time ___ Int. ___	AN-42D25	Tape	Bupivacaine 0.25% 2.5 mg/ml Local Anesthetics (Color: Gray)	LAN-39	Label	LAN-42D25	Label	ATRAcurium 10 mg/ml Date ___ Time ___ Int. ___	AN-39D10	Tape	ATRAcurium 10 mg/ml Muscle Relaxants (Color: Fl. Red)	Calcium Chloride ___ mg/ml Date ___ Time ___ Int. ___	AN-110	Tape	Calcium Chloride Miscellaneous (Color: White)	LAN-39D10	Label	LAN-110	Label	Atropine ___ mg/ml Date ___ Time ___ Int. ___	AN-1	Tape	Atropine Anticholinergic Agents (Color: Green)	Cisatracurium ___ mg/ml Date ___ Time ___ Int. ___	AN-162	Tape	Cisatracurium Muscle Relaxants (Color: Fl. Red)	LAN-1	Label	LAN-162	Label	Atropine 0.4 mg/ml Date ___ Time ___ Int. ___	AN-1D04	Tape	Atropine 0.4 mg/ml Anticholinergic Agents (Color: Green)	Cisatracurium 2 mg/ml Date ___ Time ___ Int. ___	AN-162D2	Tape	Cisatracurium 2 mg/ml Muscle Relaxants (Color: Fl. Red)	LAN-1D04	Label	LAN-162D2	Label	Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)	Curare ___ mg/ml Date ___ Time ___ Int. ___	AN-2	Tape	Curare Muscle Relaxants (Color: Fl. Red)																																																								
ATRAcurium ___ mg/ml Date ___ Time ___ Int. ___	AN-39	Tape	ATRAcurium Muscle Relaxants (Color: Fl. Red)	Bupivacaine 0.25% 2.5 mg/ml Date ___ Time ___ Int. ___	AN-42D25	Tape	Bupivacaine 0.25% 2.5 mg/ml Local Anesthetics (Color: Gray)																																																																																																																
	LAN-39	Label			LAN-42D25	Label		ATRAcurium 10 mg/ml Date ___ Time ___ Int. ___	AN-39D10	Tape	ATRAcurium 10 mg/ml Muscle Relaxants (Color: Fl. Red)	Calcium Chloride ___ mg/ml Date ___ Time ___ Int. ___	AN-110	Tape	Calcium Chloride Miscellaneous (Color: White)	LAN-39D10	Label	LAN-110	Label	Atropine ___ mg/ml Date ___ Time ___ Int. ___	AN-1	Tape	Atropine Anticholinergic Agents (Color: Green)	Cisatracurium ___ mg/ml Date ___ Time ___ Int. ___	AN-162	Tape	Cisatracurium Muscle Relaxants (Color: Fl. Red)	LAN-1	Label	LAN-162	Label	Atropine 0.4 mg/ml Date ___ Time ___ Int. ___	AN-1D04	Tape	Atropine 0.4 mg/ml Anticholinergic Agents (Color: Green)	Cisatracurium 2 mg/ml Date ___ Time ___ Int. ___	AN-162D2	Tape	Cisatracurium 2 mg/ml Muscle Relaxants (Color: Fl. Red)	LAN-1D04	Label	LAN-162D2	Label	Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)	Curare ___ mg/ml Date ___ Time ___ Int. ___	AN-2	Tape	Curare Muscle Relaxants (Color: Fl. Red)																																																																				
ATRAcurium 10 mg/ml Date ___ Time ___ Int. ___	AN-39D10	Tape	ATRAcurium 10 mg/ml Muscle Relaxants (Color: Fl. Red)	Calcium Chloride ___ mg/ml Date ___ Time ___ Int. ___	AN-110	Tape	Calcium Chloride Miscellaneous (Color: White)																																																																																																																
	LAN-39D10	Label			LAN-110	Label		Atropine ___ mg/ml Date ___ Time ___ Int. ___	AN-1	Tape	Atropine Anticholinergic Agents (Color: Green)	Cisatracurium ___ mg/ml Date ___ Time ___ Int. ___	AN-162	Tape	Cisatracurium Muscle Relaxants (Color: Fl. Red)	LAN-1	Label	LAN-162	Label	Atropine 0.4 mg/ml Date ___ Time ___ Int. ___	AN-1D04	Tape	Atropine 0.4 mg/ml Anticholinergic Agents (Color: Green)	Cisatracurium 2 mg/ml Date ___ Time ___ Int. ___	AN-162D2	Tape	Cisatracurium 2 mg/ml Muscle Relaxants (Color: Fl. Red)	LAN-1D04	Label	LAN-162D2	Label	Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)	Curare ___ mg/ml Date ___ Time ___ Int. ___	AN-2	Tape	Curare Muscle Relaxants (Color: Fl. Red)																																																																																
Atropine ___ mg/ml Date ___ Time ___ Int. ___	AN-1	Tape	Atropine Anticholinergic Agents (Color: Green)	Cisatracurium ___ mg/ml Date ___ Time ___ Int. ___	AN-162	Tape	Cisatracurium Muscle Relaxants (Color: Fl. Red)																																																																																																																
	LAN-1	Label			LAN-162	Label		Atropine 0.4 mg/ml Date ___ Time ___ Int. ___	AN-1D04	Tape	Atropine 0.4 mg/ml Anticholinergic Agents (Color: Green)	Cisatracurium 2 mg/ml Date ___ Time ___ Int. ___	AN-162D2	Tape	Cisatracurium 2 mg/ml Muscle Relaxants (Color: Fl. Red)	LAN-1D04	Label	LAN-162D2	Label	Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)	Curare ___ mg/ml Date ___ Time ___ Int. ___	AN-2	Tape	Curare Muscle Relaxants (Color: Fl. Red)																																																																																												
Atropine 0.4 mg/ml Date ___ Time ___ Int. ___	AN-1D04	Tape	Atropine 0.4 mg/ml Anticholinergic Agents (Color: Green)	Cisatracurium 2 mg/ml Date ___ Time ___ Int. ___	AN-162D2	Tape	Cisatracurium 2 mg/ml Muscle Relaxants (Color: Fl. Red)																																																																																																																
	LAN-1D04	Label			LAN-162D2	Label		Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)	Curare ___ mg/ml Date ___ Time ___ Int. ___	AN-2	Tape	Curare Muscle Relaxants (Color: Fl. Red)																																																																																																								
Atropine 1 mg/ml Date ___ Time ___ Int. ___	LAN-1D1	Label	Atropine 1 mg/ml Anticholinergic Agents (Color: Green)	Curare ___ mg/ml Date ___ Time ___ Int. ___	AN-2	Tape	Curare Muscle Relaxants (Color: Fl. Red)																																																																																																																

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

ANESTHESIA

ANESTHESIA LABELS & TAPE ALPHABETICALLY LISTED

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description	Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description
Decadron _____ mg/ml Date _____ Time _____ Int. _____	AN-203	Tape	Decadron Miscellaneous (Color: White)	 Edrophonium _____ mg/ml Date _____ Time _____ Int. _____	ANS-50	Tape	Edrophonium Relaxant Antagonists (Color: Red/White)
	LAN-203	Label		 Edrophonium 10 mg/ml Date _____ Time _____ Int. _____	LANS-50D10	Label	
Demerol _____ mg/ml Date _____ Time _____ Int. _____	AN-46	Tape	Demerol Narcotics (Color: Lt. Blue)	EPHEDrine _____ mg/ml Date _____ Time _____ Int. _____	AN-5	Tape	EPHEDrine Vasopressors (Color: Violet)
	LAN-46	Label		EPHEDrine 5 mg/ml Date _____ Time _____ Int. _____	AN-5D5	Tape	
Dexamethasone _____ mg/ml Date _____ Time _____ Int. _____	AN-132	Tape	Dexamethasone Miscellaneous (Color: White)	EPHEDrine 10 mg/ml Date _____ Time _____ Int. _____	AN-5D10	Tape	EPHEDrine 10 mg/ml Vasopressors (Color: Violet)
	LAN-132	Label		EPINEPHrine _____ mcg/ml Date _____ Time _____ Int. _____	AN-6	Tape	
Dexamethasone 4 mg/ml Date _____ Time _____ Int. _____	AN-132D4	Tape	Dexamethasone 4 mg/ml Miscellaneous (Color: White)	EPINEPHrine 1 mg/ml Date _____ Time _____ Int. _____	AN-6D1	Tape	EPINEPHrine 1 mg/ml Vasopressors (Color: Violet/Black)
	LAN-132D4	Label		EPINEPHrine 10 mcg/ml Date _____ Time _____ Int. _____	AN-6D10	Tape	
Diazepam _____ mg/ml Date _____ Time _____ Int. _____	AN-3	Tape	Diazepam Tranquilizers (Color: Orange)	EPINEPHrine 1:200,000 Date _____ Time _____ Int. _____	LAN-6D12	Label	EPINEPHrine 1:200,000 Vasopressors (Color: Violet/Black)
	LAN-3	Label		EPINEPHrine 1:400,000 Date _____ Time _____ Int. _____	LAN-6D14	Label	
Dilaudid _____ mg/ml Date _____ Time _____ Int. _____	AN-200	Tape	Dilaudid Narcotics (Color: Lt. Blue)	 Esmolol _____ mg/ml Date _____ Time _____ Int. _____	ANS-20	Tape	Esmolol Hypotensive Agents (Color: Violet/White)
	LAN-200	Label		 Esmolol 10 mg/ml Date _____ Time _____ Int. _____	ANS-20D10	Tape	
Dilaudid 0.5 mg/ml Date _____ Time _____ Int. _____	AN-200D05	Tape	Dilaudid 0.5 mg/ml Narcotics (Color: Lt. Blue)	Etomidate _____ mg/ml Date _____ Time _____ Int. _____	AN-51	Tape	Etomidate Induction Agents (Color: Yellow)
	LAN-200D05	Label		Etomidate 2 mg/ml Date _____ Time _____ Int. _____	AN-51D2	Tape	
Diphenhydramine 50 mg/ml Date _____ Time _____ Int. _____	AN-210D50	Tape	Diphenhydramine 50 mg/ml Miscellaneous (Color: White)	Fentanyl _____ mcg/ml Date _____ Time _____ Int. _____	AN-7	Tape	Fentanyl Narcotics (Color: Lt. Blue)
	LAN-210D50	Label		Fentanyl 10 mcg/ml Date _____ Time _____ Int. _____	AN-7D10	Tape	
INJECTION DIPRIVAN® propofol 10 mg/ml Date _____ Time _____ Int. _____	AN-155	Tape	Diprivan Propofol Induction Agents (Color: Yellow)				
	LAN-155	Label					
DOBUTamine _____ mg/ml Date _____ Time _____ Int. _____	AN-113	Tape	DOBUTamine Vasopressors (Color: Violet)				
	LAN-113	Label					
DOPamine _____ mg/ml Date _____ Time _____ Int. _____	AN-47	Tape	DOPamine Vasopressors (Color: Violet)				
	LAN-47	Label					
DOXAcurium _____ mg/ml Date _____ Time _____ Int. _____	AN-99	Tape	DOXAcurium Muscle Relaxants (Color: Fl. Red)				
Droperidol _____ mg/ml Date _____ Time _____ Int. _____	AN-4	Tape	Droperidol Major Tranquilizers (Color: Salmon)				
	LAN-4	Label					
Droperidol 2.5 mg/ml Date _____ Time _____ Int. _____	AN-4D25	Tape	Droperidol 2.5 mg/ml Major Tranquilizers (Color: Salmon)				
	LAN-4D25	Label					
Drug _____ mg/ml Date _____ Time _____ Int. _____	AN-14	Tape	Drug Miscellaneous (Color: White)				
	LAN-14	Label					
Duramorph 0.5 mg/ml Date _____ Time _____ Int. _____	LAN-211D05	Label	Duramorph 0.5 mg/ml Narcotics (Color: Lt. Blue)				

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

ANESTHESIA

ANESTHESIA

ANESTHESIA LABELS & TAPE ALPHABETICALLY LISTED

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description
Fentanyl 50 mcg/ml Date _____ Time _____ Int. _____	AN-7D50	Tape	Fentanyl 50 mcg/ml Narcotics (Color: Lt. Blue)
	LAN-7D50	Label	
Flumazenil _____ mg/ml Date _____ Time _____ Int. _____	ANS-18	Tape	Flumazenil Relaxant Antagonists (Color: Red/White)
	LANS-18	Label	
Flumazenil 0.1 mg/ml Date _____ Time _____ Int. _____	ANS-18D01	Tape	Flumazenil 0.1 mg/ml Relaxant Antagonists (Color: Red/White)
	LANS-18D01	Label	
Furosemide _____ mg/ml Date _____ Time _____ Int. _____	AN-134	Tape	Furosemide Miscellaneous (Color: White)
	LAN-134	Label	
Furosemide 10 mg/ml Date _____ Time _____ Int. _____	AN-134D10	Tape	Furosemide 10 mg/ml Miscellaneous (Color: White)
	LAN-134D10	Label	
Glycopyrrolate _____ mg/ml Date _____ Time _____ Int. _____	AN-9	Tape	Glycopyrrolate Anticholinergic Agents (Color: Green)
	LAN-9	Label	
Glycopyrrolate 0.2 mg/ml Date _____ Time _____ Int. _____	AN-9D02	Tape	Glycopyrrolate 0.2 mg/ml Anticholinergic Agents (Color: Green)
	LAN-9D02	Label	
Heparin _____ mg/ml Date _____ Time _____ Int. _____	AN-10	Tape	Heparin Miscellaneous (Color: White)
	LAN-10	Label	
Heparin _____ Units/ml Date _____ Time _____ Int. _____	AN-10U	Tape	Heparin (Units) Miscellaneous (Color: White)
	LAN-10U	Label	
Heparin 1000 units/ml Date _____ Time _____ Int. _____	AN-10D1000	Tape	Heparin 1000 units/ml Miscellaneous (Color: White)
	LAN-10D1000	Label	
HydrALAZine _____ mg/ml Date _____ Time _____ Int. _____	ANS-118	Tape	HydrALAZine Hypotensive Agents (Color: Violet/White)
	LANS-118	Label	
HydrALAZine 20 mg/ml Date _____ Time _____ Int. _____	ANS-118D20	Label	HydrALAZine 20 mg/ml Hypotensive Agents (Color: Violet/White)
	LANS-118D20	Label	
HydrALAZine 4 mg/ml Date _____ Time _____ Int. _____	LANS-118D4	Label	HydrALAZine 4 mg/ml Hypotensive Agents (Color: Violet/White)
HYDRomorphone _____ mg/ml Date _____ Time _____ Int. _____	AN-16	Tape	HYDRomorphone Narcotics (Color: Lt. Blue)
	LAN-16	Label	
HYDRomorphone 0.2 mg/ml Date _____ Time _____ Int. _____	AN-16D02	Tape	HYDRomorphone 0.2mg/ml Narcotics (Color: Lt. Blue)
	LAN-16D02	Label	
Inapsine _____ mg/ml Date _____ Time _____ Int. _____	AN-54	Tape	Inapsine Major Tranquilizers (Color: Salmon)
	LAN-54	Label	

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description
Insulin 100 units/ml Date _____ Time _____ Int. _____	AN-207D100	Tape	Insulin 100 units/ml Miscellaneous (Color: White)
	LAN-207D100	Label	
Isuprel _____ mg/ml Date _____ Time _____ Int. _____	AN-57	Tape	Isuprel Vasopressors (Color: Violet)
	LAN-57	Label	
Ketalar _____ mg/ml Date _____ Time _____ Int. _____	AN-59	Tape	Ketalar Induction Agents (Color: Yellow)
Ketamine _____ mg/ml Date _____ Time _____ Int. _____	AN-60	Tape	Ketamine Induction Agents (Color: Yellow)
	LAN-60	Label	
Ketamine 10 mg/ml Date _____ Time _____ Int. _____	AN-60D10	Tape	Ketamine 10 mg/ml Induction Agents (Color: Yellow)
	LAN-60D10	Label	
Ketamine 50 mg/ml Date _____ Time _____ Int. _____	LAN-60D50	Label	Ketamine 50 mg/ml Induction Agents (Color: Yellow)
Labetalol _____ mg/ml Date _____ Time _____ Int. _____	ANS-117	Tape	Labetalol Hypotensive Agents (Color: Violet/White)
	LANS-117	Label	
Labetalol 5 mg/ml Date _____ Time _____ Int. _____	ANS-117D5	Tape	Labetalol 5 mg/ml Hypotensive Agents (Color: Violet/White)
	LANS-117D5	Label	
Levophed _____ mg/ml Date _____ Time _____ Int. _____	AN-62	Tape	Levophed Vasopressors (Color: Violet)
	LAN-62	Label	
Lidocaine _____ mg/ml Date _____ Time _____ Int. _____	AN-11	Tape	Lidocaine Local Anesthetics (Color: Gray)
	LAN-11	Label	
Lidocaine _____ % Date _____ Time _____ Int. _____	AN-11PC	Tape	Lidocaine _____ % Local Anesthetics (Color: Gray)
	LAN-11PC	Label	
Lidocaine 1% Date _____ Time _____ Int. _____	AN-11D1	Tape	Lidocaine 1% Local Anesthetics (Color: Gray)
Lidocaine 1% 10 mg/ml Date _____ Time _____ Int. _____	AN-11D10	Tape	Lidocaine 1% 10 mg/ml Local Anesthetics (Color: Gray)
	LAN-11D10	Label	
Lidocaine 2% 20 mg/ml Date _____ Time _____ Int. _____	AN-11D20	Tape	Lidocaine 2% 20 mg/ml Local Anesthetics (Color: Gray)
	LAN-11D20	Label	
Lispro 100 units/ml Date _____ Time _____ Int. _____	AN-208D100	Tape	Lispro 100 units/ml Miscellaneous (Color: White)
	LAN-208D100	Label	
Lorazepam _____ mg/ml Date _____ Time _____ Int. _____	AN-26	Tape	Lorazepam Tranquilizers (Color: Orange)
	LAN-26	Label	

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

ANESTHESIA

ANESTHESIA LABELS & TAPE ALPHABETICALLY LISTED

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description
Magnesium 0.5 mg/ml Date _____ Time _____ Int. _____	AN-209D05	Tape	Magnesium 0.5 mg/ml Miscellaneous (Color: White)
	LAN-209D05	Label	
Marcaine _____ mg/ml Date _____ Time _____ Int. _____	AN-64	Tape	Marcaine Local Anesthetics (Color: Gray)
	LAN-64	Label	
Meperidine _____ mg/ml Date _____ Time _____ Int. _____	AN-12	Tape	Meperidine Narcotics (Color: Lt. Blue)
	LAN-12	Label	
Meperidine 10 mg/ml Date _____ Time _____ Int. _____	AN-12D10	Tape	Meperidine 10 mg/ml Narcotics (Color: Lt. Blue)
	LAN-12D10	Label	
Mepivacaine _____ mg/ml Date _____ Time _____ Int. _____	AN-34	Tape	Mepivacaine Local Anesthetics (Color: Gray)
	LAN-34	Label	
Methohexital 10 mg/ml Date _____ Time _____ Int. _____	AN-68	Tape	Methohexital 10 mg/ml Inductions Agents (Color: Yellow)
	LAN-68	Label	
Metoclopramide _____ mg/ml Date _____ Time _____ Int. _____	AN-17	Tape	Metoclopramide Miscellaneous (Color: White)
	LAN-17	Label	
Metoclopramide 5 mg/ml Date _____ Time _____ Int. _____	LAN-17D5	Label	Metoclopramide 5 mg/ml Miscellaneous (Color: White)
Metoclopramide 10 mg/ml Date _____ Time _____ Int. _____	AN-17D10	Tape	Metoclopramide 10 mg/ml Miscellaneous (Color: White)
	LAN-17D10	Label	
Metocurine _____ mg/ml Date _____ Time _____ Int. _____	AN-30	Tape	Metocurine Muscle Relaxants (Color: Fl. Red)
Metoprolol _____ mg/ml Date _____ Time _____ Int. _____	ANS-49	Tape	Metoprolol Hypotensive Agents (Color: Violet/White)
	LANS-49	Label	
Metoprolol 1 mg/ml Date _____ Time _____ Int. _____	ANS-49D1	Tape	Metoprolol 1 mg/ml Hypotensive Agents (Color: Violet/White)
	LANS-49D1	Label	
Midazolam _____ mg/ml Date _____ Time _____ Int. _____	AN-120	Tape	Midazolam Tranquilizers (Color: Orange)
	LAN-120	Label	
Midazolam 1 mg/ml Date _____ Time _____ Int. _____	AN-120D1	Tape	Midazolam 1 mg/ml Tranquilizers (Color: Orange)
	LAN-120D1	Label	
Milrinone _____ mg/ml Date _____ Time _____ Int. _____	AN-215	Tape	Milrinone Miscellaneous (Color: White)
	LAN-215	Label	
MIVAcurium _____ mg/ml Date _____ Time _____ Int. _____	AN-160	Tape	MIVAcurium Muscle Relaxants (Color: Fl. Red)
	LAN-160	Label	

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description
Morphine _____ mg/ml Date _____ Time _____ Int. _____	AN-13	Tape	Morphine Narcotics (Color: Lt. Blue)
	LAN-13	Label	
Morphine 1 mg/ml Date _____ Time _____ Int. _____	AN-13D1	Tape	Morphine 1 mg/ml Narcotics (Color: Lt. Blue)
	LAN-13D1	Label	
Morphine PF 0.5 mg/ml Date _____ Time _____ Int. _____	AN-13D05	Tape	Morphine PF 0.5 mg/ml Narcotics (Color: Lt. Blue)
	LAN-13D05	Label	
Morphine PF 1 mg/ml Date _____ Time _____ Int. _____	AN-13PFD1	Tape	Morphine PF 1 mg/ml Narcotics (Color: Lt. Blue)
	LAN-13PFD1	Label	
Morphine S 1 mg/ml Date _____ Time _____ Int. _____	LAN-13SD1	Label	Morphine Sulfate 1 mg/ml Narcotics (Color: Lt. Blue)
Naloxone _____ mg/ml Date _____ Time _____ Int. _____	ANS-14	Tape	Naloxone Narcotics Antagonists (Color: Lt. Blue/White)
	LANS-14	Label	
Naloxone 0.04 mg/ml Date _____ Time _____ Int. _____	ANS-14D004	Tape	Naloxone 0.04 mg/ml Narcotics Antagonists (Color: Lt. Blue/White)
	LANS-14D004	Label	
Narcan _____ mg/ml Date _____ Time _____ Int. _____	ANS-71	Tape	Narcan Narcotics Antagonists (Color: Lt. Blue/White)
	LANS-71	Label	
Neostigmine _____ mg/ml Date _____ Time _____ Int. _____	ANS-15	Tape	Neostigmine Relaxant Antagonists (Color: Red/White)
	LANS-15	Label	
Neostigmine 1 mg/ml Date _____ Time _____ Int. _____	ANS-15D1	Tape	Neostigmine 1 mg/ml Relaxant Antagonists (Color: Red/White)
	LANS-15D1	Label	
Neosyneprine _____ mg/ml Date _____ Time _____ Int. _____	AN-72	Tape	Neosyneprine Vasopressors (Color: Violet)
	LAN-72	Label	
Neosyneprine 100 mcg/ml Date _____ Time _____ Int. _____	AN-72D100	Tape	Neosyneprine 100 mcg/ml Vasopressors (Color: Violet)
	LAN-72D100	Label	
Nimbex _____ mg/ml Date _____ Time _____ Int. _____	AN-25	Tape	Nimbex Muscle Relaxants (Color: Fl. Red)
	LAN-25	Label	
Nitroglycerine _____ mg/ml Date _____ Time _____ Int. _____	ANS-16	Tape	Nitroglycerine Hypotensive Agents (Color: Violet/White)
	LANS-16	Label	
Nitroprusside _____ mg/ml Date _____ Time _____ Int. _____	ANS-17	Tape	Nitroprusside Hypotensive Agents (Color: Violet/White)
	LANS-17	Label	
Norepinephrine _____ mg/ml Date _____ Time _____ Int. _____	AN-75	Tape	Norepinephrine Vasopressors (Color: Violet)
	LAN-75	Label	

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

ANESTHESIA

ANESTHESIA LABELS & TAPE ALPHABETICALLY LISTED

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description	Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description
Norcuron _____ mg/ml Date _____ Time _____ Int. _____	AN-150	Tape	Norcuron Muscle Relaxants (Color: Fl. Red)	Pentothal 25 mg/ml Date _____ Time _____ Int. _____	AN-79	Tape	Pentothal 25 mg/ml Inactive Agents (Color: Yellow)
	LAN-150	Label			LAN-79	Label	
Normal Saline _____ mg/ml Date _____ Time _____ Int. _____	AN-22	Tape	Normal Saline Miscellaneous (Color: White)	Pitocin _____ units/ml Date _____ Time _____ Int. _____	AN-23	Tape	Pitocin Miscellaneous (Color: White)
	LAN-22	Label			LAN-23	Label	
Nubain _____ mcg/ml Date _____ Time _____ Int. _____	AN-28	Tape	Nubain Narcotics (Color: Lt. Blue)	Precedex _____ mg/ml Date _____ Time _____ Int. _____	AN-24	Tape	Precedex Induction Agents (Color: Yellow)
	LAN-28	Label			LAN-24	Label	
Ondansetron _____ mg/ml Date _____ Time _____ Int. _____	AN-55	Tape	Ondansetron Miscellaneous (Color: White)	Promethazine 2.5 mg/ml Date _____ Time _____ Int. _____	AN-206D25	Tape	Promethazine 2.5 mg/ml Miscellaneous (Color: White)
	LAN-55	Label			LAN-206D25	Label	
Ondansetron 2 mg/ml Date _____ Time _____ Int. _____	AN-55D2	Tape	Ondansetron 2 mg/ml Miscellaneous (Color: White)	Propofol _____ mg/ml Date _____ Time _____ Int. _____	AN-27	Tape	Propofol Induction Agents (Color: Yellow)
	LAN-55D2	Label			LAN-27	Label	
Ondansetron 4 mg/ml Date _____ Time _____ Int. _____	LAN-55D4	Label	Ondansetron 4 mg/ml Miscellaneous (Color: White)	Propofol 10 mg/ml Date _____ Time _____ Int. _____	AN-27D10	Tape	Propofol 10 mg/ml Induction Agents (Color: Yellow)
				LAN-27D10	Label		
Oxytocin _____ mg/ml Date _____ Time _____ Int. _____	AN-76	Tape	Oxytocin Miscellaneous (Color: White)	Prostigmine _____ mg/ml Date _____ Time _____ Int. _____	ANS-84	Tape	Prostigmine Relaxant Antagonists (Color: Red/White)
	LAN-76	Label					
Oxytocin 10 units/ml Date _____ Time _____ Int. _____	AN-76D10	Tape	Oxytocin 10 units/ml Miscellaneous (Color: White)	Protamine _____ mg/ml Date _____ Time _____ Int. _____	AN-19	Tape	Protamine Miscellaneous (Color: White)
	LAN-76D10	Label			LAN-19	Label	
PANcuronium _____ mg/ml Date _____ Time _____ Int. _____	AN-18	Tape	PANcuronium Muscle Relaxants (Color: Fl. Red)	Protamine 10 mg/ml Date _____ Time _____ Int. _____	AN-19D10	Tape	Protamine 10 mg/ml Miscellaneous (Color: White)
	LAN-18	Label			LAN-19D10	Label	
PANcuronium 1 mg/ml Date _____ Time _____ Int. _____	AN-18D1	Tape	PANcuronium 1 mg/ml Muscle Relaxants (Color: Fl. Red)	Reglan _____ mg/ml Date _____ Time _____ Int. _____	AN-65	Tape	Reglan Miscellaneous (Color: White)
	LAN-18D1	Label			LAN-65	Label	
Pavulon _____ mg/ml Date _____ Time _____ Int. _____	AN-77	Tape	Pavulon Muscle Relaxants (Color: Fl. Red)	Remifentanil _____ mcg/ml Date _____ Time _____ Int. _____	AN-38	Tape	Remifentanil Narcotics (Color: Lt. Blue)
				LAN-38	Label		
Phenylephrine _____ mg/ml Date _____ Time _____ Int. _____	AN-81	Tape	Phenylephrine Vasopressors (Color: Violet)	Remifentanil 50 mcg/ml Date _____ Time _____ Int. _____	AN-38D50	Tape	Remifentanil 50 mcg/ml Narcotics (Color: Lt. Blue)
	LAN-81	Label			LAN-38D50	Label	
Phenylephrine _____ mcg/ml Date _____ Time _____ Int. _____	AN-81MCG	Tape	Phenylephrine _____ mcg/ml Vasopressors (Color: Violet)	Robinul _____ mg/ml Date _____ Time _____ Int. _____	AN-88	Tape	Robinul Anticholinergic Agents (Color: Green)
	LAN-81MCG	Label			LAN-88	Label	
Phenylephrine 40 mcg/ml Date _____ Time _____ Int. _____	AN-81D40	Tape	Phenylephrine 40 mcg/ml Vasopressors (Color: Violet)	ROcuronium _____ mg/ml Date _____ Time _____ Int. _____	AN-161	Tape	ROcuronium Muscle Relaxants (Color: Fl. Red)
	LAN-81D40	Label			LAN-161	Label	
Phenylephrine 100 mcg/ml Date _____ Time _____ Int. _____	AN-81D100	Tape	Phenylephrine 100 mcg/ml Vasopressors (Color: Violet)	ROcuronium 10 mg/ml Date _____ Time _____ Int. _____	AN-161D10	Tape	ROcuronium 10 mg/ml Muscle Relaxants (Color: Fl. Red)
	LAN-81D100	Label			LAN-161D10	Label	
Pentobarbital _____ mg/ml Date _____ Time _____ Int. _____	AN-78	Tape	Pentobarbital Induction Agents (Color: Yellow)	Romazicon _____ mg/ml Date _____ Time _____ Int. _____	ANS-19	Tape	Romazicon Relaxant Antagonists (Color: Red/White)
				LAN-19	Label		

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

ANESTHESIA

ANESTHESIA LABELS & TAPE ALPHABETICALLY LISTED

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description
Ropivacaine _____ mg/ml Date _____ Time _____ Int. _____	AN-32	Tape	Ropivacaine Local Anesthetics (Color: Gray)
	LAN-32	Label	
Ropivacaine _____ % Date _____ Time _____ Int. _____	AN-32PC	Tape	Ropivacaine _____ % Local Anesthetics (Color: Gray)
	LAN-32PC	Label	
Ropivacaine 0.25% Date _____ Time _____ Int. _____	LAN-32D25	Label	Ropivacaine 0.25% Local Anesthetics (Color: Gray)
Ropivacaine 0.5% Date _____ Time _____ Int. _____	LAN-32D05	Label	Ropivacaine 0.5% Local Anesthetics (Color: Gray)
Sodium Bicarbonate _____ mg/ml Date _____ Time _____ Int. _____	AN-89	Tape	Sodium Bicarbonate Miscellaneous (Color: White)
	LAN-89	Label	
Sodium Chloride _____ mg/ml Date _____ Time _____ Int. _____	AN-90	Tape	Sodium Chloride Miscellaneous (Color: White)
	LAN-90	Label	
Stadol _____ mg/ml Date _____ Time _____ Int. _____	AN-205	Tape	Stadol Narcotics (Color: Lt. Blue)
Sublimaze _____ mcg/ml Date _____ Time _____ Int. _____	AN-93	Tape	Sublimaze Narcotics (Color: Lt. Blue)
	LAN-93	Label	
Succinylcholine _____ mg/ml Date _____ Time _____ Int. _____	AN-20	Tape	Succinylcholine Muscle Relaxants (Color: Fl. Red/Black)
	LAN-20	Label	
Succinylcholine 20 mg/ml Date _____ Time _____ Int. _____	AN-20D20	Tape	Succinylcholine 20 mg/ml Muscle Relaxants (Color: Fl. Red/Black)
	LAN-20D20	Label	
Sufenta _____ mcg/ml Date _____ Time _____ Int. _____	AN-147	Tape	Sufenta Narcotics (Color: Lt. Blue)
	LAN-147	Label	
Sufentanil _____ mcg/ml Date _____ Time _____ Int. _____	AN-122	Tape	Sufentanil Narcotics (Color: Lt. Blue)
	LAN-122	Label	
Sufentanil 50 mcg/ml Date _____ Time _____ Int. _____	AN-122D50	Tape	Sufentanil 50 mcg/ml Narcotics (Color: Lt. Blue)
	LAN-122D50	Label	
Syntocinon _____ mg/ml Date _____ Time _____ Int. _____	AN-44	Tape	Syntocinon Miscellaneous (Color: White)
Tensilon _____ mg/ml Date _____ Time _____ Int. _____	ANS-96	Tape	Tensilon Relaxant Antagonists (Color: Red/White)
Terbutaline 1 mg/ml Date _____ Time _____ Int. _____	AN-213D1	Tape	Terbutaline 1 mg/ml Miscellaneous (Color: White)
	LAN-213D1	Label	

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Description
Thiopental 2.5% Date _____ Time _____ Int. _____	AN-21	Tape	Thiopental 2.5% Induction Agents (Color: Yellow)
	LAN-21	Label	
Toradol _____ mg/ml Date _____ Time _____ Int. _____	AN-201	Tape	Toradol Miscellaneous (Color: White)
	LAN-201	Label	
TRAcrium _____ mg/ml Date _____ Time _____ Int. _____	AN-95	Tape	TRAcrium Muscle Relaxants (Color: Fl. Red)
Ultiva _____ mcg/ml Date _____ Time _____ Int. _____	AN-45	Tape	Ultiva Narcotics (Color: Lt. Blue)
Valium _____ mg/ml Date _____ Time _____ Int. _____	AN-101	Tape	Valium Tranquilizers (Color: Orange)
	LAN-101	Label	
Vasopressin _____ Units/ml Date _____ Time _____ Int. _____	AN-216	Tape	Vasopressin Miscellaneous (Color: White)
	LAN-216	Label	
VEcuronium _____ mg/ml Date _____ Time _____ Int. _____	AN-103	Tape	VEcuronium Muscle Relaxants (Color: Fl. Red)
	LAN-103	Label	
VEcuronium 1 mg/ml Date _____ Time _____ Int. _____	AN-103D1	Tape	VEcuronium 1 mg/ml Muscle Relaxants (Color: Fl. Red)
	LAN-103D1	Label	
Verapamil _____ mg/ml Date _____ Time _____ Int. _____	AN-50	Tape	Verapamil Miscellaneous (Color: White)
	LAN-50	Label	
Versed _____ mg/ml Date _____ Time _____ Int. _____	AN-149	Tape	Versed Tranquilizers (Color: Orange)
	LAN-149	Label	
Xylocaine _____ mg/ml Date _____ Time _____ Int. _____	AN-105	Tape	Xylocaine Local Anesthetics (Color: Gray)
	LAN-105	Label	
Zemuron _____ mg/ml Date _____ Time _____ Int. _____	AN-52	Tape	Zemuron Muscle Relaxants (Color: Fl. Red)
	LAN-52	Label	
Zofran _____ mg/ml Date _____ Time _____ Int. _____	AN-58	Tape	Zofran Miscellaneous (Color: White)
	LAN-58	Label	



For Anesthesia
Label Dispensers,
See page 86.

DPM-16

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

ANESTHESIA

ANESTHESIA LABELS & TAPE GUIDELINES

Anesthesia Labels & Tapes Custom


Don't see a label you need? Custom imprinted drug tapes and labels are also available. Each conforms to the ASTM Standard D4774-11 (2017) unless otherwise noted. This standard color codes each label by drug class.

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Minimum & Incremental Order Qty.
Anticholinergic Agents	ANAA	Tape	12 RLS (Color: Green)
	LANAA	Label	
Hypotensive Agents	ANHA	Tape	12 RLS (Color: Violet/White)
	LANHA	Label	
Induction Agents	ANIA	Tape	12 RLS (Color: Yellow)
	LANIA	Label	
Local Anesthetics	ANLA	Tape	12 RLS (Color: Gray)
	LANLA	Label	
Major Tranquilizers	ANMT	Tape	12 RLS (Color: Salmon)
	LAMNT	Label	
Miscellaneous	ANW	Tape	12 RLS (Color: White)
	LANW	Label	
Muscle Relaxants	ANMR	Tape	12 RLS (Color: Fl. Red)
	LANMR	Label	

Label Size: 1-1/2" x 1/2"	Product Number	Tape/Label	Minimum & Incremental Order Qty.
Narcotics	ANN	Tape	12 RLS (Color: Lt. Blue)
	LANN	Label	
Narcotics Antagonists	ANNA	Tape	12 RLS (Color: Lt. Blue/White)
	LANNA	Label	
Narcotics	ANCNT	Tape	12 RLS (Color: Lt. Blue/Salmon)
Major Tranquilizers	LANCNT	Label	
Relaxant Antagonists	ANRA	Tape	12 RLS (Color: Red/White)
	LANRA	Label	
Tranquilizers	ANT	Tape	12 RLS (Color: Orange)
	LANT	Label	
Vasopressors	ANV	Tape	12 RLS (Color: Violet)
	LANV	Label	

Anesthesia Tape and Label Dispenser

Designed specifically to dispense PDC Healthcare Anesthesia tapes and labels. These dispensers are compatible with most anesthesia crash carts. For more information regarding other tape and label dispensers, please see page 156.



Product Number	Description
DSP-HC10	Clear plastic dispenser can hold 10 label rolls. Equipped with 9 clear removable dividers that provide equally separated compartments. Size: (W) 8-1/4" x (H) 5" x (D) 5-1/2"
DPM-16	White metal dispenser can hold up to 16 rolls of 1/2" wide tape. Divided into 4 sections that are 2" wide. Size: (W) 12-1/4" x (H) 3-5/16" x (D) 5"
500P-LAN	Plastic dispenser can hold 15 rolls of Anesthesia tape or label products. Includes inserts, hanging bars and bottom suction cups for secure placement. Size: (W) 12-1/2" x (H) 3-1/2" x (D) 4-1/2"

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.