

Bar Code ID Solutions

PDC Healthcare developed the first Bar Code Wristband System in 1984. Today, more hospitals choose our Bar Code Wristbands than any other system in the world. Before choosing a system, consider PDC Healthcare's experience, expertise, and proven advantages. We offer the highest quality, most reliable Bar Code ID Solutions in healthcare today that meet important patient safety regulatory requirements, as listed below.

Choose from the following products:

- Thermal Printer Wristbands & Tags print in mere seconds through direct thermal printers.
- Laser Printer Wristbands & Labels provide excellent printability and scannability, while protecting printed information from smearing or bleeding.
- Insert Style Wristbands feature convenient pockets to protect patient data.

Our Bar Code ID Solutions provide medical-grade quality for maximum durability, safety, and security.

Bar Code ID Solutions help healthcare organizations meet applicable patient safety regulatory guidelines, including:

- **The Joint Commission National Patient Safety Goals:** Identify Patients Correctly, Improve Staff Communication, Use Medicines Safely, Prevent Infection, Identify Patient Safety Risks, and Prevent Mistakes in Surgery.
- **Centers for Medicare & Medicaid Services (CMS):** Stage 2 meaningful use requirements to use assistive technologies with eMAR.
- **Health Insurance Portability and Accountability Act (HIPAA):** to protect patients' medical records and other health information.
- **American Health Association (AHA):** to prevent the chance of lost or transferred data.
- **Five Rights of Patient Safety:** to ensure right patient, medication, dose, time, and method.



Table of Contents

Bar Code Direct Thermal Printer Wristbands	7
Mother/Infant Thermal Bar Code Tags	13
Bar Code Label Wristbands	14
Direct Thermal Label Patient ID Wristbands	15
Bar Code Super Soft Direct Thermal Wristbands	16
Inkjet Wristbands	17
Bar Code Color Laser Wristbands	18
Bar Code Color Laser Wristbands & Labels	19
Bar Code Laser Printable Labels	28
Bar Code Insert Wristbands	31



Easily verify printer and wristband compatibility with PDC Healthcare's Product-to-Printer Reference Guide. Visit: www.pdchealthcare.com/en-us/product-to-printer-guide.html

BAR CODE ID SOLUTIONS

BAR CODE DIRECT THERMAL PRINTER WRISTBANDS

CompuBand® Direct Thermal with Snap Closure

- Patented PermaPrint® surface provides best imprint quality, maximum resistance to solvents, and seamless printing
- Prints in seconds through direct thermal printers, which activate a heat-sensitive thermal material
- Durable, strong, lightweight and comfortable
- CompuBand® features a softer material. Ideal for patients with sensitive skin, including babies and the elderly
- Provides high-quality bar code for scanning and accurate patient ID
- Tamperproof, secure snap closure

Precision® CompuBand® Snap



Product Number	Description	Snaps/Bag
8615-00-PDJ	Single-Post Snaps	250

CompuBand® Adult Direct Thermal

- Soft CompuBand® features softer material on 7041S, 7043S, and 7049S
- Includes Single-Post Snaps (Also sold separately. See top right.)



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Box
7041-**-PDM	11-1/2" L x 1-1/8" W	4" L x 1-1/8" H	1.5"	In	6.48"	250/RL-2RLS
7041S-**-PDM			1.5"	In	5.85"	250/RL-2RLS
7043-**-PDH			1"	In	4.15"	100/RL-2RLS
7043S-**-PDH			1"	In	3.80"	100/RL-2RLS
7049S-**-PDO			1"	In	4.98"	175/RL-2RLS
7055-**-PDO††			1"	Out	3.80"	85/RL-4RLS
7340-11-PDM†			3"	Out	7.20"	250/RL-2RLS
7341-**-PDJ†††			3"	In	5.5"	125/RL-2RLS
7342-**-PDM			3"	Out	7.09"	250/RL-2RLS

Substitute ** for numbered color:



†7340 Style wristband is available in White ONLY. Indexing Bar located at perf is compatible with all Zebra Printers
 ††7055 Style wristband is not available in 13-Blue
 †††7341 Style wristband is available in White and Green ONLY.

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.

CompuBand® Adult/Pediatric Direct Thermal

- Soft CompuBand® features softer material on 7023S, 7024S and 7029S
- Includes Single-Post Snaps (Also sold separately. See top right.)

For Alert Snap options, see pages 59-60.



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Box
7020-11-PDM	11-1/2" L x 1-1/8" W	3-1/8" L x 1-1/8" H	3"	In	7.09"	250/RL-2RLS
7023-**-PDM	11-1/2" L x 3/4" W	3-1/8" L x 3/4" H	1.5"		6.58"	250/RL-2RLS
7023S-**-PDM	11-1/2" L x 3/4" W	3-1/8" L x 3/4" H	1.5"		5.85"	250/RL-2RLS
7024-**-PDH	11-1/2" L x 3/4" W	3-1/8" L x 3/4" H	1"		4.10"	100/RL-2RLS
7024S-**-PDH	11-1/2" L x 3/4" W	3-1/8" L x 3/4" H	1"		3.80"	100/RL-2RLS
7029S-**-PDO	11-1/2" L x 3/4" W	3-1/8" L x 3/4" H	1"		4.96"	175/RL-2RLS

Substitute ** for numbered color:



ALL WRISTBANDS ARE LATEX-FREE.

BAR CODE ID SOLUTIONS

BAR CODE DIRECT THERMAL PRINTER WRISTBANDS

CompuBand® Extra Wide Direct Thermal Wristbands

- Patented PermaPrint® surface provides best imprint quality, maximum resistance to solvents, and seamless printing
- Includes Single-Post Snaps (Also sold separately. See page 7.)



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Box
7050-11-PDM	11-1/2" L x 1 3/8" W	4" L x 1 3/8" H	1.5"	In	6.25"	250/RL-2RLS
7056-11-PDH			1"		4.00"	100/RL-2RLS

CompuBand® X Wristband

- Patented PermaPrint® surface provides best imprint quality, maximum resistance to solvents, and seamless printing
- "X-Hole" pattern allows you to print along the entire length of the band
- Includes Single Post Snap (Also sold separately. See page 7.)



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Box
7550-11-PDK	11-1/2" L x 1-1/8" W	11-1/2" L x 1-1/8" H	1.5"	In	5.23"	150/RL-2RLS

Conf-ID-ent™ Adult Direct Thermal Wristband with Snap Closure

- Durable - resistant to tearing
- Supports text, bar codes and graphics
- Excellent first time scan rates
- Flexible and comfortable design
- Available in 3 core sizes: 1", 1.5" and 3"
- Available wind direction: wound in or wound out
- Snaps included
- Snap closure

For Alert Snap options, see pages 59-60.



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Case
WBCA22T-1†	11-1/2"	3" x 1"	1.5"	Out	4.4"	80/RL-4RLS
WBCA23N-**	11-1/2"	3" x 1"	1.5"	In	5"	125/RL-2RLS
WBCA23T-**	11-1/2"	3" x 1"	1.5"	Out	6.5"	250/RL-2RLS
WBCA33N-**	11-1/2"	3" x 1"	3"	In	3.80"	125/RL-2 RLS
WBZC412T-1-***†	11"	3" x 1"	1"	Out	5"	152/RL-6RLS
WBZC435T†	11"	3" x 1"	3"	Out	8"	327/RL-6RLS
WBCA11T-**	11-1/2"	3.75" x 1.125"	1"	Out	4"	85/RL-4RLS
WBCA12N-**	11-1/2"	3.75" x 1.125"	1"	In	4.2"	100/RL-2RLS
WBCA25N-**	11-1/2"	3.75" x 1.125"	1.5"	In	6.5"	250/RL-2RLS

Substitute ** for numbered color: 1 2 3 4 5 6
 White Blue Yellow Green Red Orange

†White (Color options available with a 15/case minimum purchase.)

Conf-ID-ent™ Pediatric Direct Thermal Wristband with Snap Closure

- Durable – resistant to tearing
- Supports text, bar codes and graphics
- Excellent first time scan rates
- Flexible and comfortable design
- Available in 2 core sizes: 1" and 1.5"
- Available wind direction: wound in or wound out
- Closure snaps included
- Snap closure

For Alert Snap options, see pages 59–60.



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Case
WBCP12T-**†	10-1/2"	3" x 1"	1"	Out	4.5"	100/RL-2RLS
WBCP12N-**		2.5" x .875"	1.5"	In	5"	125/RL-2RLS
WBCP25N-**		2.5" x .875"	1.5"	In	6.5"	250/RL-2RLS

Substitute ** for numbered color:

1	2	3	4	5	6
White	Blue	Yellow	Green	Red	Orange

†White (Color options available with a 15/case minimum purchase.)

Precision® Infant Direct Thermal Wristbands

Precision® Infant Direct Thermal Wristbands provide positive patient identification for newborns using bar code technology to help ensure patient safety.

- Patented PermaPrint® surface provides best imprint quality, maximum resistance to solvents, and seamless printing
- Provides a flat surface for best bar code scanning results, eliminating curvature issues
- Improved direct thermal surface protects print & bar code from water, soap, and hand sanitizer
- Accommodates most thermal printers
- Easy to apply adhesive closure ensures safety and security
- Available in White or Snuggly Bears® pattern



Product Number	Size of Carrier	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Box
7701-11-PDL	8" L x 1-1/8" W	1-1/2" L x 3/4" H	1.5"	In	4.95"	200/RL-2RLS
7701-49-PDL						
7702-49-PDL						
7701-11-PDL†						
7701-49-PDL†						

7701-11 White	7701-49 Blue Snuggly Bears	7702-49 Pink Snuggly Bears	†7710-49 Snuggly Bears with Black Bar at Perf	†7710-11 White with Black Bar at Perf

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.

ALL WRISTBANDS ARE LATEX-FREE.

BAR CODE ID SOLUTIONS

BAR CODE ID SOLUTIONS

BAR CODE DIRECT THERMAL PRINTER WRISTBANDS

ScanBand® Direct Thermal with Adhesive Closure

- Patented PermaPrint® surface provides best imprint quality, maximum resistance to solvents, and seamless printing
- ScanBand® provides a high-quality bar code for scanning and accurate patient ID
- Supports text, linear and 2D bar codes, photos and graphics
- Tamper-evident adhesive closure

ScanBand® Adult Direct Thermal



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Box
7042-**-PDM*	11-1/2" L x 1-1/8" W	4-1/2" L x 1-1/8" H	1.5"	In	5.32"	250/RL-2RLS
7045-**-PDL†			1"	In	4.71"	200/RL-2RLS
7046-**-PDL††			1"	In	4.76"	200/RL-2RLS
7048-**-PDL†			1"	In	4.76"	200/RL-2RLS
7054-**-PDO*			1"	Out	3.15"	85/RL-4RLS
7057-11-PDK			1"	Out	4.28"	150/RL-2RLS
7058-11-PDO			1.5"	Out	5.02"	220/RL-2RLS
7240-**-PDM			1"	In	5.25"	250/RL-2RLS
7343-**-PDM			3"	Out	6.05"	250/RL-2RLS
7344-**-PDR†			3"	Out	7.88"	500/RL-2RLS
7142H-**-PDM††			1.5"	In	5.32"	250/RL-2RLS
7051-11-PDM***			3"	Out		250/RL-2RLS

Substitute ** for numbered color:



*7042, 7048 and 7054 Style wristbands are not available in 17-Orange

†Indexing Bar located at perf is compatible with all Zebra Printers

††With sensor hole

ScanBand® Adult/Pediatric Direct Thermal



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Box
7022-**-PDM	11-1/2" L x 3/4" W	4-1/2" L x 3/4" H	1.5"	In	5.23"	250/RL-2RLS
7052-**-PDH†			1"	In	3.46"	100/RL-2RLS
7028-**-PDM			3"	Out	6"	250/RL-2RLS
7030-11-PDL††*			1.5"	In	5"	200/RL-2RLS

Substitute ** for numbered color:



†7052 Style wristband is available in White, Blue and Red ONLY

††Indexing Bar located at perf is compatible with all Zebra Printers

*Without sensor hole

Conf-ID-ent™ Adult Direct Thermal Wristband with Adhesive Closure

- Durable – resistant to tearing
- Supports text, bar codes and graphics
- Excellent first-time scan rates
- Flexible and comfortable design
- Available in 3 core sizes: 1", 1.5" and 3"
- Available wind direction: wound in or wound out
- Adhesive closure



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Box
WBAA16N-**	11-1/2" L x 1-1/8" W	4-1/2" x 1-1/8"	1"	In	4.75"	200/RL-2RLS
WBAA23N-**			1.5"	In	5.32	200/RL-2RLS
WBAA35T-**			3"	Out	6.05"	250/RL-2RLS
WBAA11T-**			1"	Out	3.2"	85/RL-4RLS
WBAA25N-**			1.5"	In	5.32	250/RL-2RLS

Substitute ** for numbered color:



Conf-ID-ent™ Adult Direct Thermal Wristband with Adhesive Closure

- Durable – resistant to tearing
- Supports text, bar codes and graphics
- Excellent first time scan rates
- Flexible and comfortable design
- Available in 2 core sizes: 1" and 3"
- Indexing bar located at perf is compatible with all Zebra Printers



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Roll
WBZA112T-1†	11" L x 1" W	5-1/2" x 3/4"	1"	Out	5"	200/RL-6RLS
WBZA134T-1†	11" L x 1" W		3"		8"	450/RL-6RLS
WBZA212T-1†	11" L x 1-1/4" W		1"		5"	200/RL-6RLS
WBZA234T-1†	11" L x 1-1/4" W		3"		8"	450/RL-6RLS
WBZA312T-1†	11" L x 3/4" W		1"		5"	200/RL-6RLS

†White (Color options available with a 15/case minimum purchase.)

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.

ALL WRISTBANDS ARE LATEX-FREE.

BAR CODE ID SOLUTIONS

BAR CODE DIRECT THERMAL PRINTER WRISTBANDS

Conf-ID-ent™ Direct Thermal Pediatric Wristbands


- Durable - resistant to tearing
- Supports text, bar codes and graphics
- Excellent first-time scan rates
- Flexible and comfortable design
- Available in 2 core sizes: 1" and 3"
- Available in Pediatric and Infant

For a complete selection of wristbands, visit us online at www.pdchealthcare.com.



Product Number	Size of Wristband	Core Size	Wound	Outer Diameter	Wristbands/Case
Pediatric - WBZP113T-1	7" L x 1" W	1"	Out	5"	1,800
Pediatric - WBZP137T-1	7" L x 1" W	3"		8"	4,200
Pediatric - WBZP139T-1	7" L x .75" W	1"		5"	1,800
Pediatric - WBZP339T-1	7" L x .75" W	3"		8"	4,800

Conf-ID-ent™ Direct Thermal Infant Wristbands



Product Number	Size of Wristband	Core Size	Wound	Outer Diameter	Wristbands/Case
Infant - WBZB113T-**	6" L x 1" W	1"	Out	5"	350/RL-6RL/CS
Infant - WBZB123T-**	6" L x 1" W	1.5"		6"	350/RL-6RL/CS
Infant - WBZB130T-**	6" L x 1" W	3"		8"	800/RL-6RL/CS
Infant - WBZB316T-**	6" L x 3/4" W	1"		5"	350/RL-6RL/CS
Infant - WBZB330T-**	6" L x 3/4" W	3"		8"	800/RL-6RL/CS

Substitute ** for numbered color:




Conf-ID-ent™ Newborn Identification Tags


Conf-ID-ent™ newborn identification tags are ideal for bar coding your smallest patients. Their unique dual-sided design and adhesive-free strap features a flat surface that provides ease of scanning in any position without disturbing the patient. Conf-ID-ent™ Newborn Identification Tags feature super soft materials, making them gentle on your most sensitive patients' skin.



- Dual-sided for ease of scanning
- Adhesive-free strap allows scanning in any position, without disturbing the patient
- Rounded corners make it soft and gentle on sensitive patients' skin
- Lightweight and comfortable design
- Non-removable liner provides excellent durability
- Supports text, bar codes and photos
- Excellent first-time scan grades
- Direct thermal



Substitute ** for numbered color:



Product Number	Core Size	Outer Diameter	Wound	Wristbands/Case
WBTAG1XN-**	1"	5"	In	500/CS-2RLS/250
WBTAG1XT-**	1"	5"	Out	500/CS-2RLS/250
WBTAG15N-**	1.5"	5"	In	500/CS-2RLS/250
WBTAG15T-**	1.5"	5"	Out	500/CS-2RLS/250
WBTAG3XN-**	3"	8"	In	1,000/CS-2RLS/500
WBTAG3XT-**	3"	8"	Out	1,000/CS-2RLS/500

ALL WRISTBANDS ARE LATEX-FREE.


BAR CODE ID SOLUTIONS

MOTHER/INFANT DIRECT THERMAL BAR CODE TAGS

Precision® Mother/Infant Thermal Bar Code Tags

- Flat surface makes scanning bar codes easier and more accurate when "curvature" is an issue
- Eliminates the need of extra wristbands for added patient information
- Supports text, linear and 2D bar codes, photo ID and graphics
- Threads easily on patient wristband for positive patient identification
- Compatible with all PDC Healthcare Mother/Father/Baby wristbands*

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.



(2 Adult Tags and 2 Infant Tags)

Product Number	Information Area	Core Size	Wound	Sets/Roll	Rolls/Box
7600-11-PDL	Adult: 1-1/8" L x 1-1/4" H	15"	In	200	2
7601-11-PDL	Infant: 1-1/8" L x 5/8" H	3"			

Available in White ONLY

*Excludes PDC Healthcare Direct Thermal Bar Code and Tyvek® wristbands



BAR CODE ID SOLUTIONS

Precision® Infant Thermal Bar Code Tags



Product Number	Information Area	Core Size	Wound	Tags/Roll	Rolls/Box
7610-11-PDR	1-1/8" L x 5/8" H	15"	In	500	2
7612-11-PDR		3"			

Available in White ONLY

Precision® Pediatric Thermal Bar Code Tag

- Flat surface makes scanning bar codes easier and more accurate when "curvature" is an issue
- Pediatric style tags are compatible with all PDC Healthcare Mother/Father/Baby wristbands*
- Threads easily on patient wristband for positive patient identification
- In addition to patient identification, tags can be used to ID bottles, breast milk, personal items, etc.



Product Number	Information Area	Core Size	Wound	Tags per Roll	Rolls/Box
7630-11-PDR	1-1/2" L x 1/2" H	15"	In	500	2

Available in White ONLY

*Excludes PDC Healthcare Direct Thermal Bar Code and Tyvek® wristbands

ALL WRISTBANDS ARE LATEX-FREE.

BAR CODE ID SOLUTIONS

BAR CODE LABEL WRISTBANDS

Sentry® Bar Code LabelBand®

- Unique slits are tamper evident
- Media supports text, linear and 2D bar codes, photos and graphics
- Clear adhesive shield provides optimal protection from solvents for accurate bar code scanning
- Large label area accommodates both laser and coated thermal labels
- Soft, smooth, non-abrasive material
- Water-resistant, durable, and strong, yet lightweight and comfortable
- Custom imprinting available upon request

Companion Labels

1" Core Direct Thermal	3" Core Direct Thermal	Laser Sheet
THERMDWB42	THERMDWB7	WBWL5080



Product Number	Size of Wristband	Information Area	Max Label	Wristbands/Box
Adult: 5080-**-PDM†	11-3/4" L	3-1/4" L x 1-1/4" H	3 x 1"	500
Adult/Pediatric: 5090-**-PDM	10-3/4" L	2-1/2" L x 1-1/4" H	2" x 1"	

Substitute ** for numbered color:



† 5080 Style band is not available in Pink and Gray

Conf-ID-ent™ Shield Wristbands

- Laminating shield protects patient information
- Supports text, bar codes and graphics
- Can be used with In-A-Snap® clasps
- Non-transferable, secure patented SecurSnap® closure provides maximum security



Product Number	Size of Wristband	Information Area	Max Label	Wristbands/Case
Adult: WBSHDA-**	11-3/4" L x 1-1/4" H	3" L x 1" H	3" x 1"	500
Adult/Pediatric: WBSHDP-**	10-3/4" L x 1-1/4" H	2" L x 1" H	2" x 1"	

Substitute ** for numbered color:



For Alert Snap options, see pages 59-60.

Short Stay® Bar Code LabelBand®

- Media supports text, linear and 2D bar codes, photos and graphics
- Clear adhesive shield provides optimal protection from solvents for accurate bar code scanning
- Tough Tyvek® material is water-resistant, durable and strong, yet lightweight and comfortable
- Soft, smooth, non-abrasive material
- Non-transferable, tamper-evident adhesive closure
- Custom imprinting available upon request



Product Number	Size of Wristband	Information Area	Max Label	Wristbands/Box
3080-**-PDR	10" L x 1" W	3" L x 1" H	2 7/8" x 7/8"	1,000

Substitute ** for numbered color:



Tyvek® is a DuPont registered trademark

BAR CODE ID SOLUTIONS

DIRECT THERMAL LABEL PATIENT ID WRISTBANDS

Sentry® DataMate® Wristbands

- Special ClearImage® surface for superior label adhesion
- Durable, strong, lightweight, comfortable and water-resistant
- Patented SecurSnap® closure provides maximum security
- Media supports text, linear and 2D bar codes, photos and graphics
- Direct thermal labels included
- Custom imprinting available upon request

For Alert Snap options,
see pages 59–60.

Sentry® DataMate® System Adult Wristband




Product Number	Size of Wristband	Information Area	Wristbands/Box
5140-**-PDM	11-1/2" L	3-9/16" L x 1" H	500 (Includes 500 Thermal Labels)

Substitute ** for numbered color:



Sentry® DataMate® System Adult/Pediatric Wristband



Product Number	Size of Wristband	Information Area	Wristbands/Box
5120-11-PDM	10-1/4" L	2-1/4" L x 1" H	500 (Includes 500 Thermal Labels)

Available in White ONLY

Sentry® DataMate® System Infant Wristband



Product Number	Size of Wristband	Information Area	Wristbands/Box
5110-11-PDJ	7" L	1-1/4" L x 11/16" H	250 (Includes 500 Thermal Labels)

Available in White ONLY

DataMate® Thermal Labels

- Accommodates both bar codes and patient data
- White, heat-sensitive thermal surface resists water and alcohol
- Tamper-evident slits



Product Number	Information Area	Labels/Box
Adult: 9130-00-PDR	3 1/4" L x 7/8" H	1,000- 2RLS/500EA
Adult/Pediatric: 9132-00-PDR	2 1/4" L x 7/8" H	1,000- 2RLS/500EA

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.

ALL WRISTBANDS ARE LATEX-FREE.

BAR CODE ID SOLUTIONS

BAR CODE SUPER SOFT DIRECT THERMAL WRISTBANDS

Maximum Comfort Meets Maximum Performance

For Alert Snap options,
see pages 59–60.

Our new ComfortBand™ is specially designed for patients with delicate skin. It provides superior comfort with softer material and a contoured shape. ComfortBand™ is available in multiple sizes to meet the needs of the entire patient population, from geriatrics to infants.

ComfortBand™ Features:

- Softer material is comfortable for even the most delicate skin
- Antimicrobial additive protects wristband surface against tested non-pathogenic bacteria*
- Tamperproof double-post snap closure for security (Adult includes double-post snaps; Infants include single-post heart-shaped snaps.)
- Contour-shaped band comfortably conforms to the wrist
- Large information area

ComfortBand™ Performance:

- Patented PermaPrint® surface provides best imprint quality, maximum resistance to solvents, and seamless printing
- Supports text, linear and 2D bar codes, photo ID and graphics
- Ensures crisp, clear bar codes for easy scannability
- Compatible with most direct thermal printers

ComfortBand™ Adult Thermal Wristbands



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Box
7800-**-PDH	10-13/16" L x 1-1/8" W	4-1/2" x 1-1/8"	1-1/2"	Out	3.7"	100/RL-2RLS/BX
7801-**-PDL			1"		4.0"	200/RL-2RLS/BX
7805-**-PDH			1-1/2"		3.7"	100/RL-2RLS/BX

Substitute ** for numbered color:



ComfortBand™ Infant Thermal Wristbands



Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Box
7850-**-PDH	7-1/4" L x 3/4" W	2-1/4" x 3/4"	1-1/2"	Out	3.7"	100/RL-2RLS/BX
7851-**-PDL			1"		4.0"	200/RL-2RLS/BX
7855-**-PDH			1-1/2"		3.7"	100/RL-2RLS/BX

Substitute ** for numbered color:



All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.



*PDC Healthcare ComfortBand™ Wristbands eliminate more than 99% of non-pathogenic bacteria within 24 hours, based on independent laboratory testing.

ALL WRISTBANDS ARE LATEX & PHTHALATE-FREE.

Exclusively approved and specially designed for use with the Epson SecurColor™ printer

Inkjet Wristband Features:

- Preferred wristband media for the Epson SecurColor™ inkjet color printers
- Extensively tested and verified for compatibility and consistent performance
- Available in snap and adhesive closure options
- Print patient identification information with bar codes, alerts, photo, and graphics directly on the wristband using the Epson SecurColor™ color printer
- Extra smooth, soft material is ideal for patients with sensitive skin
- Latex-free, Phthalate-free, and PVC-free for patient safety
- Wristbands are 100% recyclable to minimize environmental impact
- Print more than 1,500 wristbands using only a single ink cartridge



Inkjet Wristband with Adhesive Closure



- Adhesive closure allows you to print entire length of the band.
- Can be used with foil tamper-evident label which leaves behind the words "VOID" when removed. (Army Branch shown at right. Also available with Navy, Air Force, Marines, Coast Guard Insignias and Blank.)



Tamper-evident labels with branch insignia.

Shows "VOID" once label is removed.

See below for tamper-evident labels.

Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Roll
8020-11-PDL	11-1/2" x 1-1/8"	Entire Band	2"	Out	3.7"	100/RL-4RLS

Inkjet Wristband with Snap Closure



- Patented "X" hole pattern allows you to print entire length of the band while also accommodating a tamper-proof clasp closure. Single-post snaps included.
- Can be used with Ident-Alert® color coded alert snaps for special precautions (Sold separately — see page 59). Ideal when new alerts are identified to prevent re-banding patient.



See Page 59 for Ident-Alert® Snaps.

Product Number	Size of Wristband	Information Area	Core Size	Wound	Outer Diameter	Wristbands/Roll
8000-11-PDL	11-1/2" x 1-1/8"	Entire Band	2"	Out	3.7"	100/RL-4RLS

Foil Tamper-evident Labels



Product Number	Description	Color	Qty./Roll
803-50-PDG	Blank Generic Void Labels	Silver	100/RL
802A-50-PDG	Army Void Labels		100/RL
802AF-50-PDG	Air Force Void Labels		100/RL
802C-50-PDG	Coast Guard Void Labels		100/RL
802M-50-PDG	Marines Void Labels		100/RL
802N-50-PDG	Navy Void Labels		100/RL
802-50-PDM	VA Void Label Set (Set includes 100 each of Army, Air Force, Coast Guard, Marines & Navy)		100/RL-5RLS/BX

BAR CODE ID SOLUTIONS

BAR CODE COLOR LASER WRISTBANDS

ClarID™ Color Laser Wristbands – Specially designed for VA Hospitals to print color laser patient ID wristbands on demand.



PDC Healthcare introduces ClarID™ color laser patient ID wristbands which print on-demand using a color laser printer to print patient photos, color code alerts, bar codes, and text directly onto the wristband. This maximizes patient safety by providing bar code automated identification and supporting color code standardization initiatives.

If you're already using a comparable laser printable ID product, you can seamlessly switch to ClarID™ without IT reformatting.



ClarID™ Wristband with Adhesive Closure

- Scans first time, every time
- Jam-Free™ printing!
- Quick and easy application streamlines the admission process
- Tamper-evident adhesive closure prevents transfer from one patient to another
- Optional single post snap closure can be ordered separately to provide an extra measure of security
- Extra smooth, soft material and contoured edges for patient comfort. Latex-free & Phthalate-free
- Tested and verified by the VA hospital in Washington, DC, ClarID™ lasts up to a month, withstanding moisture, alcohol, and hand sanitizer to preserve patient data and bar codes
- Use Ident-Alert® color coded alert snaps for special precautions (sold separately). Ideal when new alerts are identified to prevent re-banding a patient
- Can be used with foil tamper-evident branch insignia label which leaves behind the words "VOID" when removed. (Army Branch shown below. Also available with Navy, Air Force, Marines and Coast Guard Insignias.)



See page 17 for tamper-evident labels.

Product Number	Size of Wristband	Information Area	Sheets/Box	Sheet Size	Wristbands/Sheet
9350-11-PDL	10-1/4" L x 1" W	Entire Band	400	3.25" W x 11" L	1

For blank and custom imprint options, call (800) 435-4242.

Precision® CompuBand® Single-Post Snap



Product Number	Snaps/Bag
8615-00-PDJ	250

Optional single-post snap closure can be ordered separately to provide an extra measure of security.

Ident-Alert® Color Coded Snaps

ALLERGY	FALL RISK	DNR	NO LATEX	LIMB ALERT
8621A-16-PDH	8621F-14-PDH	8621D-33-PDH	8621LT-22-PDH	8621L-37-PDH

For more Ident-Alert® Color Coded Snaps, see page 59.

Snaps/Box
200

Introducing DuraSoft®

Inspired Softness.

Discover DuraSoft®, a softer laser ID solution inspired by customer demand for a lightweight laser wristband that's more comfortable for patients and soft enough for even the most delicate skin.

Practical Convenience.

Without a stiff, cumbersome laminate overlay, applying DuraSoft® to patients is fast and easy, allowing caregivers more time for patient care. DuraSoft® streamlines steps in the admissions process to reduce wait time and increase workflow efficiency. And because it's paper-based, DuraSoft® performs like a saint in virtually any laser printer.



DuraSoft® Laser Patient ID System: Easy as 1-2-3

1. Print

Easy drop-in replacement works with your existing laser printer and format so there's little to no change needed. Jam Free™ printing in virtually any laser printer.



2. Remove

- + No laminate shields to fold and stick
- + No assembly needed

= Quick & Easy



3. Apply

Fast and easy application saves nursing and admission staff time to improve workflow efficiency and allow more time for patient care.



Meet Stage 2 Meaningful Use Requirements with eMAR*

Using DuraSoft® as part of your electronic medication administration record (eMAR) system can help you meet meaningful use requirements established by the Centers of Medicare & Medicaid Services (CMS) for Electronic Tracking of Medications.



BAR CODE ID SOLUTIONS

*To meet CMS Stage 2 Meaningful Use requirements, hospitals must automatically track more than 10 percent of medication orders using assistive technologies, such as bar code or RFID, with an electronic medication administration record (eMAR). Failure to meet regulations can result in significant incentive payment penalties in fiscal year 2015.



PDC Healthcare DuraSoft® Wristbands eliminate more than 99% of non-pathogenic bacteria within 24 hours, based on independent laboratory testing.

BAR CODE ID SOLUTIONS

BAR CODE LASER PRINTER WRISTBANDS & LABELS


DuraSoft® Features.

- Ultra soft and smooth wristband material provides patient comfort and durability for up to 15 days
- Quick and easy application saves nursing and admission staff time to improve workflow efficiencies
- Easy drop-in replacement works with your existing laser printer system and format so there's little to no change required.
- Secure, tamper-evident adhesive closure
- Patented design featuring unique cross-cuts provides additional options:
 - Optional color-coded alert snaps to identify special-risk patients (See page 59)
 - Optional plastic snap closure for additional security (See page 60)
- Antimicrobial additive protects wristband surface against tested non-pathogenic bacteria*
- Superior Bar Code Scanning. Excellent scan rates for linear and 2D bar codes
- Moisture-resistant. Specially designed material resists moisture, alcohol, and hand sanitizer to preserve patient data and bar codes. Latex-free and Phthalate-free.
- Jam Free™ printing with laser printers
- Standard file holes make it easy to add the laser sheet to a chart or binder



DuraSoft® Wristband and Labels


- Drop in replacement for LaserBand® (PLS-102)



9500-11-PDR

DuraSoft® Mother/Father/Baby


- Drop in replacement for 9329-11-PDM



9510-11-PDM

DuraSoft® Infant/Pediatric

- Drop in replacement for 9333-11-PDM



9511-11-PDM

Product Number	9500-11-PDR	9510-11-PDM	9511-11-PDM
Size of Sheet	8-1/2" W x 11" L	8-1/2" W x 11" L	8-1/2" W x 11" L
Sheets/Box	1,000	500	500
Labels & Wristbands/Sheet	20 Labels & 1 Wristband	1 Crib Card & 4 Wristbands (2 Adult, 2 Infant)	16 Labels & 2 Wristbands
Label & Wristband Sizes	Wristband: 10" L x 3/4" H Labels: 2-1/2" L x 1" H	Adult Wristband: 10" L x 1" H Pediatric Wristband: 7" L x 1" H Label: 3-1/4" L x 2" H	Adult Wristband: 10" L x 1" H Pediatric Wristband: 7" L x 1" H Label: 2-1/2" L x 1" H




*PDC Healthcare DuraSoft® Wristbands eliminate more than 99% of non-pathogenic bacteria within 24 hours, based on independent laboratory testing.

BAR CODE ID SOLUTIONS

BAR CODE LASER PRINTER WRISTBANDS & LABELS

DuraSoft® Wristband and Labels


- Drop in replacement for Standard Register (A7038)



9503-11-PDR


Product Number	9503-11-PDR
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Labels & Wristbands/Sheet	24 Labels & 1 Wristband
Label & Wristband Sizes	Wristband: 10" L x 7/8" H Labels: 2-3/5" L x 1" H

- Drop in replacement for Standard Register (A7015)



9504-11-PDR

Product Number	9504-11-PDR
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Labels & Wristbands/Sheet	32 Labels & 1 Wristband
Label & Wristband Sizes	Wristband: 10" L x 3/4" H Labels: 2-5/7" L x 3/4" H



9508-11-PDR

Product Number	9508-11-PDR
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Labels & Wristbands/Sheet	20 Labels & 1 Wristband
Label & Wristband Sizes	Wristband: 10" L x 3/4" H Labels: 2-1/2" L x 1" H

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.


ALL WRISTBANDS ARE LATEX-FREE.

BAR CODE ID SOLUTIONS

BAR CODE LASER PRINTER WRISTBANDS & LABELS


DuraSoft® Wristband and Labels

- Drop in replacement for LaserBand® (LB2-ADULT-L2)




9505-11-PDR

- Drop in replacement for LaserBand® (PLS-102A)




9509-11-PDR

- Drop in replacement for LaserBand® (LB2-ADULT-L3)



9506-11-PDR

- Drop in replacement for LaserBand® (PLS-101A)



9507-11-PDR



*PDC Healthcare DuraSoft® Wristbands eliminate more than 99% of non-pathogenic bacteria within 24 hours, based on independent laboratory testing.

Product Number	9505-11-PDR	9506-11-PDR	9507-11-PDR	9509-11-PDR
Size of Sheet	8-1/2" W x 11" L	8-1/2" W x 11" L	3" W x 11" L	8-1/2" W x 11" L
Sheets/Box	1,000	1,000	1,000	1,000
Labels & Wristbands/Sheet	20 Labels & 1 Wristband	20 Labels & 1 Wristband	1 Wristband	20 Labels & 1 Wristband
Label & Wristband Sizes	Wristband: 10" L x 7/8" H Labels: 2-1/2" L x 15/16" H	Wristband: 10" L x 7/8" H Labels: 2-1/2" L x 31/32" H	Wristband: 10" L x 3/4" H	Wristband: 10" L x 3/4" H Labels: 2-1/2" L x 1" H

ALL PRODUCTS ON THIS PAGE ARE LATEX-FREE.


DataMate® Laser Printer Wristbands and Labels

The DataMate® system consists of a laser-printable wristband with adhesive closure and labels that adhere to patient records. The wristband's clear, fold-over protective window is designed to prevent printed information from smearing and bleeding to assure error-free reading/scanning. Single-piece construction allows you to print wristbands and labels from the same laser tray.

DataMate® Laser Wristband and Labels

- Wristband's clear, fold-over protective window is designed to protect information from smearing and bleeding
- Accepts linear and 2D bar codes
- One wristband designed for adults with additional labels for duplicating patient information
- Available with tamper-Kimevident slits to promote patient safety

For a complete selection of wristbands & labels, visit us online at www.pdchealthcare.com.




Product Number	9200-**-PDR	TM-PLS-102
Size of Sheet	8-1/2" W x 11" L	8-1/2" W x 11" L
Sheets/Box	1,000	1,000
Labels & Wristbands/Sheet	20 Labels & 1 Wristband (with TE Slits)	20 Labels & 1 Wristband
Label & Wristband Sizes	Wristband: 8-5/8" L x 3/4" H Labels: 2-1/2" L x 1" H	Wristband: 8-5/8" L x 3/4" H Labels: 2-1/2" L x 1" H

Substitute ** for numbered color: 00 White 15 Green

Available with colored borders – 25/cs. minimum (1,000 shts./cs.)

DataMate® Laser Wristband and Labels with File Holes

- Wristband's clear, fold-over protective window is designed to protect information from smearing and bleeding
- Accepts linear and 2D bar codes
- One wristband designed for adults with additional labels for duplicating patient information



Product Number	TM-PLS-WH
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Labels & Wristbands /Sheet	20 Labels & 1 Wristband
Label & Wristband Sizes	Wristband: 8-5/8" L x 3/4" H Labels: 2-1/2" L x 1" H

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.


ALL PRODUCTS ON THIS PAGE ARE LATEX-FREE.

BAR CODE ID SOLUTIONS

BAR CODE LASER PRINTER WRISTBANDS & LABELS

DataMate® Laser Wristband and Labels


- Wristband's clear, fold-over protective window is designed to protect information from smearing and bleeding
- Accepts linear and 2D bar codes
- One wristband designed for adults with additional labels for duplicating patient information



Product Number	TM2-FF-102
Size of Sheet	8-1/2" L x 11-5/8" W
Sheets/Box	1,000
Labels & Wristbands /Sheet	20 Labels & 1 Wristband
Label & Wristband Sizes	Wristband: 8-3/8" L x 1" H Label: 2-11/16" L x 3/4" H

DataMate® Laser Wristband and Labels

- Wristband's clear, fold-over protective window is designed to protect information from smearing and bleeding
- Accepts linear and 2D bar codes
- One wristband designed for adults with additional labels for duplicating patient information



Product Number	TM-PLS-102A
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Labels & Wristbands /Sheet	20 Labels & 1 Wristband
Label & Wristband Sizes	Wristband: 8-5/8" L x 3/4" H Label: 2-1/2" L x 1" H


ALL PRODUCTS ON THIS PAGE ARE LATEX-FREE.

BAR CODE ID SOLUTIONS

BAR CODE LASER PRINTER WRISTBANDS & LABELS

DataMate® Laser Wristband and Labels

- Wristband's clear fold-over protective window is designed to protect information from smearing and bleeding
- One wristband designed for adults with additional labels for duplicating patient information



Product Number	TM-SRC-02
Size of Sheet	8 1/2" W x 11-5/8" L
Sheets/Box	1,000
Labels & Wristband /Sheet	28 labels & 1 wristband
Wristband & Label Sizes	Wristband: 8-5/8" L x 3/4" H Label: 2-1/16" L x 3/4" H

Individual Wristbands

- Economical solution to patient identification
- A single band for easy replacement



Product Number	TM-IND-02A	TM-IND-02-R	TM-IND-02-B
Size of Sheet	3" L x 11"W	3" L x 11" W	3" L x 11" W
Bands/Box	500	500	500
Wristband Size	8-3/8" L x 3/4" W	8-3/8" L x 3/4" W	8-3/8" L x 3/4" W

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.


ALL PRODUCTS ON THIS PAGE ARE LATEX-FREE.

BAR CODE ID SOLUTIONS

BAR CODE LASER PRINTER WRISTBANDS & LABELS

Mother/Father/Baby Wristbands and Crib Card

- Wristband's clear fold-over protective window is designed to protect information from smearing and bleeding
- Accepts linear and 2D bar codes
- Tamper-evident slits on wristband promote patient safety
- Includes ID/Basinette Card




For a complete selection of wristbands & labels, visit us online at www.pdchealthcare.com.

Product Number	TM-MFB-02A
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	500
Wristbands & ID Cards/Sheet	2 Adult Wristbands 2 Baby Wristbands & 1 Crib Card
Wristband & ID Card Sizes	Adult Wristband: 8-3/8" L x 3/4" H Baby Wristband: 4-1/4" L x 1/2" H ID Card: 3-1/4" L x 2" H


Infant Wristbands and Labels

- Accepts linear and 2D bar codes
- 2 infant wristbands
- Also available in individual sheet format



Product Number	TM-INF-02
Size of Sheet	8-1/2" L x 11" W
Sheets/Box	1,000
Labels & Wristbands/Sheet	20 Labels & 2 Wristbands
Wristband & Label Sizes	Wristband: 4-1/4" L x 1/2" H Label: 2-1/2" L x 1" H

TM-INF-02




TM-INDINF-2A

Product Number	TM-INDINF-2A
Size of Sheet	3" L x 11" W
Bands/Box	600
Wristbands/Sheet	2 Wristbands
Wristband Size	4-1/4" L x 1/2" W

ALL PRODUCTS ON THIS PAGE ARE LATEX-FREE.

DataMate® Pediatric Laser Labels

- Wristband's clear, fold-over protective window is designed to protect information from smearing and bleeding
- Accepts linear and 2D bar codes
- One wristband designed for pediatric patients with additional labels for duplicating patient information
- Standard healthcare file hole locations to fit patient charts




Product Number	TM-PEDS-CLR
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Label & Wristband /Sheet	20 labels & 1 wristband
Wristband & Label Sizes	Wristband: 6" L x 1/2" H Label: 2-1/2" L x 1" H

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.

DataMate® Laser Labels with File Holes

- Accepts linear and 2D bar codes
- Standard Healthcare file hole locations to fit patient charts



TM-ADMIT-02
9101X-00-PDR

TM-ADMIT-10

Product Number	TM-ADMIT-02	TM-ADMIT-10
Size of Sheet	8-1/2" W x 11" L	3-3/4" W x 11" L
Sheets/Box	1,000	1,000
Labels/Sheet	30 labels	10 labels
Labels Sizes	2-1/2" L x 1" H	2-1/2" L x 1" H

DataMate® Companion Laser Chart Labels guarantee consistent Jam Free™ Printing of patient information.

- Support text, bar codes and graphics
- Non-adhesive area around the edge prevents printer jams
- Increase work flow and prevent printer maintenance costs

Product Number	9101X-00-PDR
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Labels/Sheet	30 labels
Labels Sizes	2-1/2" L x 1" H


ALL PRODUCTS ON THIS PAGE ARE LATEX-FREE.

BAR CODE ID SOLUTIONS


BAR CODE LASER PRINTER LABELS

DataMate® Laser Labels


Media supports text, linear and 2D bar codes, photos and graphics



Product Number	TM-ADMIT-42
Size of Sheet	8-1/2"W x 11" L
Sheets/Box	1,000
Labels/Sheet	42 labels
Labels Sizes	2-5/8" L x 3/4" H



Product Number	9105-00-PDR
Size of Sheet	8-1/2"W x 11" L
Sheets/Box	1,000
Labels/Sheet	30 labels
Labels Sizes	2" L x 1" H




Product Number	TM-ADMIT-EID
Size of Sheet	8-1/2"W x 11" L
Sheets/Box	1,000
Labels/Sheet	33 labels
Labels Sizes	(27) 2-1/2" L x 1" H (3) 2-1/2" L x 7/16" H (2) 3" L x 3/4" H (1) 1-1/2" L x 3/4" H

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.

ALL PRODUCTS ON THIS PAGE ARE LATEX-FREE.

DataMate® Laser Labels


- Media supports text, linear and 2D bar codes, photos and graphics



For a complete selection of wristbands & labels, visit us online at www.pdchealthcare.com.

Product Number	TM-ADMIT-32
Size of Sheet	8-1/2" L x 11" W
Sheets/Box	1,000
Labels/Sheet	32 labels
Labels Sizes	2-19/32" L x 1" H


Laser Labels



DPSSL1-258
DPSSL1-258-MI

DPSSL1-256

Product Number	DPSSL1-258	DPSSL1-258-MI	DPSSL1-256
Size of Sheet	8-1/2" W x 11" L	8-1/2" W x 11" L	8-1/2" W x 11" L
Sheets/Box	100	1000	250
Labels/Sheet	30 labels	30 labels	33 labels
Labels Sizes	2-5/8" L x 1" H	2-5/8" L x 1" H	2-5/6" L x 1" H



DataMate® Companion Laser Chart Labels guarantee consistent Jam Free™ printing of patient information.

- Media supports text, linear and 2D bar codes, photos and graphics
- Non-adhesive area around the edge prevents printer jams
- Increase workflow and prevent printer maintenance costs

Product Number	DPX-258-MI
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1000
Labels/Sheet	30 labels
Labels Sizes	2-5/8" L x 1" H


ALL PRODUCTS ON THIS PAGE ARE LATEX-FREE.

BAR CODE ID SOLUTIONS

BAR CODE LASER PRINTER LABELS


Laser Labels

- Media supports text, linear and 2D bar codes, photos and graphics
- Laser labels



Product Number	WBWL5080
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Labels/Sheet	24 labels
Labels Sizes	3 L x 1" H


- Media supports text, linear and 2D bar codes, photos and graphics
- Multi-size laser labels



Product Number	WBLBL8
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Labels/Sheet	31 Labels
Labels Sizes	(27) 2-5/8" L x 1" H, (2) Infant 1-7/16" L x 5/8" H, (1) Pediatric 2" L x 1" H, (1) Adult 3" L x 1" H

DataMate® Poly Laser Labels


- Poly labels provide additional strength and durability
- Media supports text, linear and 2D bar codes, photos and graphics
- Tamper-evident slits on 4 wristband labels promote patient safety



Product Number	9145-00-PDR
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Labels/Sheet	32 labels
Labels Sizes	2-5/8" L x 1" H

All Bar Code Wristbands comply with The Joint Commission National Patient Safety Goals, Centers for Medicare & Medicaid Services (CMS), Health Insurance Portability and Accountability (HIPAA), and American Hospital Association (AHA) industry standards.

- Poly labels provide additional strength and durability
- Media supports text, linear and 2D bar codes, photos and graphics
- Tamper-evident slits on 3 wristband labels promote patient safety



Product Number	9140-00-PDR
Size of Sheet	8-1/2" W x 11" L
Sheets/Box	1,000
Labels/Sheet	42 labels
Labels Sizes	2-5/8" L x 3/4" H

ALL PRODUCTS ON THIS PAGE ARE LATEX-FREE.

Insert Style Bar Code Wristbands

- Snap or adhesive closure
- Media supports text, linear and 2D bar codes, photos and graphics
- Insert cards are included
- Custom imprinting available upon request

For a complete selection of wristbands, visit us online at www.pdchealthcare.com.

Safeguard® Sealident® Super Wide Wristband

- Tri-laminate plastic material
- Adhesive closure
- Water-resistant Sealident® pocket
- Will accommodate a 1-1/8" wide bar code label

Companion Labels

1" Core Direct Thermal	3" Core Direct Thermal	Laser Sheet
THERMDWB42	THERMDWB7	WBWL5080



Product Number	Size of Wristband	Information Area	Wristbands/Box
881-**-PDJ	13" L	3-1/2" L x 1-1/8" H	250

Substitute ** for numbered color: 11 13 14 15 16 17
White Blue Yellow Green Red Orange

Securline® Bar Code Wristband

- Strong soft vinyl
- Patented SecurSnap® closure provides maximum security
- Will accommodate a 7/8" wide bar code label



Product Number	Size of Wristband	Information Area	Wristbands/Box
940-**-PDJ	12" L	3-1/2" L x 7/8" H	250

Substitute ** for numbered color: 10 11 15 16
Clear White Green Red

PrimeBand®

- Strong, soft vinyl
- Patented SecurSnap® closure provides maximum security
- Will accommodate a 1-1/4" wide bar code label



Product Number	Size of Wristband	Information Area	Wristbands/Box
860-**-PDJ	12" L	4" L x 1-1/4" H	250

Substitute ** for numbered color: 10 11 15 16 17
Clear White Green Red Orange

ALL WRISTBANDS ARE LATEX-FREE.

BAR CODE ID SOLUTIONS