

Central Service Labels

PDC Healthcare provides an extensive selection of Central Service Labels that communicate important messages among caregivers and staff. Choose from thousands of stock, pre-printed Central Service Labels that communicate pertinent information to identify sterile items, expiration dates, and much more.

PDC Healthcare offers a wide variety of patient charge labels for use with laser, EDP, and thermal printers. Our Keystone Labeling System, which is a stand-alone label printing solution, is the perfect complement to our thermal labels.

Central Service Labels help healthcare organizations meet applicable patient safety regulatory guidelines, including:

- **The Joint Commission National Patient Safety Goals:** Identify Patients Correctly, Improve Staff Communication, Use Medicines Safely, Prevent Infection, Identify Patient Safety Risks, and Prevent Mistakes in Surgery.





Table of Contents

| | |
|--|----|
| Patient Charge Labeling Solutions..... | 88 |
| Keystone™/Healthcoder Labeling Solutions..... | 89 |
| Healthcoder Labeling Solutions/ADS/ESI Nova Systems..... | 91 |
| Patient Charge Cards & Labels..... | 92 |
| Hand-held Labeling Solutions | 93 |
| PromoTouch™ Gun/ Instrument Marking Tape..... | 95 |
| Lot Control Labels/Sterility Labels/ Clean Labels..... | 96 |
| Tray & Return Labels/Spree-D-Remove™ Adhesive Remover..... | 97 |
| Expiration/Implant Labels..... | 98 |

CENTRAL SERVICE

PATIENT CHARGE LABELING SOLUTIONS


Patient Charge Piggyback Laser Labels



DPSSL-PC2-2

DPSSL-PC2-6

DPSSL-PC2-7



DPSSL-PC3-7


DPSSL-PC5-2

DPSSL-PC5-6


DPSSL-PC5-7

| Product Number | Label Size | Sheet Size | Labels/Sheet | Material/Adhesive | Box Qty. |
|----------------|-----------------|--------------|--------------|---|----------------|
| DPSSL-PC2-2 | 2" x 7/8" | 8-1/2" x 11" | 40 | Laser Piggyback/ Glove-Friendly Permanent | 250 SHTS/BX |
| DPSSL-PC2-6 | 2" x 7/8" | | 40 | | |
| DPSSL-PC2-7 | 2" x 7/8" | | 40 | | |
| DPSSL-PC3-7 | 2-5/8" x 13/16" | | 36 | | |
| DPSSL-PC5-2 | 1-3/4" x 3/4" | | 44 | | |
| DPSSL-PC5-6 | 1-3/4" x 3/4" | | 44 | | |
| DPSSL-PC5-7 | 1-3/4" x 3/4" | | 44 | | |


For Monarch 9400 Thermal Printers




DPS-MPC-2




DPS-MPC-6



DPS-MPC-7



DPS-MPC2-2



DPS-MPC2-6

NOTE: There is a black sensor bar on the back of these liners.

| Product Number | Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|----------------|---------------|---------------|-----------|--|----------|
| DPS-MPC-2 | 1-3/8" x 3/4" | 1-5/8" | 4" | Direct Thermal Piggyback/ Extra Permanent | 2,000/RL |
| DPS-MPC-6 | 1-3/8" x 3/4" | 1-5/8" | | | |
| DPS-MPC-7 | 1-3/8" x 3/4" | 1-5/8" | | | |
| DPS-MPC2-2 | 2" x 1" | 2-1/8" | | | |
| DPS-MPC2-6 | 2" x 1" | 2-1/8" | | | |

Keystone™ Labeling System

Stand-alone, On-demand Label Printing

- Facilitates label printing for patient charging and product identification
- Not dependent on IS system - labels can still be printed during an IT system failure and do not require a dedicated computer

Simple Plug & Play

- Pre-loaded with label templates and labels

Small Footprint

- Ultra compact keyboard and heavy-duty industrial printer

Large Memory Capacity

- Supports up to 10,000 unique label records


Flexibility

- Supports multiple bar code formats, including 1-Dimensional (Code 39, Code 128) and 2-Dimensional (Datamatrix/PDF 417)



| Product Number | Description |
|----------------|--------------------------|
| KEYSTN-1 | Keystone Labeling System |

For Keystone™ & HealthCoder® Plus Printers



HCP-TD-115-Y



HCP-TD-115-O




HCP-TD-115-W




HCP-TD-115-B

Keystone™
Compatible!




HCP-TD-115-G



HCP-TD-115-V

| Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|-------------|---------------|-----------|---|------------------|
| 1.575" x 1" | 1-21/32" | 3" | Direct Thermal Piggyback/ Glove-Friendly Permanent | 3,500/RL-6RLS/BX |




HCP-TD-W
HCP-TD-W2

| Product Number | Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|----------------|-------------|---------------|-----------|---|------------------|
| HCP-TD-W | 1-1/2" x 1" | 1-5/8" | 3" | Direct Thermal/Glove-Friendly Permanent | 2,250/RL-6RLS/BX |
| HCP-TD-W2 | | | | Direct Thermal/Glove-Friendly Extra Permanent | |


For Keystone™ and all HealthCoder® Printers



DPTH-212112



DPTH-21SY



DPTH-OR21


Keystone™
Compatible!

| Product Number | Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|----------------|-----------------|---------------|-----------|--|------------------|
| DPTH-212112 | 1-1/2" x 2-1/2" | 1-5/8" | 1-5/8" | Direct Thermal/ Glove-Friendly Permanent | 800/RL-6RLS/BX |
| DPTH-21SY | 1" x 2" | 1-1/8" | | | 1,000/RL-6RLS/BX |
| DPTH-OR21 | 1" x 2" | 1-1/8" | | | |

CENTRAL SERVICE

KEYSTONE™ / HEALTHCODER® LABELING SOLUTIONS

For Keystone™ and all HealthCoder® Printers



DPTH-378178



DPTH-31




DPTH-OR31


Keystone™
Compatible!

| Product Number | Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|----------------|-----------------|---------------|-----------|--|----------------|
| DPTH-378178 | 1-7/8" x 3-7/8" | 2" | 1-5/8" | Direct Thermal/ Glove-Friendly Permanent | 500/RL-6RLS/BX |
| DPTH-31 | 1" x 3" | 1-1/8" | | | 700/RL-6RLS/BX |
| DPTH-OR31 | 1" x 3" | 1-1/8" | | | 700/RL-6RLS/BX |


For all HealthCoder® 5000 Printers



DPSTH5-Y




DPSTH5-B




DPSTH5-G

| Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|----------------|---------------|-----------|--|------------------|
| 1.575" x .945" | 1-21/32" | 1-5/8" | Direct Thermal Piggyback/Glove- Friendly Permanent | 1,500/RL-6RLS/BX |


For all HealthCoder® 8200, 8000 and 5000 Printers




DPSTH-W




DPSTH-Y




DPSTH-O



DPSTH-B



DPSTH-G




DPSTH-V




DPSTH-YNX

| Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|----------------|---------------|-----------|---|------------------|
| 1.575" x .945" | 1-21/32" | 1-5/8" | Direct Thermal Piggyback/ Glove-Friendly Permanent | 3,000/RL-6RLS/BX |

DPSTH-YNX is recommended for hand application.



DPTH-916W
DPTH-916W2




DPTH-916SY
DPTH-916Y2

| Product Number | Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|----------------|----------------|---------------|-----------|--|------------------|
| DPTH-916W | 1.575" x .945" | 1-7/10" | 1-5/8" | Direct Thermal/Glove Friendly Permanent | 2,250/RL-6RLS/BX |
| DPTH-916W2 | | | | Direct Thermal/Glove Friendly Extra Permanent | 2,250/RL-6RLS/BX |
| DPTH-916SY | | | | Direct Thermal/Glove Friendly Permanent | 2,200/RL-6RLS/BX |
| DPTH-916Y2 | | | | Direct Thermal/ Extra Permanent | 2,250/RL-6RLS/BX |


For all HealthCoder® 8x00 Printers

| Product Number | Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|----------------|-----------------|---------------|-----------|---|-------------------|
| HC-1 | 1-6/10" x 9/10" | 1-11/16" | 1-5/8" | Direct Thermal Piggyback/ Glove-Friendly Permanent | 3,000/RL-10RLS/BX |
| HC-3 | | | 3" | | |



HC-1
HC-3


For Industrial Size Thermal Printers



DPS-ALL-2



DPS-HMS-2



DPS-JIT-Y

| Product Number | Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|----------------|------------------|---------------|-----------|---|------------------|
| DPS-ALL-2 | 2" x 1" | 2-1/8" | 3" | Direct Thermal Piggyback/Extra Permanent | 1,500/RL-6RLS/BX |
| DPS-HMS-2 | 1-3/4" x 3/4" | 2" | 3" | Direct Thermal Piggyback/Glove-Friendly Permanent | 3,500/RL |
| DPS-JIT-Y | 1-9/16" x 13/16" | 1-5/8" | 1-5/8" | Direct Thermal Piggyback/Extra Permanent | 3,000/RL-6RLS/BX |

For ADS and Other Patient Charge Systems



DPS-PCBAR2-6

| Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|------------------|---------------|-------------|--|-----------|
| 2-4/10" x 15/16" | 11" | Fanfold/EDP | Piggyback Litho/Glove-Friendly Permanent | 25,000/BX |

For ESI Nova Systems



DPS-WDB4-2




DPS-WDB4-6

| Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|----------------|---------------|-------------|---------------------------|-----------|
| 1-7/10" x 5/6" | 9-1/2" | Fanfold/EDP | Flexo Piggyback/Permanent | 25,000/BX |

CENTRAL SERVICE

PATIENT CHARGE CARDS & LABELS

For Multi-Purpose Patient Charge Labeling





| Label Size | Carrier Width | Core Size | Material/Adhesive | Qty. |
|---------------|---------------|-------------|---------------------------|-----------|
| 1-3/4" x 3/4" | 10-7/8" | Fanfold/EDP | Flexo Piggyback/Permanent | 25,000/BX |

Patient Charge Cards

100# paper cards for all HealthCoder® piggyback labels. For compatible HealthCoder® Labels, see pages 89-91.

| Product Number | Card Size | Card Stock | Qty. |
|----------------|-----------|------------|--------|
| DPSTH-PCC | 7" x 5" | 100# | 250/BX |




DPSTH-PCC

Our patient charge cards are 100# paper cards for piggyback labels. We carry a large selection of piggyback labels. See pages 88-93.

Piggyback labels are used whenever a label must be applied to one surface and re-applied to another surface. The two-ply construction has a bottom label, which adheres to the chargeable item, and a top label that can be removed from the bottom label and transferred to the patient's charge card.

| Product Number | Card Size | Card Stock | Qty. |
|----------------|-----------|------------|--------|
| DPS-57 | 7" x 5" | 100# | 250/BX |



DPS-57



CS-2® Hand-Held Labeling System

The CS-2® Labeling System has been designed primarily for labeling any items in Central Service/Sterile that need to be autoclaved. It is also great for labeling in Dietary, Pharmacy and Laboratory.

With the lightweight CS-2® labeler, you can print two lines of information with up to eight characters per line. The available characters include two Gregorian dates for expiration and autoclave dating with additional bands for load and sterilizer numbers. The CS-2® System single ply labels are available with a specially developed autoclavable adhesive in a variety of colors and styles.




CS-2



CS-2INK

| Product Number | Description | Qty. |
|----------------|---|--------|
| CS-2 | Two-line hand-held labeling device | 1EA/BX |
| CS-2INK | Ink roller for CS-2 hand-held labeling device | |

CS-2® Autoclavable Labels



| Product Number | Label Size | Core Size | Material /Adhesive | Qty. |
|----------------|---------------|-----------|-------------------------------|-------------------|
| CS-W | 15/16" x 5/8" | 1" | Litho/ Permanent | 1,000/RL-24RLS/BX |
| CS-Y | | | | |
| CS-O | | | | |
| CS-BL | | | | |
| CS-G | | | | |
| CS-R | | | | |
| CS-E | | | | |
| CS-RO | | | | |
| CS-VS | | | | |
| CS-YP | | | | |
| CS-YP | | | Litho Piggyback/ Permanent | 700/RL-24RLS/BX |

PC-6® Hand-Held Labeling Device

The PC-6® System provides you with the ability to print three lines of information with up to 12 characters per line – including alpha, numeric and symbols. Great for Patient Charging!



| Label Size | Core Size | Material/Adhesive | Qty. |
|------------------|-----------|---|-----------------|
| 7/8" x 1-22/100" | 1" | Smudgeproof Litho Piggyback/ Permanent | 750/RL-12RLS/BX |



| Product Number | Description | Qty. |
|----------------|--|------|
| PC6-LD | Three-line hand-held labeling device | 1EA |
| PC6-INK | PC6® ink roller can also be used in the PC5-LD labeling device | |

CENTRAL SERVICE


HAND-HELD LABELING SOLUTIONS

TXII™ Hand-Held Labeling System

Ergonomically designed to offer optimum ease of operation, the TXII™ Hand-Held Labeling System allows you to easily date and code items in the Laboratory, Central Service, Warehouse and Pharmacy departments, or wherever a dating or coding system is required.

The TXII™ device's rugged, high-impact plastic provides for exceptional ease of handling and consistent printing on all labels. The TXII™ system provides you with:


- Imprinted labels that display one line of information
- 10 digit dial-in capabilities, including Gregorian dating or large numbers
- Labels in a variety of colors, construction and adhesives



| Product Number | Description | Qty. |
|----------------|---|--------|
| TXII-2 | TXII Labeler with Gregorian Dating in Positions 1 & 8 | 1EA/BX |
| TXII-3 | TXII Labeler with all Large Numbers | 1EA/BX |
| TX-IR5 | Ink roller for TXII hand-held labeling devices | 5/PKG |

Visit us online at
www.pdchealthcare.com
 for more labeling systems.

TXII™ Hand-Held Labeling System Labels



| Product Number | Label Size | Core Size | Material/Adhesive | Qty. |
|----------------|------------|-----------|---------------------------------------|----------|
| TX-B | 1" x 1/2" | 1" | Smudgeproof Litho/Permanent | 1,500/RL |
| TXP-2 | | | Smudgeproof Litho/Permanent | 1,500/RL |
| TX-E | | | Smudgeproof Litho/Permanent | 1,500/RL |
| TX-PY | | | Smudgeproof Litho Piggyback/Permanent | 1,000/RL |
| TX-RO | | | Smudgeproof Litho/Permanent | 1,500/RL |

PromoTouch™ Label Gun Saves Time!

- Dispenses Promotouch™ compatible labels (see below)
- Quickly & easily label your supplies



| Product Number | Qty. |
|----------------|------|
| PROGUN | 1 EA |

PromoTouch™ Compatible Labels



59707980 — Yellow (R)



59705684 — Yellow (P)



59705679 — Fl. Orange (P)

| Product Number | Size | Roll Qty. |
|----------------|-----------------|-----------|
| 59707980 | 1-3/4" X 1-1/8" | 250 |
| 59705684 | | |
| 59705679 | | |


(P) = Permanent Adhesive
(R) = Removable Adhesive

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

Autoclavable Labels and Ribbons

| Synthetic Label with Autoclavable Adhesive | | | | |
|--|-------------|-----------|---|----------------------|
| Product Number | Label Size | Core Size | Material/Adhesive | Qty. |
| TTRS-1124A | 4" x 1-1/2" | 3" | Thermal Transfer Synthetic/ Glove-Friendly Permanent | 3,600/RL- 2RLS/BX |

| Thermal Transfer Ribbon | | | | |
|-------------------------|----------------|---------------|-------|------|
| Product Number | Ribbon Size | Type | Color | Qty. |
| TTR-FRD10233 | 4.02" x 1,476' | Premium Resin | Black | 1/BX |



TTR-FRD10233

TTRS-1124A

Instrument Marking Tape

Instromark® tape instantly color codes all surgical instruments for easy tracking and identification.

- Durable, vinyl tape
- Available in 9 colors, plus transparent
- Autoclavable



IN-600-0
Transparent

IN-600-1

IN-600-2

IN-600-3

IN-600-5

IN-600-6

IN-600-7

IN-600-12

IN-600-18

IN-600-19

| Tape Width | Tape Length | Core Size | Material/Adhesive | Roll Size |
|------------|-------------|-----------|----------------------|-------------|
| 3/8" | 200" | 1" | Vinyl Tape/Removable | 3/8" x 200" |



CENTRAL SERVICE

LOT CONTROL LABELS / STERILITY LABELS / CLEAN LABELS

Lot Control Labels



59704156 — Fl. Red (P)



59702206 — Fl. Orange (P)




59704160 — Green (P)

For more Lot Control Labels, see page 137.


| Label Size | Roll Qty. |
|---------------|-----------|
| 1-5/8" x 7/8" | 1,000 |

Sterility Labels




MV04FR0465 — Fl. Red (P)

| Label Size | Roll Qty. |
|-------------|-----------|
| 2-1/4" x 1" | 420 |




MV04SY4849 — Yellow (P)



59705151 — Red (P)


| Label Size | Roll Qty. |
|---------------|-----------|
| 2-1/4" x 7/8" | 1,000 |

Clean Labels




MV08FG2692 — Fl. Green (R)

| Product Number | Size | Roll Qty. |
|----------------|---------------|-----------|
| MV08FG2692 | 2-15/16" x 2" | 333 |
| 59705743 | 4" x 2-1/8" | 500 |
| 59712529 | 4" x 2-1/8" | 500 |



59705743 — Fl. Green (R)



59712529 — Fl. Yellow (P)



MV04FC0460 — Fl. Chartreuse (R)

| Label Size | Roll Qty. |
|-------------|-----------|
| 2-1/4" x 1" | 420 |


See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

(P) = Permanent Adhesive
(R) = Removable Adhesive


Tray Labels



59705785 — Fl. Yellow (R)




MV07FO4885 — Fl. Orange (R)



59705783 — Fl. Orange (R)
59705784 — Orange (P)



59705377 — Fl. Green (P)



59705587 — Fl. Red (P)

| Product Number | Label Size | Roll Qty. |
|----------------|-------------------|-----------|
| 59705785 | 2" x 5" | 500 |
| MV07FO4885 | 2-15/16" x 1-1/2" | 333 |
| 59705783 | 2-3/8" x 1-3/4" | 1,000 |
| 59705784 | 2-3/8" x 1-3/4" | 1,000 |
| 59705377 | 1-1/4" x 3/8" | 1,000 |
| 59705587 | 1-5/8" x 7/8" | 1,000 |

Return To Labels



MV04FO7481 — Fl. Orange (R)

| Label Size | Roll Qty. |
|-------------|-----------|
| 2-1/4" x 1" | 420 |

(P) = Permanent Adhesive
(R) = Removable Adhesive



Spee-D-Remove™ Adhesive Remover

Quickly and painlessly remove adhesive residue with Spee-D-Remove™ adhesive remover.

- A few drops is all you need — simply apply and wipe with cloth!
- Works on metal, glass, plastic, ceramic, chrome, and most other hard surfaces
- Good for clean up jobs both big and small

| Product Number | Description | Qty. |
|----------------|--|------|
| SPD-RM2 | Spee-D-Remove™ Adhesive Remover — 8 oz. Bottle | 1 EA |

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

CENTRAL SERVICE

EXPIRATION / IMPLANT LABELS


Expiration Labels



59705156 — Fl. Orange (P)



MV01FC1588 — Fl. Chartreuse (R)



5970719181000 — Fl. Red (P)
5970719182000 — Fl. Red (R)



MV04FC1588 — Fl. Chartreuse (R)



59705704 — Pink (R)

| Product Number | Label Size | Roll Qty. |
|----------------|-----------------|-----------|
| 59705156 | 2-3/8" x 1-3/4" | 1,000 |
| MV04FC1588 | 2-1/4" x 1" | 420 |
| MV01FC1588 | 1-1/4" x 5/16" | 760 |
| 5970719181000 | 1-5/8" x 7/8" | 1,000 |
| 5970719182000 | 1-5/8" x 7/8" | 1,000 |
| 59705704 | 2-1/4" x 7/8" | 1,000 |


Implant Labels



MV06FP0468 — Fl. Pink (R)



59705722 — Fl. Red (R)



MV03FP1228P — Fl. Pink (P)



59705726 — Fl. Pink (R)
59705729 — Fl. Pink (P)

| Product Number | Label Size | Roll Qty. |
|----------------|---------------|-----------|
| MV06FP0468 | 2-15/16" x 1" | 333 |
| 59705722 | 2-7/8" x 7/8" | 1,000 |
| MV03FP1228 | 1-7/16" x 1" | 666 |
| 59705726 | 2-1/4" x 7/8" | 1,000 |
| 59705729 | 2-1/4" x 7/8" | 1,000 |

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

(P) = Permanent Adhesive
(R) = Removable Adhesive