


NEW Specialty Cryogenic Labels



Superior Performance in Temperatures as Low as -196°C / -320°F

These specialty labels are specifically engineered for clear identification during cryogenic transportation and storage, withstanding the harshest, low-temperature environments. They are durable, with top performance even when exposed to liquid nitrogen or ultra-low temperature freezers, thaw cycles, chemicals and solvents.



- Qualified through rigorous durability testing in extreme environments to ensure top performance
- Greater adhesion to small rounded surfaces such as test tubes
- Endure extended storage in low temperature exposures such as liquid nitrogen, and low temperature freezers
- Prevents edge-lifting even when repeatedly subjected to drastic temperature changes
- Polyester coated to ensure durability and crisp, clean printing
- Available in a wide variety of sizes and configurations for thermal transfer, direct thermal and laser printers

Cryogenic Label Options

Product Number/SKU	Printer Type	Label Size	Quantity
LASCR1	Laser	2 ⁵ / ₈ " x 1"	30 Per Sheet / 100 Sheets per Box
LASCR2	Laser	1 ¹ / ₂ " x ³ / ₄ "	60 Per Sheet / 100 Sheets per Box
LASCR3	Laser	1 ¹ / ₄ " x ¹ / ₂ "	80 Per Sheet / 100 Sheets per Box
LASCR4	Laser	1" Circle	70 Per Sheet / 100 Sheets per Box
LASCR5	Laser	¹ / ₂ " Circle	90 Per Sheet / 100 Sheets per Box
LASCR6	Laser	³ / ₄ " Circle	90 Per Sheet / 100 Sheets per Box
LASCR7	Laser	1 ¹ / ₂ " x 1"	30 Per Sheet / 100 Sheets per Box
TTCR1	Thermal Transfer	2" x ¹ / ₄ " – 1" Core	2000 Labels per Roll
TTCR2	Thermal Transfer	1 ¹ / ₂ " x ³ / ₄ " – 1" Core	2000 Labels per Roll
TTCR3	Thermal Transfer	1" x 2 ¹ / ₂ " – 1" Core	2000 Labels per Roll
TTCR4	Thermal Transfer	1" x 2 ¹ / ₂ " – 3" Core	3000 Labels per Roll
TTCR5	Thermal Transfer	1" x 2 ¹ / ₄ " – 3" Core	3000 Labels per Roll
TTCR6	Thermal Transfer	2" x ¹ / ₂ " – 3" Core	3000 Labels per Roll
TTCR7	Thermal Transfer	1 ¹ / ₂ " x ³ / ₄ " – 3" Core	3000 Labels per Roll
DTCR1	Direct Thermal	1 ¹ / ₂ " x ¹ / ₂ " – 1" Core	2000 Labels per Roll
DTCR2	Direct Thermal	1 ¹ / ₂ " x ³ / ₄ " – 1" Core	2000 Labels per Roll
DTCR3	Direct Thermal	1" x ¹ / ₂ " – 1" Core	2000 Labels per Roll
DTCR4	Direct Thermal	1" x 2 ¹ / ₄ " – 3" Core	3000 Labels per Roll
DTCR5	Direct Thermal	2" x ¹ / ₂ " – 3" Core	3000 Labels per Roll
DTCR6	Direct Thermal	1 ¹ / ₂ " x ³ / ₄ " – 3" Core	3000 Labels per Roll
DTCR7	Direct Thermal	1" x 2 ¹ / ₂ " – 3" Core	3000 Labels per Roll

Contact us today for samples and more information.
 Call **800-435-4242** | pdchealthcare.com