

Disposable Pillowcases

- Control contamination and eliminate laundering
- Moisture resistant and comfortable
- Choice of fabric, and color size

Product Number	Material Type	Color	Size	Qty.
SP-365	Tissue/Poly	White	21" x 30"	100/Case
SP-363	Tissue/Poly	Blue	21" x 30"	100/Case
SP-355	Fabrice ^l *	White	21" x 30"	100/Case
SP-353	Fabrice ^l *	Blue	21" x 30"	100/Case

*A super-strong, fluid-proof and tear-resistant material with superior draping and softness.



Drape Sheets

Provide comfort and modesty with our drape sheet products.

Product Number	Material Type	Color	Size	Qty.
SP-305	2-Ply Tissue	White	36" x 48"	100/Case
SP-022**	2-Ply Tissue	White	40" x 48"	100/Case
SP-302	2-Ply Tissue	White	40" x 48"	100/Case
SP-310	3-Ply Tissue	White	40" x 48"	100/Case
SP-311	3-Ply Tissue	White	40" x 60"	100/Case
SP-309	3-Ply Tissue	White	40" x 72"	50/Case
SP-374	T/P/T	Blue	40" x 48"	50/Case
SP-375	T/P/T	Blue	40" x 60"	50/Case
SP-370	T/P/T	Blue	40" x 90"	50/Case

T/P/T is Tissue/Paper/Tissue

**Space Saver - 4 pack of 25 drapes per case



Smooth Table Paper - White

Choice of narrow or wide rolls and long or short paper.

Product Number	Style	Size	Qty.
PP-182	Premium Smooth	18" x 225"	12 Rolls/Case
PP-183	Premium Smooth	18" x 260"	12 Rolls/Case
PP-202	Premium Smooth	20" x 225"	12 Rolls/Case
PP-212	Premium Smooth	21" x 225"	12 Rolls/Case
PP-213	Premium Smooth	21" x 260"	12 Rolls/Case
PP-S243	Premium Smooth	24" x 260"	8 Rolls/Case



EXAM ROOM SUPPLIES

SHORTS/GOWNS/CAPES/PONCHOS



Sewn elastic waistband for secure fit.

Generous cut for ample room.

Disposable Exam Shorts

- Non-woven S/M/S (Spunbond/Meltblown/Spunbond)
- Repels moisture — ideal for whirlpools, hydrotherapy and ice packs
- Security and modesty, even when wet
- Latex-free

Product Number	Material Type	Color	Size	Qty.
PP-919	S/M/S Non-woven	Blue	Medium (28"-34" waist)	100/Case
PP-917	S/M/S Non-woven	Blue	Large (32"-40" waist)	100/Case
PP-918	S/M/S Non-woven	Blue	X-Large (38"-48" waist)	100/Case
PP-920	S/M/S Non-woven	Blue	Plus Size (60" waist)	50/Case

Available in 3 Colors




Front/Back Opening T/P/T

PP-420

Exam Gowns


- Knee length coverage for hygiene, comfort and modesty
- Choice of colors
- Available in 3-Ply T/P/T or 2-Ply Tissue/Poly

Colors:



Product Number	Material	Color	Size	Qty.
PP-420	3-Ply T/P/T	White	30" x 42"	50/Case
PP-834	3-Ply T/P/T	Mauve	30" x 42"	
PP-520	3-Ply T/P/T	Blue	30" x 42"	
PP-525	Tissue/Poly	White	30" x 42"	
PP-526	Tissue/Poly	Blue	30" x 42"	

Available in 4 Colors



MAM-516

Exam Capes

- Access to upper torso is easy
- 3-Ply T/P/T (Tissue/Poly/Tissue)
- Choice of colors

Colors:



Product Number	Material	Color	Size	Qty.
PP-C315	3-Ply Tissue	White	30" x 21"	100/Case
PP-916	3-Ply T/P/T	Mauve	30" x 21"	
PP-415	3-Ply T/P/T	White	30" x 21"	
MAM-516	3-Ply T/P/T*	Mauve	30" x 21"	
MAM-515	3-Ply T/P/T	Blue	30" x 21"	
MAM-518	3-Ply T/P/T	Peach	30" x 21"	
MAM-510	3-Ply T/P/T*	Blue	30" x 21"	

*Indicates recycled tissue/poly/tissue.



Snap closure

Oversized Mammography Ponchos

- Accommodates full-figured patients

Product Number	Material	Color	Closure	Size	Qty.
PP-981	Non-woven	Blue	Snap Closure	45" x 48"	50/Case

Our disposables are available in a wide selection of materials!



2-Ply & 3-Ply Tissue
An economical fabric that features fluid absorbency and supple comfort. Recommended for shorter procedures.



3-Ply Tissue/Poly/Tissue
Fluid-proof poly film is laminated between two layers of soft, absorbent tissue. More durable than standard tissue.



Non-woven Fabric
Bonded man-made fibers, offering strength, softness and durability. Ideally suited for lengthy procedures.