

Laboratory Labels

PDC Healthcare knows that when it comes to patient safety in laboratory operations, positive identification is mission critical. We offer labeling systems designed specifically for LIS systems, including Sunquest, Cerner, SoftLab, Meditech, and others.

- **V-Notch™ Secondary Labels**, along with BD Vacutainer® Plus Plastic Tubes, provide correct label placement the first time and every time! So you can ensure accurate identification and an optimal healthcare experience for patients and staff alike.
- **ISBT Compliant Bank Labels** enhance patient safety, while supporting ISBT 128 standards.
- **Point-of-Care Labeling Systems** use bar code scanning and wireless mobility for convenient, positive patient identification.

Laboratory Labels help healthcare organizations meet applicable patient safety regulatory guidelines, including:

- **The Joint Commission National Patient Safety Goals:** Identify Patients Correctly, Improve Staff Communication, Use Medicines Safely, Prevent Infection, Identify Patient Safety Risks, and Prevent Mistakes in Surgery.



Table of Contents

V-Notch™ Labels.....	106
System Labels.....	107
Point-of-Care Labeling Solutions.....	126
ISBT Blood Bank Labels.....	127
Microscope Slide Labeling Solutions.....	128
Consecutive Number Labels/ Tamper-Evident Labels & Tape.....	133
Color Circles & Specialty Tape	134
Time Tape®	135
Stat/Specimen Labels.....	136
Inventory Management/ Lot Control Labels.....	137
Refrigeration/Blood Labels.....	138
Hazard Communication Labels	139
General Purpose Labels/Tape.....	140
Custom Labeling.....	142



View the V-Notch™ Video Online

- Learn how you and your lab can benefit from standardized blood collection tube labeling for accuracy, efficiency and patient safety.
- Discover how V-Notch™ Secondary Labels also help save you time and money!
- Watch simple, step-by-step instructions on how to apply V-Notch™ Secondary Labels to BD Vacutainer® Plus Plastic Tubes.

www.pdchealthcare.com/vnotch

LABORATORY

V-NOTCH™ LABELS

V-Notch™ Secondary Labels:

The Simple Solution for Accurate Labeling & ID. V-Notch™ Secondary Labels provide quick, easy and accurate identification of BD Vacutainer® Plus Plastic Tubes. The V-shaped notch aligns with the color-coded notch and sidebar to provide:

- Consistent & Easy Labeling Guides
- Accurate & Reliable Bar Code Scanning
- Quicker Test Results
- Operational Efficiencies That Save Time & Money
- Increased Clinician Satisfaction
- Patient Safety & Satisfaction
- Accurate ID



before



after

Look for labels with the  symbol.

To Apply the V-Notch™ Labels:


STEP 1: Peel your V-Notch™ Secondary Label from the liner.




STEP 2: Visually align the V-Notch™ with the color-coded notch on the BD Vacutainer® Plus plastic tube.



STEP 3: Before you apply the V-Notch™ Secondary label, make sure the color-coded notch and the color-coded sidebar on the BD Vacutainer® Plus Plastic Tube are clearly visible. Smoothly wrap the label around the tube.



Meditech Labeling Solutions (For Direct Thermal Printers)



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
LDWHP58	4" x 1-1/4"	4-1/8"	1"	In	Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-8RLS/CS
LDWHP61			1"	In	Direct Thermal/Glove-Friendly Permanent	Yellow	1,000/RL-8RLS/CS
LDWHP59			1"	In	Direct Thermal/Glove-Friendly Permanent	Blue Border	1,000/RL-8RLS/CS
LDWHP62			1"	In	Direct Thermal/Glove-Friendly Permanent	Orange Border	1,000/RL-8RLS/CS
LDWHP63			1"	In	Direct Thermal/Glove-Friendly Permanent	Red Border	1,000/RL-8RLS/CS
LDWHP60			1"	In	Direct Thermal/Glove-Friendly Permanent	Green Border	1,000/RL-8RLS/CS
LDWHP64			3"	Out	Direct Thermal/Glove-Friendly Permanent	White	3,600/RL
LMEDY91			3"	Out	Direct Thermal/Glove-Friendly Permanent	Yellow	3,600/RL
LDWHP65			3"	Out	Direct Thermal/Glove-Friendly Permanent	Blue Border	3,600/RL
LDWHP66			3"	Out	Direct Thermal/Glove-Friendly Permanent	Orange Border	3,600/RL
LDWHP67			3"	Out	Direct Thermal/Glove-Friendly Permanent	Red Border	3,600/RL
LMEDG91			3"	Out	Direct Thermal/Glove-Friendly Permanent	Green Border	3,600/RL
LMEDW10			1"	In	Direct Thermal IR/Glove-Friendly Permanent	White	1,000/RL-8RLS/CS
LMEDY10			1"	In	Direct Thermal IR/Glove-Friendly Permanent	Yellow	1,000/RL-8RLS/CS
LMEDB10			1"	In	Direct Thermal IR/Glove-Friendly Permanent	Blue Border	1,000/RL-8RLS/CS
LMEDO10			1"	In	Direct Thermal IR/Glove-Friendly Permanent	Orange Border	1,000/RL-8RLS/CS
LMEDR10			1"	In	Direct Thermal IR/Glove-Friendly Permanent	Red Border	1,000/RL-8RLS/CS
LMEDG10			1"	In	Direct Thermal IR/Glove-Friendly Permanent	Green Border	1,000/RL-8RLS/CS
LMED51			3"	Out	Direct Thermal IR/Glove-Friendly Permanent	White	3,600/RL
LMEDY51			3"	Out	Direct Thermal IR/Glove-Friendly Permanent	Yellow	3,600/RL
LMEDB51			3"	Out	Direct Thermal IR/Glove-Friendly Permanent	Blue Border	3,600/RL
LMEDO91			3"	Out	Direct Thermal IR/Glove-Friendly Permanent	Orange Border	3,600/RL
LMEDR91			3"	Out	Direct Thermal IR/Glove-Friendly Permanent	Red Border	3,600/RL
LMEDG51			3"	Out	Direct Thermal IR/Glove-Friendly Permanent	Green Border	3,600/RL
LDWHP155			1"	Out	Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-8RLS/CS
LDWHP160			1"	Out	Direct Thermal/Glove-Friendly Permanent	Red Border	1,000/RL-8RLS/CS



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TD1-MIT-3BD	4" x 1-1/4"	4-1/8"	1"	Out	Direct Thermal/Glove-Friendly Permanent	White	1,000/RL - 6RLS/BX
TD-MIT-3BD			3"		Direct Thermal/Glove-Friendly Permanent	White	3,600/RL - 1RL/BX
TD-MIT-3BBD			3"		Direct Thermal/Glove-Friendly Permanent	Blue Border	3,600/RL - 1RL/BX
TD-MIT-3GBD			3"		Direct Thermal/Glove-Friendly Permanent	Green Border	3,600/RL - 1RL/BX
TD-MIT-3RBD			3"		Direct Thermal/Glove-Friendly Permanent	Red Border	3,600/RL - 1RL/BX
DPS-MIT-3BD			3"		Direct Thermal IR/Glove-Friendly Permanent	White	3,600/RL - 1RL/BX

BD, BD Logo, and V-Notch™ trademarks are the property of Becton, Dickinson and Company. ©2005 BD.

LABORATORY

MEDITECH SYSTEM LABELS

Meditech Labeling Solutions

LDWHP146 LDWHP125 LDWHP123 LDWHP165 3-3/8" x 2-3/8"				LDWHP126 LMED52		LDWHP154 LDWHP153		LDWHP164		LDWHP124		LDWHP177	
1-1/8" x 2"													
9/16" x 7/8"													

Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
LDWHP146	3-3/8" x 5-1/4"	3-1/2"	3'	Out	Direct Thermal/Glove-Friendly Permanent	White	500/RL-2RLS/CS
LDWHP125			3'		Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-2RLS/CS
LDWHP126			3'		Direct Thermal/Glove-Friendly Permanent	Pink Border	1,000/RL-2RL/CS
LDWHP154			3'		Direct Thermal/Glove-Friendly Permanent	Purple Border	1,000/RL-2RLS/CS
LDWHP164			3'		Direct Thermal/Glove-Friendly Permanent	Red Border	1,000/RL-2RLS/CS
LDWHP123			3'		Direct Thermal IR/Glove-Friendly Permanent	White	1,000/RL-2RLS/CS
LDWHP124			3'		Direct Thermal IR/Glove-Friendly Permanent	Pink Border	1,000/RL-2RLS/CS
LDWHP153			3'		Direct Thermal IR/Glove-Friendly Permanent	Purple Border	1,000/RL-2RLS/CS
LDWHP165			1'		Direct Thermal/Glove-Friendly Permanent	White	500/RL
LMED52			1'		Direct Thermal/Glove-Friendly Permanent	Pink Border	500/RL
LDWHP177			3'		Direct Thermal/Glove-Friendly Permanent	Yellow Border	1,000/RL-2RLS/CS

LDWHP712 3-3/8" x 2-1/4"			
1-1/8" x 2-1/4"			
9/16" x 3/4"			

LDWHP127 LDWHG01			
1-1/8" x 2-1/4"			
9/16" x 3/4"			

LDWHP128 LDWHP71 LMED61 LDWHP711			
1-1/8" x 2-1/4"			
9/16" x 3/4"			

Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
LDWHP712	3-3/8" x 5-1/4"	3-1/2"	1'	Out	Direct Thermal/Glove-Friendly Permanent	Red Border	500/RL-2RLS/CS
LDWHP127			3'		Direct Thermal/Glove-Friendly Permanent	Pink Border	500/RL-2RLS/CS
LDWHP128			3'		Direct Thermal/Glove-Friendly Permanent	White	500/RL-2RLS/CS
LDWHG01			3'		Direct Thermal/Glove-Friendly Permanent	Pink Border	1,000/RL
LDWHP71			3'		Direct Thermal/Glove-Friendly Permanent	White	1,000/RL
LMED61			3'		Direct Thermal IR/Glove-Friendly Permanent	White	1,000/RL-2RLS/CS
LDWHP711			1'		Direct Thermal/Glove-Friendly Permanent	White	500/RL-2RLS/CS

Cerner Labeling Solutions



For a complete selection of Cerner Labeling Solutions, visit us online at www.pdchealthcare.com.



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TD-PN6-SL	4" x 4-1/2"	4-1/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-2RLS/BX
DPS-PN6-SL			3"		Direct Thermal IR/Glove-Friendly Permanent	White	1,000/RL-2RLS/BX
TD-PN6-SLZ			1-1/2"		Direct Thermal/Glove-Friendly Permanent	White	600/RL-4RLS/BX
TD1-PN6-SL			1"		Direct Thermal/Glove-Friendly Permanent	White	500/RL-4RLS/BX
TD-PN6-SLBD			3"		Direct Thermal IR/Glove-Friendly Permanent	White	1,000/RL-2RLS/BX
DPS-PN6-SLBD			3"		Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-2RLS/BX
TD-PN6-SLZBD			1-1/2"		Direct Thermal/Glove-Friendly Permanent	White	600/RL-4RLS/BX
TD1-PN6-SLBD			1"		Direct Thermal/Glove-Friendly Permanent	White	500/RL-4RLS/BX
TD-PN6-SLY			3"		Direct Thermal/Glove-Friendly Permanent	Yellow	1000/RL-2RLS/BX
LPNSY1			1"		Direct Thermal/Glove-Friendly Permanent	Yellow	350/RL-2RLS/BX
TD-PN6-SLYBD			3"		Direct Thermal/Glove-Friendly Permanent	Yellow	1,000/RL-2RLS/BX
LDWHP113			3"		Direct Thermal/Glove-Friendly Permanent	Green	1,000/RL-2RLS/CS
LDWHP137			3"		Direct Thermal/Glove-Friendly Permanent	Light Pink	1,000/RL-2RLS/CS
LPN6Y			3"		Direct Thermal IR/Glove-Friendly Permanent	Yellow	1,000/RL-2RLS/CS
TD-PN6A-SL			3"		Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-2RLS/BX
TD-PN6A-SLBD			3"		Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-2RLS/BX
TD-PN6B-SL			3"		Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-2RLS/BX
DPS-PN6			3"		Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-2RLS/BX

LABORATORY

LABORATORY


CERNER SYSTEM LABELS

Cerner Labeling Solutions




For more
Slide Labels, see
pages 128-129.

Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
LDWHP181	2" x 1-1/8"	2-1/8"	1"	Out	Direct Thermal/Glove-Friendly Permanent	2,500/RL
TD-PN7-S	2" x 1-1/8"		3"		Direct Thermal Synthetic/Extra Permanent	5,000/RL-6RLS/BX









Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TD-PN8	2" x 1-1/4"	2-1/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	3,000/RL-6RLS/BX
TD-PN8-S			3"		Direct Thermal Synthetic/Extra Permanent	3,000/RL-6RLS/BX
TD1-PN8			1"		Direct Thermal/Glove-Friendly Permanent	1,000/RL-6RLS/BX
TD-PN8-BD			3"		Direct Thermal/Glove-Friendly Permanent	3,000/RL-6RLS/BX
TD1-PN8-BD			1"		Direct Thermal/Glove-Friendly Permanent	1,000/RL-6RLS/BX
LPN8S			3"		Direct Thermal IR/Glove-Friendly Permanent	3,000/RL-6RLS/CS
LDWHP148			3"		Direct Thermal Synthetic/Glove-Friendly Permanent	5,000/RL

BD, BD Logo, and V-Notch™ trademarks are the property of Becton, Dickinson and Company. ©2005 BD.





Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TD-PN12-SL	2" x 2-3/16"	2-1/8"	3"	Out	Direct Thermal/ Glove-Friendly Permanent	1,000/RL-4RLS/BX
DPS-PN12-SL					Direct Thermal IR/ Glove-Friendly Permanent	

Cerner Labeling Solutions



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TD-PN4-BD	4" x 2-3/16"	4-1/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	White	2,600/RL-2RLS/BX
LDTSYP1	4" x 2-3/16"	4-1/8"	3"		Direct Thermal/Glove-Friendly Permanent	Yellow	2,850/RL-2RLS/CS
TD-PN4-BDS	4" x 2-3/16"	4-1/8"	3"		Direct Thermal Synthetic/Extra Permanent	White	2,600/RL-2RLS/BX
DPS-PN4-BD	4" x 2-3/16"	4-1/8"	3"		Direct Thermal IR/Glove-Friendly Permanent	White	2,600/RL-2RLS/BX
TD-PN4-NTCH	4" x 2-3/16"	4-1/8"	3"		Direct Thermal/Glove-Friendly Permanent	White	2,600/RL-2RLS/BX
LDTSYP84	4" x 2-3/16"	4-1/8"	1"		Direct Thermal/Glove-Friendly Permanent	Pink Border	800/RL-2RLS/CS
MCCG-PN4W	4" x 2-3/16"	4-1/8"	1"		Direct Thermal/ Glove-Friendly Permanent	White	1,000/RL
TD-PN4-S	5-1/4" x 2-1/2"	5-1/2"	Fanfold/EDP		Direct Thermal/Glove-Friendly Permanent	White	5,000/BX






Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TD-PN5	2-1/4" x 2-1/2"	2-3/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	White	2,000/RL-4RLS/BX
LDYHP75					Direct Thermal/Glove-Friendly Permanent	Yellow	2,000/RL-4RLS/CS
LDRHP75					Direct Thermal/Glove-Friendly Permanent	Red Border	2,000/RL-4RLS/CS
DPS-PN5					Direct Thermal IR/Glove-Friendly Permanent	White	2,000/RL-4RLS/BX
TD-PN5-BD					Direct Thermal/Glove-Friendly Permanent	White	2,000/RL-4RLS/BX




LABORATORY

CERNER SYSTEM LABELS

Cerner Labeling Solutions





Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TD-PN-115162	2" x 1-15/16"	2-1/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	2,700/RL-6RLS/BX
TD-PN115162A			3"		Direct Thermal/Glove-Friendly Permanent	2,700/RL-6RLS/BX
TD1PN115162A			1"		Direct Thermal/Glove-Friendly Permanent	1,000/RL-6RLS/BX
LDWHPA166			3"		Direct Thermal/Glove-Friendly Permanent	3,000/RL




Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TD-PN3-BD	7-1/2" x 2-3/16"	7-5/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	2,400/RL
TD-PN3-NTCH	7-1/2" x 2-3/16"	7-5/8"				
LCERP29	8-1/2" x 2-1/2"	8-3/4"				


Cerner Labeling Solutions





Product Number	Label Size	Carrier Width	Core Size	Material/Adhesive	Qty.
TD-MPN-122	8-15/32" x 1-3/8"	8-1/2"	Fanfold	Direct Thermal/Glove-Friendly Permanent	4,000/BX
DPS-MPN-01	8-1/2" x 1-3/4"			Direct Thermal IR/Glove-Friendly Permanent	5,000/BX



For a complete selection of Lab Labels, visit us online at www.pdchealthcare.com.



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
LDWHP553	4" X 2-3/4"	4-1/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	2,000/RL
TDS-1142	2" x 1-1/4"	2-1/8"			Direct Thermal Synthetic/Extra Permanent	2,000/RL-4RLS/BX








Product Number	Label Size	Carrier Width	Core Size	Material/Adhesive	Qty.
DPS-L1A	4-1/4" x 11/12"	9-1/2"	Fanfold/EDP	Smudge Proof Flexo/Glove-Friendly Permanent	5,000/BX
DPS-L1TWIN				Smudge Proof Litho/Glove-Friendly Permanent	

LABORATORY

SUNQUEST SYSTEM LABELS


Sunquest Labeling Solutions





Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TD-SQ4	4-1/8" x 1-3/16"	4-3/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	White	4,300/RL-2RLS/BX
TD1-SQ4			1"		Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-6RLS/BX
TD1-SQ4-BD			1"		Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-6RLS/BX
TD-SQ4-XP			3"		Direct Thermal/ Glove-Friendly Extra Permanent	White	4,300/RL-2RLS/BX
TD-SQ4-XPBD			3"		Direct Thermal/ Glove-Friendly Extra Permanent	White	4,300/RL-2RLS/BX
DPS-SQ4			3"		Direct Thermal IR/Glove-Friendly Permanent	White	4,300/RL-2RLS/BX
TD-SQ4Y			3"		Direct Thermal/Glove-Friendly Permanent	Yellow	4,300/RL-2RLS/BX
TD-SQ4Y-XP			3"		Direct Thermal/ Glove-Friendly Extra Permanent	Yellow	4,300/RL-2RLS/BX
DPS-SQ4Y			3"		Direct Thermal IR/Permanent	Yellow	4,300/RL-2RLS/BX
TD-SQ4-BD			3"		Direct Thermal/Glove-Friendly Permanent	White	4,300/RL-2RLS/BX
TD-SQ4Y-BD			3"		Direct Thermal/Glove-Friendly Permanent	Yellow	4,300/RL-2RLS/BX
TD-SQ4YB			3"		Direct Thermal/Glove-Friendly Permanent	Yellow Border	4,300/RL-2RLS/BX
TD-SQ4YB-BD			3"		Direct Thermal/Glove-Friendly Permanent	Yellow Border	4,300/RL-2RLS/BX
TD-SQ4B			3"		Direct Thermal/Glove-Friendly Permanent	Blue Border	4,300/RL-2RLS/BX
TD-SQ4B-BD			3"		Direct Thermal/Glove-Friendly Permanent	Blue Border	4,300/RL-2RLS/BX
DPS-SQ4B			3"		Direct Thermal IR/Glove-Friendly Permanent	Blue Border	4,300/RL-2RLS/BX
TD-SQ4R			3"		Direct Thermal/Glove-Friendly Permanent	Red Border	4,300/RL-2RLS/BX
TD-SQ4R-BD			3"		Direct Thermal/Glove-Friendly Permanent	Red Border	4,300/RL-2RLS/BX
TD-SQ4R-XP			3"		Direct Thermal/ Glove-Friendly Extra Permanent	Red Border	4,300/RL-2RLS/BX
DPS-SQ4R			3"		Direct Thermal/Glove-Friendly Permanent	Red Border	4,300/RL-2RLS/BX
TD-SQ4-SL			3"		Direct Thermal/Glove-Friendly Permanent	White	4,300/RL-2RLS/BX
TD-SQ4-SLBD	3"	Direct Thermal/Glove-Friendly Permanent	White	4,300/RL-2RLS/BX			

Sunquest Labeling Solutions



The following labels are designed specifically for Sunquest Information Systems. The BD symbol in this section indicates a V-Notch™ label, which are developed in conjunction with BD to prevent misidentification errors from unaligned labels, while facilitating accurate bar code scanning.



V-Notch™



V-Notch™



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TD-SQBLSTR	4-1/8" x 1-3/16"	4-3/16"	1-1/2"	In	Direct Thermal/Glove-Friendly Permanent	White	1,200/RL-8RLS/BX
TD-SQBLTR-XP		4-3/16"	1-1/2"	In	Direct Thermal/ Glove-Friendly Extra Permanent	White	1,200/RL-8RLS/BX
DPS-SQBLSTR		4-3/16"	1-1/2"	In	Direct Thermal IR/Glove-Friendly Permanent	White	1,200/RL-8RLS/BX
TD-SQBLSTRBD		4-3/16"	1-1/2"	In	Direct Thermal/Glove-Friendly Permanent	White	1,200/RL-8RLS/BX
DPS-SQBLSRED		4-3/16"	1-1/2"	In	Direct Thermal IR/Glove-Friendly Permanent	Red Border	1,200/RL-8RLS/BX
TD-SQBLSRED		4-3/16"	1-1/2"	In	Direct Thermal/Glove-Friendly Permanent	Red Border	1,200/RL-8RLS/BX
TD-SQBLSR-BD		4-3/16"	1-1/2"	In	Direct Thermal/Glove-Friendly Permanent	Red Border	1,200/RL-8RLS/BX
TD1-SQ4A		4-3/8"	1"	Out	Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-6RLS/BX



"Each year more than 7 billion laboratory tests are performed in North America, influencing an estimated 70% of medical decisions."



Implementation Guide for the use of Bar Code Technology in Healthcare published by HIMSS

PDC Healthcare Laboratory Labels provide positive identification for superior patient safety.

LABORATORY

SUNQUEST SYSTEM LABELS


Sunquest Labeling Solutions





V-Notch™




Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TD-SQBL	4-1/8" x 1-3/16"	4-3/8"	1-1/2"	In	Direct Thermal/Glove-Friendly Permanent	1,200/RL-8RLS/BX
DPS-SQBLS-SL		4-5/16"			Direct Thermal IR/Glove-Friendly Permanent	
TD-SQBLS-SL		4-5/16"			Direct Thermal/Glove-Friendly Permanent	
TD-SQBLSL-BD		4-5/16"			Direct Thermal/Glove-Friendly Permanent	




Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
DPS-MSQ-22	3-1/16" x 1-7/16"	3-1/8"	3"	Out	Direct Thermal IR/Glove-Friendly Permanent	1,200/RL-5RLS/BX
TD-MSQ-22					Direct Thermal/Glove-Friendly Permanent	



V-Notch™






Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TD-MSQ-23	3-1/16" x 1-7/16"	3-1/8"	1"	Out	Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-5RLS/BX
TD-MSQ-23BD					Direct Thermal/Glove-Friendly Permanent	White	
TD-MSQ23-R					Direct Thermal/Glove-Friendly Extra Permanent	Red Border	



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TTRS-23	3" x 2"	3-1/8"	1"	Out	Thermal Transfer Synthetic/Glove-Friendly Permanent	1,000/RL

Sunquest Labeling Solutions



Product Number	Sheet Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
DPS-SQ6	8-7/8" x 6-7/8"	11"	Fanfold/EDP	—	Smudge Proof Flexo/Glove-Friendly Permanent	1,000/BX
DPS-SQ7	8" x 5-3/8"	9"	Fanfold/EDP	—	Smudge Proof Flexo/Glove-Friendly Permanent	1,000/BX
TD1A-SQ7	3-3/4" x 6-5/8"	3-25/32"	1-1/2"	In	Direct Thermal/Glove-Friendly Permanent	400/RL-8RLS/BX

LABORATORY

Save Time and Money!




ValuePlus™ is the only standardization service of its kind that combines best practices for patient safety, regulatory compliance, and fiscal control. With unmatched personal attention, PDC Healthcare can guide you through a quick and convenient procurement process.

For details, visit us online at www.pdchealthcare.com


LABORATORY

SUNQUEST SYSTEM LABELS


Sunquest Labeling Solutions



Product Number	Label Size	Carrier Width	Core Size	Material/Adhesive	Qty.
DPS-SQ2	2-5/8" x 13/16"	4-3/4"	Fanfold/EDP	Smudge Proof Flexo/ Permanent	5,000/BX




Product Number	Sheet Size	Carrier Width	Core Size	Material/Adhesive	Qty.
DPS-SQ3	12" x 5"	12"	Fanfold/EDP	Smudge Proof Flexo/Glove-Friendly Permanent	1,000/BX



For a complete selection of Lab Labels, visit us online at www.pdchealthcare.com.


Product Number	Sheet Size	Carrier Width	Core Size	Material/Adhesive	Qty.
DPS-SQ5	8" x 3-7/16"	9"	Fanfold/EDP	Smudge Proof Flexo/Glove-Friendly Permanen	2,000/BX

CPSI Labeling




Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TDWHP01	3-1/2" x 15/16"	3-5/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	White	5,000/RL
TDWXP01				Out	Direct Thermal/Extra Permanent	White	5,000/RL
TDRED01				Out	Direct Thermal/Glove-Friendly Permanent	Red Border	5,000/RL
TDYLP01				Out	Direct Thermal IR/Permanent	Yellow	5,000/RL
TDWHP03				Out	Thermal Transfer/Glove-Friendly Permanent	White	5,000/RL
TD-CPSI-312A				In	Direct Thermal/Glove-Friendly Permanent	White	5,000/RL-2RLS/BX
TD-CPSI312BD	Out	Direct Thermal/Glove-Friendly Permanent	White	5,000/RL-2RLS/BX			

VA Labeling Solutions



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
DPS-CHCS	1-1/2" x 3-3/4"	1-5/8"	3"	Out	Direct Thermal IR/Glove-Friendly Permanent	1,000/RL-4RLS/BX





Product Number	Label Size	Carrier Width	Core Size	Material/Adhesive	Qty.
DPS-VA1	8-3/16" x 1-5/16"	9-1/4"	Fanfold	Smudgeproof Flexo/Permanent	5,000/BX

LABORATORY


VA SYSTEM & DYMO PRINTER LABELS

VA Labeling Solutions





Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
DPS-VA2	2-3/8" x 4"	2-1/2"	3"	Out	Direct Thermal IR/ Glove-Friendly Permanent	800/RL-4RLS/BX
TD-VA2	2-3/8" x 4"	2-1/2"			Direct Thermal/ Glove-Friendly Permanent	800/RL-4RLS/BX
TD-VA2-BD	2-3/8" x 4"	2-1/2"			Direct Thermal/ Glove-Friendly Permanent	800/RL-4RLS/BX
LVA43	2-3/8" x 4"	2-1/2"			Direct Thermal/ Glove-Friendly Permanent	1,400/RL-10RLS/CS
LVA55	3-3/8" x 4"	3-1/2"			Direct Thermal/ Glove-Friendly Permanent	800/RL
LVA46	3-1/2" x 2-7/8"	3-3/4"			Direct Thermal/ Glove-Friendly Permanent	1,700/RL-4RLS/CS

For Dymo/Costar Printers



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
LDWHP103	3/4" x 2"	1-7/32"	1"	Out	Direct Thermal IR/Glove-Friendly Permanent	500/RL
SJCSLP10	1" x 2-1/8"	1-7/32"			Direct Thermal/Glove-Friendly Permanent	500/RL
SJCSLP11	2-1/4" x 1-1/4"	2-7/16"			Direct Thermal/Glove-Friendly Permanent	1,000/RL

SMS Labeling Solutions









LABORATORY

Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TD-MSMS-82	4-1/4" x 6"	4-5/16"	3"	In	Direct Thermal/Glove-Friendly Permanent	White	500/RL-4RLS/BX
TD-MSMS-82BD		4-5/16"		In	Direct Thermal/Glove-Friendly Permanent	White	500/RL-4RLS/BX
TD-MSMS-82P		4-5/16"		In	Direct Thermal/Glove-Friendly Permanent	Pink	500/RL-4RLS/BX
TD-MSMS-82Y		4-5/16"		In	Direct Thermal/Glove-Friendly Permanent	Yellow	500/RL-4RLS/BX
LSMSMP	4-1/4" x 6"	4-3/8"	3"	In	Direct Thermal IR/Glove-Friendly Permanent	White	500/RL-4RLS/CS
TD-SMS-2A		4-3/8"		Out	Direct Thermal/Glove-Friendly Permanent	White	870/BX-2RLS/BX
LDWHP12		4-3/8"		Out	Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-2RLS/CS
LDWHP12B		4-3/8"		Out	Direct Thermal/Glove-Friendly Permanent	Blue Border	1,000/RL-2RLS/CS
LDWHP12R		4-3/8"		Out	Direct Thermal/Glove-Friendly Permanent	Red Boder	1,000/RL-2RLS/CS
LDWHP12Y		4-3/8"		Out	Direct Thermal/Glove-Friendly Permanent	Yellow Border	1,000/RL-2RLS/CS

LABORATORY

SOFTLAB / HEALTH SCIENCE SYSTEM LABELS



SoftLab Labeling Solutions



For a complete selection of Lab Labels, visit us online at www.pdchealthcare.com.




Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TD-SL1212	2-1/2" x 1"	2-5/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	White	4,000/RL-2RLS/BX
TD-SL1212Y	2-1/2" x 1"	2-5/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	Yellow	4,000/RL-2RLS/BX
TD1B-1212XP	2-1/2" x 1"	2-5/8"	1-5/8"	In	Direct Thermal/ Glove-Friendly Extra Permanent	White	2,000/RL-12RLS/BX
TD1B1212XPBD	2-1/2" x 1"	2-5/8"	1-5/8"	In	Direct Thermal/ Glove-Friendly Extra Permanent	White	2,000/RL-12RLS/BX
TD-SL3	3" x 1-3/8"	3-1/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	White	2,500/RL
TD1-SL3	3" x 1-3/8"	3-1/8"	1"	Out	Direct Thermal/Glove-Friendly Permanent	White	1,500/RL - 3RLS/BX
TD1-SL3-BD	3" x 1-3/8"	3-1/8"	1"	Out	Direct Thermal/Glove-Friendly Permanent	White	1,500/RL - 3RLS/BX
TD1-SL3RB	3" x 1-3/8"	3-1/8"	1"	Out	Direct Thermal/Glove-Friendly Permanent	Red Border	1,500/RL - 3RLS/BX
TD1-SL3YB	3" x 1-3/8"	3-1/8"	1"	Out	Direct Thermal/Glove-Friendly Permanent	Yellow Border	1,500/RL - 3RLS/BX

For Health Science Systems




Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
LDWHP141	4-1/8" x 3"	4-1/4"	1-1/2"	In	Direct Thermal/Glove-Friendly Permanent	White	500/RL-20RLS/CS
LDWHP83						Pink	

Citation Labeling Solutions




Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
LTWHP52	3-1/8" x 1-1/4"	3-1/4"	3"	Out	Thermal Transfer/Glove-Friendly Permanent	White	4,200/RL-2RLS/CS
LTYEP52					Thermal Transfer/Glove-Friendly Permanent	Yellow	
LTREP52					Thermal Transfer/Glove-Friendly Permanent	Red Border	
LCIT18					Direct Thermal/Glove-Friendly Permanent	White	
LCYIT18					Direct Thermal/Glove-Friendly Permanent	Red Border	

For MedLab Systems



Product Number	Label Size	Carrier Width	Core Size	Material/Adhesive	Qty.
DPS-MLA	8-3/16" x 1-5/16"	9-1/4"	Fanfold/ EDP	Smudgeproof Flexo/Glove-Friendly Permanent	5,000/BX
DPS-ML1				Light Weight/Glove-Friendly Permanent	

For Other Laboratory Information Systems






Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
LDWHP111	1-7/8" x 1-1/8"	2"	3"	Out	Direct Thermal/Glove-Friendly Permanent	5,000/RL
TD-12-XP	2" x 1"	2-1/8"	3"	Out	Direct Thermal/ Glove-Friendly Extra Permanent	4,000/RL-3RLS/BX
TD1-12-XP	2" x 1"	2-1/8"	1"	Out	Direct Thermal/ Glove-Friendly Extra Permanent	2,500/RL-4RLS/BX
TD1-114214	2-1/4" x 1-1/4"	2-3/8"	1"	Out	Direct Thermal/Glove-Friendly Permanent	1,135/RL-12RLS/BX
TD-114214-XP	2-1/4" x 1-1/4"	2-3/8"	1"	Out	Direct Thermal/ Glove-Friendly Extra Permanent	1,135/RL-12RLS/BX
GP-1-1212	2-1/2" x 1"	2-5/8"	1"	In	Direct Thermal/Glove-Friendly Permanent	1,300/RL-8RLS/BX
LDWHP161	3-1/2" x 1"	3-5/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	5000/RL




LABORATORY

OTHER LABORATORY INFORMATION SYSTEM LABELS

For Other Laboratory Information Systems




Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TI-138178	1-7/8" x 1-3/8"	2-3/8"	1-1/2"	Out	Direct Thermal IR/ Glove-Friendly Permanent	1,000/RL-6RLS/BX
TI-138218	2-1/8" x 1-3/8"	2-1/4"	3"	In		4,000/RL-5RLS/BX
TI-1212	2-1/2" X 1"	2-5/8"	3"	Out		4,000/RL-4RLS/BX



LDWHP765

Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
LDTSLT01	4" x 5-1/2"	4-1/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	Orange Border	1,000/RL
LDWHP138			3"		Direct Thermal/Glove-Friendly Permanent	White	1,000/RL-6RLS/CS
LDWHP368			3"		Thermal Transfer/Glove-Friendly Permanent	White	1,000/RL-6RLS/CS
LDWHP764			1"		Direct Thermal/Glove-Friendly Permanent	White	500/RL
LDWXP137			3"		Thermal Transfer/ Glove-Friendly Extra Permanent	White	1,000/RL-6RLS/CS
LDWHP765			1"		Direct Thermal/Glove-Friendly Permanent	Orange	500/RL



For Other Laboratory Information Systems



For a complete selection of Lab Labels, visit us online at www.pdchealthcare.com.



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
TD-14A	4" x 1"	4-1/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	White	5,000/RL-2RLS/BX
TD1-14A	4" x 1"	4-1/8"	1"		Direct Thermal/Glove-Friendly Permanent	White	1,500/RL-4RLS/BX
TD-14BR-P	4" x 1"	4-1/8"	3"		Direct Thermal/Glove-Friendly Permanent	Pink Border	5,000/RL
TD-14BR	4" x 1"	4-1/8"	3"		Direct Thermal/Glove-Friendly Permanent	White	5,000/RL
TD1-14	4" x 1"	4-1/8"	1"		Direct Thermal/Glove-Friendly Permanent	White	1,300/RL-4RLS/BX
TD3-14	4" x 1"	4-1/8"	3"		Direct Thermal/Glove-Friendly Permanent	White	3,300/RL-4RLS/BX




Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Color	Qty.
LDWHP162	4-1/4" x 4"	4-3/8"	3"	Out	Direct Thermal/Glove-Friendly Permanent	White	1,000/RL
LDWHPK013	2-1/2" x 2-3/4"	2-5/8"	3"		Thermal Transfer Synthetic/Glove-Friendly Permanent	Red	1,800/RL
LTPB622	1-3/4" X 3/4"	4"	3"		Thermal Transfer Piggyback/Glove-Friendly Permanent	White/Yellow	8,000/RL
LDWHP01	3" X 15/16"	3-1/8"	Fanfold	—	Direct Thermal/Glove-Friendly Permanent	White	7,000/CS

LABORATORY

POINT-OF-CARE LABELING SOLUTIONS

Labeling Solutions for Point-of-Care Safety




Increase patient safety and assure accuracy of identification at the point-of-care. Whether it's medication administration or specimen collection, PDC Healthcare labels help hospitals utilize the latest bar code scanning and wireless mobility products for medical care.

Labeling at the point-of-collection

- Provides an immediate method of positive identification
- Increases accuracy of labeling specimens and containers
- Can help reduce costs associated with retesting

When accuracy and speed of diagnosis count, you can rely on PDC Healthcare's point-of-care labeling solutions. For additional information, please contact your local PDC Healthcare sales representative.



Product Number	Label Size	Carrier Width	Core Size	Wound	Max OD	Material/Adhesive	Color	Qty.
TD75-1142	2" x 1-1/4"	2-1/8"	3/4"	Out	2"	Direct Thermal/Glove-Friendly Permanent	White	250/RL-8RLS/BX
TD75F-1142	2" x 1-1/4"					Direct Thermal Synthetic/Extra Permanent	White	250/RL-8RLS/BX
TD75-1142A	2" x 1-1/4"					Direct Thermal/Glove-Friendly Permanent	White	250/RL-8RLS/BX
TD75-1142-BD	2" x 1-1/4"					Direct Thermal/Glove-Friendly Permanent	White	250/RL-8RLS/BX
TD75-17162BD	2" x 1-7/16"					Direct Thermal/Glove-Friendly Permanent	White	250/RL-8RLS/BX
TD75-23162BD	2" x 2-3/16"					Direct Thermal/Glove-Friendly Permanent	White	200/RL-8RLS/BX
TD75-SQ-BD	2" x 2-3/16"					Direct Thermal/Glove-Friendly Permanent	White	200/RL-8RLS/BX
LDQL220	2" x 1"					Direct Thermal/Glove-Friendly Permanent	White	350/RL-36RLS/CS
LDQLP220	2" x 1"	Direct Thermal/Glove-Friendly Permanent	White	350/RL-36RLS/CS				
LDQLPY220	2" x 1"	Direct Thermal/Glove-Friendly Permanent	Yellow Border	350/RL-36RLS/CS				

ISBT Blood Labels


Reduce Errors and Improve Patient Safety with ISBT Compliant Labels

We can help you meet the ISBT Blood Bank standards with our complete line of blank and pre-printed labels.

Available in single quadrant, double quadrant and full-face configurations, our labels are precisely sized to run easily through your thermal printer while producing accurate bar codes and eye-readable information.



All our ISBT 128 Blood Bag Labels on this page are for direct blood bag application.




We also carry ribbon!
Call us at (800) 435-4242

Product Number	Description	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
59715118	Without Cut-out	4" x 4"	4-1/8"	3"	Out	Synthetic Thermal Transfer/ Glove-Friendly Permanent	1,000/RL
59715119	Without Cut-out, Kit	4" x 4"	4-1/8"	3"			1,000/RL-4RLS/CS, 1 RL RIBBON
59715040	With Cut-out	4" x 4"	4-1/8"	3"			1,000/RL
59715003	With Cut-out, Kit	4" x 4"	4-1/8"	3"			1,000/RL-4RLS/CS, 1 RL RIBBON
59715306	With Cut-out	4" x 2"	4-1/8"	3"			1,750/RL
59715307	Kit	4" x 2"	4-1/8"	3"			1,750/RL-4RLS/CS, 1 RL RIBBON
59715042	Without Cut-out	2" x 2"	2-1/8"	3"			2,700/RL
59715001	Kit	2" x 2"	2-1/8"	3"			2,700/RL-3RLS/CS, 1 RL RIBBON
59715288	With Cut-out	4" x 4"	4-1/8"	1"			500/RL-4RLS/BX

ABO/Rh Blood Group

Product Number	Description	ISBT 128 Bar code	Label Size	Material/Adhesive	Qty.
59715013	A Rh POSITIVE	6200	2" x 2"	Synthetic/ Glove-Friendly Permanent	500/RL
59715012	A Rh NEGATIVE	600			
59715039	B Rh POSITIVE	7300			
59715038	B Rh NEGATIVE	1700			
59715058	O Rh POSITIVE	5100			
59715057	O Rh NEGATIVE	9500			
59715023	AB Rh POSITIVE	8400			
59715022	AB Rh NEGATIVE	2800			




LABORATORY

MICROSCOPE SLIDE LABELING SOLUTIONS




PDC Healthcare microscope slide solutions are available in tape and label form. Choose from white paper or tag stock, with permanent or Glove-Friendly adhesive in roll or sheet forms. They can be provided blank, imprinted with your facility and department name, or color coded in the size and thickness indicated. Accession numbering is also available. Our StainerShield® (page 130) and MicroTrax® (page 131) thermal transfer label/ribbon combination products are durable and dependable solutions for identifying patient tests.

Paper Slide Labels







For a complete selection of Lab Labels, visit us online at www.pdchealthcare.com.



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TD3-7878	7/8" x 7/8"	4"	3"	Out	Direct Thermal/Glove-Friendly Permanent	5,000/RL-3RLS/BX
TTR-7878	7/8" x 7/8"	4"	1-1/2"		Thermal Transfer/Glove-Friendly Permanent	7,400/RL-3RLS/BX
TD-MM1	3-1/4" x 6"	3-3/8"	3"		Direct Thermal/Glove-Friendly Permanent	1,000/RL-8RLS/BX
DPS-MM1	3-1/4" x 6"	3-3/8"	3"		Direct Thermal IR/Glove-Friendly Permanent	




Product Number	Label Size	Carrier Width	Core Size	Material/Adhesive	Qty.
DPS-AP1	15/16" x 15/16"	7"	Fanfold/EDP	Smudgeproof Litho/Glove-Friendly Permanent	25,000/BX



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
SLIDE24	15/16" x 15/16"	4-1/4"	3"	Out	Direct Thermal/Glove-Friendly Permanent	20,000/RL
TD1-PN-1516		4"	1"	Out	Direct Thermal/Glove-Friendly Permanent	1,800/RL-4RLS/BX
TTR-PN-1516		4"	3"	In	Thermal Transfer/Glove-Friendly Permanent	7,400/RL-3RLS/BX
TTR-1516		4"	1-1/2"	Out	Thermal Transfer/Glove-Friendly Permanent	7,400/RL-3RLS/BX



Synthetic Slide Labels



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
SLIDE64	7/8" x 7/8"	4"	3"	Out	Thermal Transfer Synthetic/Permanent	20,000/RL
TTRS-3-1516	15/16" x 15/16"	3-1/16"	1"		Thermal Transfer Synthetic/Permanent	6,000/RL
SLIDE68	3-3/4" x 15/16"	3-7/8"	1"		Direct Thermal Synthetic/Permanent	1,800/RL

Stock & Custom Microscope Labels – Roll Form


Product Number	Description	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
MRS-78	Stock – Roll Form	7/8" x 7/8"	5"	1"	Out	Smudgeproof Litho/ Glove-Friendly Permanent	1,000/RL
MRSC-78	Custom – Roll Form*						



*Imprint can be on the top or bottom of the label.


Stock & Custom Microscope Labels – Sheet Form

Stock



MSS-38 MST-38 (Tag)

Custom



MSSC-38 MSTC-38 (Tag)

Product Number	Sheet Form	Label Size	Sheet Size	Labels/ Sheet	Material/Adhesive	Qty.
MSS-38	Stock	7/8" x 7/16"	5" x 7-9/10"	80	Smudgeproof Litho/Glove-Friendly Permanent	13 SHTS/BX
MST-38	Stock	7/8" x 7/16"		80	Tag/Glove-Friendly Permanent	13 SHTS/BX
MSS-78	Stock	7/8" x 7/8"		40	Smudgeproof Litho/Glove-Friendly Permanent	25 SHTS/BX
MST-78	Stock	7/8" x 7/8"		40	Tag/Glove-Friendly Permanent	25 SHTS/BX
MSSC-38	Custom*	7/8" x 7/16"		80	Smudgeproof Litho/Glove-Friendly Permanent	13 SHTS/BX
MSTC-38	Custom*	7/8" x 7/16"		80	Tag/Glove-Friendly Permanent	13 SHTS/BX
MSSC-78	Custom*	7/8" x 7/8"		40	Smudgeproof Litho/Glove-Friendly Permanent	25 SHTS/BX
MSTC-78	Custom*	7/8" x 7/8"		40	Tag/Glove-Friendly Permanent	25 SHTS/BX

*Imprint can be on the top or bottom of the label.

LABORATORY


MICROSCOPE SLIDE LABELING SOLUTIONS

StainerShield® Slide Labeling Solutions

StainerShield®
Slide Labels are Xylene
and Stain Resistant!

Tired of spending valuable time and money on misidentified specimens, redraw samples and repeat testing? StainerShield® slide labels are a simplified solution to the complex problem of properly identifying and tracking specimen slides in a laboratory environment. Lab technicians can now print on-demand patient-specific labels and apply them to tissue and specimen slides as the slides are produced, before they go through the stainer. No over laminate flap, UV ink, or etching system is needed.

These labels can be printed on-demand using standard desktop thermal printers. In one simple printing step, technicians can print slide labels with both text and bar codes – including 2D bar codes – that contain data directly from the lab's information system.




Product Number	Label Size	Carrier Width	Core Size	Wound	Max OD	Material/Adhesive	Qty.
TDSS1-1-7878	7/8" x 7/8"	1"	1"	In	4-13/32"	StainerShield® Direct Thermal/Glove-Friendly Permanent	1,800/RL-4RLS/BX
TDSS15-17878			1-1/2"		4-13/32"		
TDSS3-1-7878			3"		5"		



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TDSS1-2-7878	7/8" x 7/8"	2"	1"	In	StainerShield® Direct Thermal/ Glove-Friendly Permanent	3,600/RL-2RLS/BX
TDSS15-27878			1-1/2"			
TDSS3-2-7878			3"			



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TDSS1-4-7878	7/8" x 7/8"	4"	1"	In	StainerShield® Direct Thermal/ Glove-Friendly Permanent	7,200/RL-1RL/BX
TDSS3-4-7878	7/8" x 7/8"		3"			
TDSS1-4-1278	7/8" x 1/2"		1"			
TDSS3-4-1278	7/8" x 1/2"		3"			





Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TDSS1-4-1516	15/16" x 15/16"	3-7/8"	1"	In	StainerShield® Direct Thermal/Glove- Friendly Permanent	7,200/RL-1RL/BX
TDSS3-4-1516			3"			

Custom StainerShield® labels are available in:
yellow, pink & light blue.




Contact your PDC Healthcare Sales Representative
for more details. (800) 435-4242

IDguard® Slide Labeling Solutions




Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TD15A-1516	15/16" x 15/16"	4"	1-1/2"	Out	Direct Thermal/ Glove-Friendly Permanent	2,500/RL-4RLS/BX
TD3-7878A	7/8" x 7/8"		3"			5,000/RL-3RLS/BX



*Recommended for
Cerner Millennium Software Systems*

Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TD1-PN-1516	15/16" x 15/16"	3-7/8"	4"	Out	Direct Thermal/ Glove-Friendly Permanent	1,800/RL-4RLS/BX

MicroTrax® Slide Labeling Solutions



Labels for MicroTrax® Printers						
Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
MICROTRXLBL4	7/8" x 7/8"	4"	1-1/2"	In	Thermal Transfer Synthetic/Glove-Friendly Permanent	7,400/RL-3RLS/BX
Thermal Transfer Ribbons						
Product Number	Ribbon Size	For Use With		Formula		Qty.
MICROTRXRBN4	4.25" x 5,496'	MICROTRXLBL4 Labels		Wax (High Mark)		1EA
Kits: Label and Ribbon Combinations						
Product Number	Box Contents					
MCRTRXLBLRBN	(3) rolls of MICROTRXLBL4 labels and (1) MICROTRXRBN4 ribbon					

LABORATORY

MICROSCOPE SLIDE LABELING SOLUTIONS

Custom Imprinted Accession Numbered Slide Labels

CMSST-38

YOUR FACILITY NAME
CITY AND STATE

99 579

YOUR FACILITY NAME
CITY AND STATE

No. 1234

CMSST-78
CMSTT-78

Can also customize
15/16" x 15/16" labels.

Available Accession Numbered Prefixes

00	No.	M	MM
The year	The abbreviation for number	One letter	Two letters

Product Number	Label Size	Sheet Size	Labels/Sheet	Material/Adhesive	Style
CMSST-38	7/8" x 7/16"	5" x 7-9/10"	80	Smudgeproof Litho/Permanent	Thin Tablet
CMSST-78	7/8" x 7/8"		40	Smudgeproof Litho/Permanent	Thin Tablet
CMSTT-78	7/8" x 7/8"		40	100# Tag/Permanent	Thick Tablet

Minimum order: 10,000 labels

TML-99
00581

PDC Healthcare
27770 N. Entertainment Dr.
Valencia, CA 91355

Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
CVR-58	5/8" x 7/8"	3/4"	3"	Out	Laminated Synthetic/Glove-Friendly Permanent	1,000/RL

Minimum order: 10,000 labels

No. 144

CNB-1024-R (Butt Cut)

No. 144

CNS-0824-R (Stripped)

Product Number	Style	Label Size	Material/Adhesive	Qty.
CNB-1024-R	Roll	1-1/2" x 5/8"	Smudgeproof Flexo/Permanent	1,000/RL
CNB-1024-T	Sheet			1,000/BX
CNS-0824-R	Roll	1-1/2" x 1/2"		1,000/RL
CNS-0824-T	Sheet			1,000/BX

Minimum Order: 10 rolls or 10,000 labels


Custom Imprinted Accession Numbered Tape

60527


Product Number	Tape Imprint Size	Core Size	Material/Adhesive	Roll Size
CNR-10	1/2" x 1-1/2"	1"	Time Tape®/Removable	1/2" x 750"

Minimum Order: 12 rolls.


Consecutive Number Labels




CNL010 (P)



CNL020 (P)




CNL030 (P)



CNL040 (P)

(P) = Permanent Adhesive
(R) = Removable Adhesive



CNL050 (P)

Product Number	Label Size	Qty.
CNL010	1-1/4" x 5/16"	1,000/RL
CNL020		2,000/RL
CNL030		3,000/RL
CNL040		4,000/RL
CNL050		5,000/RL

Tamper-Evident Labels & Tape



Ideal for labeling various specimens and containers. These labels provide a cost-effective means of security where tamper-evident seals are necessary.




TRL-61 — Red

For a complete selection of Tamper-Evident Labels & Tape, visit us online at www.pdchealthcare.com.



TRL-2N — Red



TRL-1 — White/Red

Product Number	Label Size	Carrier Width	Core Size	Material/Adhesive	Qty.
TRL-61	6-1/2" x 1-1/2"	1-1/2"	1"	Acetate/Permanent	200/RL-1 RL/BX
TRL-2N	6" x 1/2"	6-1/2"	Sheeted	Acetate/Permanent	1,000/BX
TRL-1	4-1/8" x 1/2"	4-1/4"	Sheeted	Light Weight/Permanent	1,000/BX

TAMPER RESISTANT TAPE

HCL8006 — Red

Product Number	Tape Imprint Size	Roll Size
HCL8006	1" x 7"	1" x 1,296"


LABORATORY

COLOR CIRCLES & SPECIALTY TAPE

Circle Labels



Circle labels are ideal for a variety of uses in the laboratory – from alerting staff to conditions needing attention, to general direction of specimen flow and urgency.



Label Size	Core Size	Carrier Width	Material/Adhesive	Qty.
3/4" Circles	1"	7/8"	Smudgeproof Litho/Permanent	1,000/RL

Hi-Lo Temp Tape



This product can withstand temperatures as high as +550°F and as low as -384.8°F. The unique self-adhesive tape may be subjected to various testing equipment such as ovens, sterilizers, and incubators. Hi-Lo temp tape combines great strength with durability and has a textured surface that may be stamped or written on with a carbon pencil.

Product Number	Tape Width	Core Size	Material/Adhesive	Roll Size
AT-312	1/2"	1"	Temp Tape/Permanent	1/2" x 300"
AT-334	3/4"			3/4" x 300"
AT-301	1"			1" x 300"

Red Border Tape



Product Number	Tape Imprint Size	Core Size	Material/Adhesive	Roll Size	Rolls/Box
RB-534	3/4" x 1"	1"	Time Tape®/ Permanent	3/4" x 500"	16
RB-501	1" x 1-5/8"			1" x 500"	12
RB-5112	1-1/2" x 2-1/2"			1-1/2" x 500"	8
RB-502	2" x 3"			2" x 500"	6

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

Time Tape® —
the toughest, most
reliable labeling tape in
the industry!

Time Tape®

Stock Time Tape® 500" Rolls

Stock Time Tape® in 500" Rolls - Available in the 18 stock colors shown below			
Product Number	Tape Width	Core Size	Qty.
T-512	1/2"	1"	500*/RL
T-534	3/4"		
T-501	1"		
T-5112	1-1/2"		
T-502	2"		

Stock Time Tape® 2160" Rolls

Stock Time Tape® in 2160" Rolls - Available in the 18 stock colors shown below			
Product Number	Tape Width	Core Size	Qty.
T-1260	1/2"	3"	2160*/RL
T-3460	3/4"		
T-160	1"		
T-11260	1-1/2"		
T-260	2"		



To order, add the "color number" to the end of the product number. For example, 1/2" x 500" in green is T-512-3.

Time Tape®, 500" Rainbow Packs

Stock Time Tape® in 500" Rolls - Rainbow Pack			
Product Number	Tape Width	Colors per Box	Box Qty.
T-512-RP	1/2"	(5) White, (4) Yellow, (3) Blue, (3) Green, (3) Orange, (3) Red, (3) Pink	24 Rolls
T-534-RP	3/4"	(4) White, (2) Yellow, (2) Blue, (2) Green, (2) Orange, (2) Red, (2) Pink	16 Rolls
T-501-RP	1"	(2) White, (2) Yellow, (2) Blue, (2) Green, (1) Orange, (2) Red, (1) Pink	12 Rolls

Time Tape® is also available in a convenient Rainbow Pack, which contains assorted rolls of 500" tape. The Rainbow Pack is available in three widths and priced per box. Order in full box quantities.

CUSTOM COLORS

Time Tape® can also be manufactured to match your specific color requirements. When ordering custom imprinted Time Tape®, please specify the following:

1) Length of the roll and quantity:

<input type="checkbox"/> 500" <input type="checkbox"/> 2160"	1/2" W	3/4" W	1" W	1-1/2" W	2" W
Minimum Qty. (500")	24	24	24	24	24
Minimum Qty. (2,160")	12	8	6	8	6

2) Width of the roll:

1/2" 3/4" 1" 1-1/2" 2"

3) Repeat length of the imprint:

Repeat length is the distance from the beginning of one imprint to the beginning of the next imprint.

REPEAT REPEAT REPE
(The repeat length is = 1-1/2")

4) Exact wording needed:

For best results, include a sample or a sketch

5) Tape Color:

Choose from one of the 18 colors shown above.

6) Color Matching:

Any color not selected from the 18 shown above is subject to a \$25.00 ink matching charge. Add 20% for each additional imprint ink color.

YOUR IMPRINT HERE

7) Reverse Plates:

When a light colored type on a dark background is needed, there is an additional 25% tinting charge.

8) Perforation:

Add an additional 10% to the tape price for perforations between imprints.

9) Special Logos:

Custom logos or unique graphics/type are subject to an additional plate charge of \$35.00 for one color and \$50.00 for two colors.

10) Proof:

Upon written request with your purchase order, we will supply a one-color proof of your label. This proof must be returned within 30 days with your approval before the order can be processed. Proof orders that are cancelled will be billed for all production charges incurred prior to cancellation.

Call your PDC Healthcare representative at
(800) 435-4242 for more details.

LABORATORY

STAT / SPECIMEN LABELS

STAT Labels



59704101 – Fl. Red (P)



59702043 – Fl. Red (P)

Label Size	Roll Qty.
1-1/4" x 3/8"	1,000



59701020 – White/Red Border (P)

Label Size	Roll Qty.
1-1/2" x 1/2"	1,000

LSTAT39 – White (P)
Bar-codes 39

LSTAT128 – White (P) Bar-codes 128

Label Size	Roll Qty.
1-1/2" x 1/2"	250

Label Size	Roll Qty.
1-5/8" x 7/8"	1,000

59702048 – Fl. Red (P)




59704098 – Fl. Red (P)



59701116 – Fl. Red (P)

Label Size	Roll Qty.
2-1/4" x 7/8"	1,000

Specimen Labels



59704131 – Fl. Red (P)



59704123 – Fl. Yellow (P)

Label Size	Roll Qty.
2-1/4" x 7/8"	1,000

59704084 – Fl. Pink (P)

59702045 – Fl. Red (P)

59704107 – Fl. Yellow (P)

59704041 – Fl. Yellow (P)


59704126 – Fl. Yellow (P)

Label Size	Roll Qty.
1-1/4" x 3/8"	1,000

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

(P) = Permanent Adhesive
(R) = Removable Adhesive


Inventory Management/Lot Control Labels



59704156 — Fl. Red (P)



59702206 — Fl. Orange (P)




59702218 — Yellow (P)




59704160 — Green (P)



59713392 — Yellow (P)



59713383 — Fl. Green (P)



59704163 — Blue (P)

Label Size	Roll Qty.
1-5/8" x 7/8"	1,000


REVISED EXPIRATION

Date _____ Time _____


Label Size	Roll Qty.
2-1/4" x 7/8"	1,000

For a complete selection of Lab Labels, visit us online at www.pdhealthcare.com.

59702193 — White (P)



59704159 — Green (P)



59704155 — Fl. Red (P)

Label Size	Roll Qty.
3" x 1-3/4"	500

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

LABORATORY

LABORATORY

REFRIGERATION / BLOOD LABELS

Storage Condition Labels

REFRIGERATE	Store at Room Temperature	REFRIGERATE	FROZEN	<table border="1"> <tr><th>Label Size</th><th>Roll Qty.</th></tr> <tr><td>1-1/4" x 3/8"</td><td>1,000</td></tr> </table>	Label Size	Roll Qty.	1-1/4" x 3/8"	1,000
Label Size	Roll Qty.							
1-1/4" x 3/8"	1,000							
59712785 — Fl. Green (P)	59704138 — Blue (P)	59704139 — Fl. Yellow (P)	59704127 — Fl. Pink (P)					

DO NOT REFRIGERATE	REFRIGERATE	<table border="1"> <tr><th>Label Size</th><th>Roll Qty.</th></tr> <tr><td>2" x 3/4"</td><td>1,000</td></tr> </table>	Label Size	Roll Qty.	2" x 3/4"	1,000
Label Size	Roll Qty.					
2" x 3/4"	1,000					
59702057 — Fl. Green (P)	59704143 — Fl. Yellow (P)					

REFRIGERATE ON ARRIVAL

59704132 — Fl. Yellow (P)

Label Size	Roll Qty.
2-7/8" x 7/8"	1,000

REFRIGERATE

59714047 — Fl. Yellow (P)

Label Size	Roll Qty.
2-1/4" x 7/8"	1,000

Blood Labels

Blood labels on this page are not intended for direct application on blood bags. They are to be applied on the blood bag base label.

<div style="border: 1px solid black; padding: 5px; background-color: #f080f0;"> <p>ABO & Rh CONFIRMED</p> <p>_____</p> <p style="text-align: center;">Date / Tech.</p> </div> <p>59702016 — Fl. Pink (P)</p>	<table border="1"> <tr><th>Label Size</th><th>Roll Qty.</th></tr> <tr><td>2" x 3/4"</td><td>1,000</td></tr> </table>	Label Size	Roll Qty.	2" x 3/4"	1,000	<div style="border: 1px solid black; padding: 5px; background-color: #f080f0;"> <p>CORD BLOOD</p> </div> <p>59702223 — Fl. Pink (P)</p>	<table border="1"> <tr><th>Label Size</th><th>Roll Qty.</th></tr> <tr><td>1-1/4" x 3/8"</td><td>1,000</td></tr> </table>	Label Size	Roll Qty.	1-1/4" x 3/8"	1,000
Label Size	Roll Qty.										
2" x 3/4"	1,000										
Label Size	Roll Qty.										
1-1/4" x 3/8"	1,000										

<div style="border: 1px solid black; padding: 5px; background-color: #ffff00;"> <p>ABO GROUP CONFIRMED Rh OF NEGATIVE UNITS CONFIRMED</p> <p>_____</p> <p style="text-align: center;">DATE/TECH</p> </div> <p>59702012 — Yellow (P)</p>	<div style="border: 1px solid black; padding: 5px; background-color: #ff0000; color: white;"> <p>UNCROSSMATCHED BLOOD</p> </div> <p>59702066 — Fl. Red (P)</p>	<table border="1"> <tr><th>Label Size</th><th>Roll Qty.</th></tr> <tr><td>2-1/4" x 7/8"</td><td>1,000</td></tr> </table>	Label Size	Roll Qty.	2-1/4" x 7/8"	1,000
Label Size	Roll Qty.					
2-1/4" x 7/8"	1,000					

<div style="border: 1px solid black; padding: 5px; background-color: #ffff00;"> <p>CORD BLOOD</p> </div> <p>59702224 — Fl. Yellow (P)</p>	<div style="border: 1px solid black; padding: 5px; background-color: #ffff00;"> <p>CORD BLOOD</p> <p>Baby's Name _____ Sex M F</p> <p>Baby's ID # _____ Birthdate _____ Time _____</p> <p>Mother _____ Type _____ Nurse Initials _____</p> </div> <p>59702225 — Fl. Yellow (P)</p>	<table border="1"> <tr><th>Label Size</th><th>Roll Qty.</th></tr> <tr><td>2-7/8" x 7/8"</td><td>1,000</td></tr> </table>	Label Size	Roll Qty.	2-7/8" x 7/8"	1,000
Label Size	Roll Qty.					
2-7/8" x 7/8"	1,000					

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.



Hazard Communication Labeling Solutions

When dealing with any biohazard materials, Hazard Communication Labels from PDC Healthcare ensure safety and efficiency. In accordance with various federal and state regulations, these Hazard Communication Labels quickly identify containers, rooms and materials. Any Hazard Communication Label can be customized in the size and design that you require.



59704555 – Fl. Red (P)



BH-405 – Fl. Red (P)



59702960 – Fl. Red (P)



59702956 – Fl. Red (P)
59701557 – Fl. Orange (P)



59704557 – Fl. Red (P)
59704971 – Fl. Orange (P)



BH-55 – Fl. Red (P)



59712724 – Fl. Red (P)



BH-11 – Fl. Red (P)

Product Number	Label Size	Qty.
59704555	2" x 3"	500/RL
BH-405	2" x 3"	500/RL
59702960	2-1/4" x 7/8"	1,000/RL
59702956	3" x 3"	500/RL
59701557	3" x 3"	500/RL

Product Number	Label Size	Qty.
59704557	6" x 6"	50/PK
59704971	6" x 6"	50/PK
BH-55	5" x 5"	250/RL
59712724	1-1/4" x 3/8"	1,000/RL
BH-11	1" x 1"	1,000/RL

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

(P) = Permanent Adhesive
(R) = Removable Adhesive

LABORATORY

GENERAL PURPOSE LABELS / TAPE

Equipment Status Labels

DEFECTIVE
MAINTENANCE
REQUIRED

DATE: _____ TECH: _____
 WO# _____
 PROBLEM: _____

 STATUS: _____

59701652 — Fl. Red (R)

Label Size	Roll Qty.
2-1/2" x 2-1/2"	500

DEFECTIVE
DO NOT USE

Comments: _____

 DATE: _____ BY: _____

59701653 — Fl. Red (P)

CLEANED — CHECKED

By: _____ Date: _____


59701756 — White (R)

Label Size	Roll Qty.
2" x 3/4"	1,000

Radioactive Warning Tape



RA-1 - Yellow



RA-4 - Yellow

Product Number	Tape Imprint Size	Core Size	Material/Adhesive	Roll Size	Qty.
RA-1	1" x 3"	1"	Time Tape®/ Removable	1" x 500"	166 Imprints/500*RL-12RLS/BX
RA-4	3/4" x 3"			3/4" x 500"	166 Imprints/500*RL-16RLS/BX

Fomalin Labels

10% FORMALIN ☠

CAUTION: CONTAINS FORMALDEHYDE

Toxic by inhalation and if swallowed, irritating to the eyes, respiratory system and skin. May cause sensitization by inhalation or skin contact. Risk of serious damage to eyes. Potential cancer hazard. Repeated or prolonged exposure increases the risk.

59707959 — Fl. Red (P)

10% FORMALIN ☠

CAUTION: CONTAINS FORMALDEHYDE

Toxic by inhalation and if swallowed, irritating to the eyes, respiratory system and skin. May cause sensitization by inhalation or skin contact. Risk of serious damage to eyes. Potential cancer hazard. Repeated or prolonged exposure increases the risk.


59707946 — Fl. Pink (P)


Label Size	Roll Qty.
2-1/4" x 7/8"	1,000

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.

(P) = Permanent Adhesive
(R) = Removable Adhesive

Reagent Label




 Reagent _____
 Prep. Date _____
 Exp. Date _____
 Storage _____
 Room Refrig.
 Initials _____

LFWHG02 – White (P)

For a complete selection of Lab Labels, visit us online at www.pdchealthcare.com.

Label Size	Roll Qty.
3" x 1-5/8"	1,000

Other Preprinted Labels

SERUM

59704022 – Fl. Red (P)

PLASMA

59704043 – Fl. Green (P)

PLASMA

59704042 – Fl. Yellow (P)

LINE DRAW

59704023 – Fl. Pink (P)

EXTRA TUBE

59713388 – Fl. Yellow (P)

Label Size	Roll Qty.
1-1/4" x 3/8"	1,000

Signature _____
 Identifier _____
 Date _____ Time _____
 MV04FO3395P

MV04FO3395P – Fl. Orange (P)

Label Size	Roll Qty.
2-1/4" x 1"	420

BEDSIDE CHECK NOT DONE

Pt Name: _____

Location: _____

MV09FP1027

MV09FP1027 – Fl. Pink (R)

Label Size	Roll Qty.
4" x 2-5/8"	375

FASTING

1/2 HOUR

1 HOUR

2 HOUR

3 HOUR

4 HOUR

5 HOUR

59707939 – Fl. Yellow (P)


Label Size	Roll Qty.
1-1/4" x 1/2"	8 LABELS/SHT-100 SHTS/PK

See Inside Back Cover for a PMS Color Guide of Fluorescent Label Colors on this Page.


PDC Healthcare is your expert in custom labels

Select Your Label Size, Color and Text on page 143 and we'll print it!

From sketch...



..to printed version.



Simplify your job by creating your own custom message labels.

- **Save Time** – eliminate time spent writing the same message
- **Reduce Errors** – clearly communicate your message
- **Get Attention** – choose from available fluorescent and brightly-colored materials to print your message on so it doesn't get overlooked. See inside back cover for a PMS color guide of Fluorescent colors.
- **Low Minimum Quantity** – minimum quantity as low as 1 roll of 500 labels

To order today, call us **(800) 435-4242** or download the order form online at www.pdchealthcare.com/customlabels.

Specialty Custom Labels

With more than 55 years of experience, PDC Healthcare is your expert in designing and producing the exact healthcare labels to meet your facility's specific needs. Give us your specifications and we'll create the labeling solution for you at the most economic cost.

Our portfolio of labeling solutions includes:

Sterile Medication Labels

Computer Generated System Labels

Label/Form Combination Sheets

Bar-codes Labels

HIPAA Compliant Self-Shred IV Labels

Self-Laminating Labels

Property Asset Tags

Consecutive Numbers

Flags and Tags

Tapes

and More . . .




Tape




Property Asset Tags




Flags




HIPAA Compliant Self-shred IV Labels



Sterile Labels



Consecutive Numbers



Computer Generated System Labels

Call us **(800) 435-4242** for a **FREE** quote today!

Create your own custom printed message labels!

Pick a size and sketch your design within the appropriate label size below (use a separate form for each unique label). Length of message will determine type size.

1) Select Your Label Size _____

1. 3" x 1-3/4"

2. 1-5/8" x 7/8"

4. 1-1/2" x 1/2"

3. 1-1/4" x 3/8"

5. 1-1/4" x 5/16"

6. 3" x 7/8"

7. 2-7/8" x 7/8"

8. 2-1/4" x 7/8"

Sketch desired imprint within the sizes
 — One time typesetting charge of \$16.95.

	Product Number	Label Size	Roll Qty.
1	<input type="checkbox"/> 59726003	3" x 1-3/4"	500
2	<input type="checkbox"/> 59726007	1-5/8" x 7/8"	1,000
3	<input type="checkbox"/> 59726001	1-1/2" x 1/2"	1,000
4	<input type="checkbox"/> 59726000	1-1/4" x 3/8"	1,000
5	<input type="checkbox"/> 59726006	1-1/4" x 5/16"	1,000
6	<input type="checkbox"/> 59726009	3" x 7/8"	1,000
7	<input type="checkbox"/> 59726002	2-7/8" x 7/8"	1,000
8	<input type="checkbox"/> 59726008	2-1/4" x 7/8"	1,000

Additional options:

- Any size
- Up to 4 ink colors
- Any logos or graphics
(please provide a clear black and white copy)
- Transparent labels

Call (800) 435-4242 for a custom quote!

2) Select Your Label Color _____

- | | |
|--|---|
| <input type="checkbox"/> Fluorescent Green | <input type="checkbox"/> Light Green |
| <input type="checkbox"/> Fluorescent Orange | <input type="checkbox"/> Light Blue |
| <input type="checkbox"/> Fluorescent Pink | <input type="checkbox"/> Pink |
| <input type="checkbox"/> Fluorescent Yellow | <input type="checkbox"/> Yellow |
| <input type="checkbox"/> Fluorescent Red | <input type="checkbox"/> White |

See Inside Back Cover for a PMS Color Guide of
Fluorescent Label Colors on this Page.

3) Select Your Ink Color _____

- | | |
|---|---|
| <input type="checkbox"/> Reflex Blue | <input type="checkbox"/> Green (PMS 342) |
| <input type="checkbox"/> Black | <input type="checkbox"/> Red (PMS 185) |

*For fluorescent material, black ink only.

4) Select Your Border Color _____

- Red**
 Green
 No Border

5) Select Your Label Color _____

- Permanent** - Adheres firmly to most surfaces.
 Removable - Removes easily, leaving little or no residue.

6) Select Your Label Color _____

Quantity of rolls _____ (minimum 1 roll)

7) Select Your Font Style _____

- | | | |
|---|---|--|
| <input type="checkbox"/> Helvetica | <input type="checkbox"/> Times New Roman | <input type="checkbox"/> <i>Brush Script</i> |
| <input type="checkbox"/> Helvetica Bold | <input type="checkbox"/> Times New Roman Bold | |

YES, I would like to start saving today!

Ordered By _____ Title _____
 Facility _____ Dept. _____
 Street _____ P.O. Box _____
 City _____ State _____ Zip _____
 Email _____

Account No. _____ PO # _____
 Phone _____ Fax _____
 Check One: BILL ME VISA MASTERCARD AMEX
 Card Number _____ "VCODE" _____
 Expiration Date _____ Print Name _____
 Signature _____

