

Comfort wrist support provides relief from wrist injury!

- Extra space below the 5h MP for palmar spread
- Generously padded, lightweight nylon for breathability and comfort
- Stockinette lining for easy one-handed application

| Product Number | Description | Qty. |
|----------------|-------------|--------|
| ORT-PWBORTX | 8" - Right | 1 Each |
| ORT-PWBOLT | 8" - Left | 1 Each |

Sizes: Determined by wrist circumference. Please specify right or left. (XXS) XX-Small, (XS) X-Small, (S) Small, (M) Medium, (L) Large, (XL) X-Large. To select your size, replace the gray X with the size code shown above.



Premium Wrist and Thumb Support

Use this support when immobilization at the MP and CMC joints is needed. Especially effective for relief in cases of painful scaphoid injuries and severe Carpal Tunnel Syndrome.

| Product Number | Description | Qty. |
|----------------|-------------|--------|
| ORT-PWTORTX | Right | 1 Each |
| ORT-PWTOLT | Left | 1 Each |

Sizes: Measure wrist circumference. (XS) X-Small up to 6", (S) Small 6" - 7", (M) Medium 7" - 8", (L) Large 8" - 9", (XL) X-Large 9" - 10".

To select your size, replace the gray X with the size code shown above.



See our new website with new features!

- Access all PDC, St. John, Pharmex and TimeMed products on one website!
- My Account 360° empowers you with instant account information
- Easier online ordering saves time
- Streamlined navigation & advanced search help you find information fast
- User-friendly tools help you explore products
- Instant online support at your fingertips



ORTHOPEDECS

CAST WALKERS



OrthoPro™ Cast Walker with rigid metal struts

- **Rigid metal struts** – for more secure immobilization and durability
- **Uprights are bowed at the ankle** – to accommodate tender and swollen ankles
- **Pronounced rocker sole** – provides a more natural gait by naturally carrying patient's motion forward
- **Super low profile** – to reduce heel height discrepancy. Helps keep hips and knees in alignment
- **Padded double wide straps** – provide comfortable immobilization from top to bottom

| Product Number | Description | Qty. |
|----------------|---|--------|
| ORT-OPW0000XX | OrthoPro™ Cast Walker with Metal Struts | 1 Each |

Sizes: (S) **Small**: fits Women's shoe size 4 - 8-1/2, Men's shoe size 3 - 7-1/2. (ML) **Medium/Large**: fits Women's shoe size 9 - 13-1/2, Men's shoe size 8 - 12-1/2. (XL) **X-Large**: fits Men's shoe size 13 - 16-1/2.

To select your size, replace the gray X with the size code shown above.



Lite-Foot™ Cast Walker is the lightest weightcast walker available

- **Lightweight** – Only 29 ounces for the small.
- **Flexible composite struts** – can be spread open to allow easy application and secure support.
- **Wider design** – more space between the struts accommodates swollen ankles and extra bandaging if necessary
- **More stable** – The wider, flatter sole has a larger footprint, resulting in a more stable gait. Non-slip tread outsole
- **Thicker cushioned heel** – absorbs impact at heel strike to increase comfort and reduce lower leg stress

| Product Number | Description | Qty. |
|----------------|----------------------------------|--------|
| ORT-LFX | Lite Foot™ Cast Walker | 1 Each |
| ORT-LTX | Lite Foot™ Cast Walker - Low Top | 1 Each |

Sizes: (S) **Small**: fits Women's shoe size 4 - 8-1/2, Men's shoe size 3 - 7-1/2. (M) **Medium/Large**: fits Women's shoe size 9 - 13-1/2, Men's shoe size 8 - 12-1/2. (X) **X-Large**: fits Men's shoe size 13 - 16-1/2.

To select your size, replace the gray X with the size code shown above.



High top all-weather cast boot

- Covers the ankle and conforms to both synthetic and plaster casting materials
- Padded angled rocker sole absorbs shock from heel strike to protect hips and knees for a more natural gait

| Product Number | Description | Qty. |
|----------------|-------------------------|--------|
| CB-PCX | High Top Poly Cast Boot | 1 Each |

Size: (PED) **Pediatric**, (XS) **X-Small** 7-1/2" L x 3-5/8" W, (S) **Small** 9-3/8" L x 4-1/4" W, (M) **Medium** 10-1/2" L x 4-7/4" W, (L) **Large** 11" L x 4-7/8" W, (XL) **X-Large** 12-5/8" L x 5-3/8" W.

To select your size, replace the gray X with the size code shown above.



Semi-rigid sole for a natural gait

- Padded nylon mesh upper for breathability
- Reinforced heel counter for a snug fit
- One buckle contact closure for easy application and fit

| Product Number | Description | Qty. |
|----------------|--------------------|--------|
| SPO-242X | Super Post-Op Shoe | 1 Each |

Size: (PED) **Pediatric**, (WS) **Women's Small**, (WM) **Women's Medium**, (WL) **Women's Large**, (MS) **Men's Small**, (MM) **Men's Medium**, (ML) **Men's Large**, (MXL) **Men's X-Large**

To select your size, replace the gray X with the size code shown above.

Synthetic Spee-D-Wrap™ Cast Padding

- Synthetic material will not absorb moisture or support bacterial growth that leads to maceration
- Tears easily — eliminating undue constriction
- Conforms easily — avoiding uncomfortable wrinkles
- Stays cool and absorbent under both synthetic and plaster casts for versatility and convenience
- Available in 4 yard lengths

| Product Number | Description | Qty. |
|----------------|-------------|--------------|
| ORT-876 | 2" Width | 12 Rolls/Bag |
| ORT-877 | 3" Width | 12 Rolls/Bag |
| ORT-878 | 4" Width | 12 Rolls/Bag |
| ORT-879 | 6" Width | 6 Rolls/Bag |



Synthetic Spee-D-Wrap™ Stockinette

Our Spee-D-Wrap™ stockinette has unusually dense stitching and ribbed knitting — assuring superior memory and support.

- Acts as the perfect cushion for preventing casts from rubbing against the skin
- "Finishes" the cast ends

| Product Number | Description | Qty. |
|----------------|-------------|--------------|
| ORT-895 | 3" Width | 25 Yards/Box |
| ORT-896 | 4" Width | |
| ORT-900 | 6" Width | |



Aquashield® — the most effective watertight cast protectors available!

Available for Half and Full Arm



AQSHA

- **Slip Resistant** — The non-slip sole on leg models provides excellent traction on wet and slippery surfaces.
- **Watertight** — A continuous seal allows the Aquashield® to be completely submersible, even when swimming!
- **Indestructible** — Tough polyurethane won't tear or rip. It's so tough it will outlast the cast or the bandaging!
- **Easy To Use** — Applies easily with one hand. No cutting or sizing required.



Available for Half and Full Leg

AQSHL

| Aquashield® for Half and Full Arm | | | |
|-----------------------------------|--------------------|------------------|--------|
| Product Number | Description | Circumference | Qty. |
| AQSHA | Small Half Arm | Less than 8-1/2" | 1 Each |
| AQRHA | Regular Half Arm | More than 8-1/2" | 1 Each |
| AQPFA | Pediatric Full Arm | - | 1 Each |
| AQSFA | Small Full Arm | Less than 10" | 1 Each |
| AQRFA | Regular Full Arm | More than 10" | 1 Each |

| Aquashield® for Half and Full Leg | | | |
|-----------------------------------|--------------------|-------------------|--------|
| Product Number | Description | Circumference | Qty. |
| AQSHL | Small Half Leg | Less than 11" | 1 Each |
| AQRHL | Regular Half Leg | More than 11" | 1 Each |
| AQPFL | Pediatric Full Leg | - | 1 Each |
| AQSFL | Small Full Leg | Less than 18-1/2" | 1 Each |
| AQRFL | Regular Full Leg | More than 18-1/2" | 1 Each |

Sizing: Measure the circumference of the limb above the cast. Pediatric sizes fit children ages 3 and under.