

Positioning Aids

- **STANDARD** – the most economical positioning aid – available in all designs. It is soft and artifact-free.
- **STAY-DRY** – a positioning aid with maximum resilience, strength and stability. It repels liquid – just wipe clean such substances as blood, barium and other agents. It is artifact-free and has a high drag coefficient which reduces slippage for easy patient management.


Chemical agents roll off Stay-Dry surface like rain on a roof. Easily wipes clean with soap and water!



6 Piece Starter Set

This set will get you started.

Qty.	Description	Size
1	Double Angle Wedge 55°	11" x 7" x 30"
1	Wedge 45°	9" x 9" x 30"
1	Circular Disc	1" x 7"
1	Rectangle	2" x 7" x 11"
2	Wedge 15°	3" x 7" x 11"




Product Number	Type	Size	Qty.
1500	Standard	Assorted	1 Set
SD2000	Stay-Dry		

Basic X-ray Set

This basic X-ray set is our most popular group saver!

Qty.	Description	Size
1	Circular Disc	1-1/2" x 7"
2	Wedge 16°	2-1/2" x 7" x 10"
2	Wedge 22°	3-1/2" x 7" x 8"
1	Finger Block	3" x 5" x 5"
2	Triangle 45°	10" x 10" x 10"
2	Triangle 45°	3" x 3" x 7-1/4"
1	Spinal Wedge 45°	7" x 7" x 24"




Product Number	Type	Size	Qty.
XR-SET	Standard	Assorted	1 Set
SDXR-SET	Stay-Dry		

Ortho Set

A perfect selection for your orthopedic needs!

Qty.	Description	Size
1	Spinal Wedge 45°	7" x 7" x 24"
2	Triangle 45°	10" x 10" x 10"
2	Wedge 18°	2-1/2" x 8" x 7"
1	Rectangle	2" x 7-1/2" x 9-1/2"
1	Rectangle	3" x 10" x 10"
1	Finger Block	3" x 5" x 5"
1	Knee Block	8" x 10" x 12"



Product Number	Type	Size	Qty.
OT-SET	Standard	Assorted	1 Set
SDOT-SET	Stay-Dry		

Cross Table Film Holder

- Size: 2-1/2" x 18" x 24"

Product Number	Type	Qty.
2121	Standard	1
SD2121	Stay-Dry	Each



Additional sizes available

Knee Block

- Size: 8" x 10" x 12"


Product Number	Type	Qty.
2033	Standard	1
SD2033	Stay-Dry	Each



Angular Bolster

- Size: 7-1/2" x 12-1/2" x 21"

Product Number	Type	Qty.
2038	Standard	1
SD2038	Stay-Dry	Each



Spinal Wedge 45°

- Size: 7" x 7" x 24"

Product Number	Type	Qty.
2037	Standard	1
SD2037	Stay-Dry	Each



Spinal Wedge 45°

- Size: 9" x 9" x 30"

Product Number	Type	Qty.
2084	Standard	1
SD2084	Stay-Dry	Each



Decubitus

- Size: 3" x 18" x 24"

Product Number	Type	Qty.
2059	Standard	1
SD2059	Stay-Dry	Each



PATIENT PROCEDURES

PROCEDURE TRAYS



Basic Procedure Tray

The paper-poly-paper fenestrated drapes guarantee impermeability, absorbency and low lint level.

Includes —

- 25G x 5/8" needle
- 18G x 1-1/2" needle
- 20cc LL syringe
- 5cc LL syringe
- 1% Lidocaine HCl, 5ml
- 3" x 3" gauze pads (3)
- PVP swabsticks (3)
- 18" x 18" fenestrated 2-3/4" window drape* with adhesive tab
- 17" x 19" towel*
- Adhesive bandage
- Sterile indicator

Product Number	Qty.
STJ-474	10/Case

*Disposable paper-poly-paper.



Arthrogram Procedure Tray

This pre-assembled sterile arthrogram tray is effectively organized to optimize procedural efficiency for accelerated patient flow.

Includes —

- 25G x 5/8" needle
- 21G x 1-1/2" needle
- 20G x 1-1/2" needle
- 18G x 1-1/2" needle
- 30cc LL syringe
- 10cc LL syringe
- 5cc LL syringe
- 1% Lidocaine HCl, 5ml
- 3" x 3" gauze pads (5)
- 18" x 18" fenestrated 2-3/4" window drape* with adhesive tab
- 17" x 19" towel*
- 20" extension set
- Applicators (3)
- Adhesive bandage
- Sterile indicator

Product Number	Qty.
STJ-471	10/Case

*Disposable paper-poly-paper.



Biopsy Tray

Specially designed to provide all the components necessary to prep and anesthetize the procedure area.

Includes —

- 25G x 5/8" needle
- 18G x 2" blunt needle
- 10cc LL syringes (2)
- 1% Lidocaine HCl, 5ml (2)
- #11 scalpel assembly
- 3" x 3" gauze pads (3)
- 8" ruler, paper
- 18" x 18" fenestrated 2-3/4" window drape* with adhesive tab
- 17" x 19" towel*
- PVP swabsticks (3)
- Alcohol swabstick
- Adhesive bandage
- Sterile indicator

Product Number	Qty.
MAM-378	10/Case

*Disposable paper-poly-paper.



Biopsy Tray

Helps speed and simplify biopsy procedures.

Includes —

- 25G x 5/8" needle
- 22G x 1-1/2" needle
- 18G x 2" blunt needle
- 10cc LL syringes (2)
- 1% Lidocaine HCl, 5ml (2)
- #11 scalpel, mini
- 3" x 3" 12 ply gauze pads (5)
- 8" ruler, plastic
- 16-1/2" x 18-1/2" towels, paper (2)
- 17-1/2" x 25" fenestrated 2-3/4" window drape* with adhesive tab
- PVP swabsticks (3)
- Alcohol swabstick
- Adhesive bandage
- Needle stop
- Foam block needle counter
- Sterile indicator

Product Number	Qty.
MAM-491	5/Case

*Disposable paper-poly-paper.



Breast Biopsy Localization Tray

Designed to facilitate breast biopsy localization — saving you time.

Includes —

- 25G x 5/8" needle
- 22G x 1-1/2" needle
- 18G x 1-1/2" needle
- 10cc LL syringe
- 5cc LL syringe
- 1% Lidocaine HCl, 5ml
- 4" x 4" gauze pads (5)
- 8" ruler, paper
- 18" x 18" fenestrated 2-3/4" window drape* with adhesive tab
- 17" x 19" towel*
- PVP swabsticks (3)
- Alcohol swabsticks (2)
- Wound closure strip 1/2" x 4"
- Sterile indicator

Product Number	Qty.
MAM-379	10/Case

*Disposable paper-poly-paper.



All Purpose Skin Prep Tray

Specially prepared to contain everything required to prep the skin prior to percutaneous procedures.

Includes —

- 25G x 5/8" needle
- 18G x 2" blunt needle
- 20cc LL syringe
- 5cc LL syringe
- 1% Lidocaine HCl, 5ml
- #11 scalpel, mini
- 3" x 3" gauze pads (5)
- 18" x 18" fenestrated 2-3/4" window drape* with adhesive tab
- 17" x 19" towels* (2)
- Applicators (3)
- Alcohol swabstick
- Adhesive bandage
- Sterile indicator

Product Number	Qty.
MAM-377	10/Case

*Disposable paper-poly-paper.