

PROTECTIVE APPAREL

APRONS/VESTS/SKIRTS

WRA-RE

Personalize It!
See p. 13

Hook and loop shoulder strap

Adjustable belt sits comfortably on your waist and evenly distributes the apron's weight over your hips

Protect-A-Tech™ Premium Back-Saver Apron

- Sewn-in tether that allows for vertical and lateral adjustments to any waist height
- Designed for even weight distribution during those long procedures
- Generous slit pocket has reinforced seams for added durability

Product Number	Weight
WRA-RE	Regular Lead
WRA-LI	Lightweight Lead
WRA-NO	No Lead

Size: Width (hip to hip) x Length (measured from sternal notch) (S) Small 20" x 24", (M) Medium 22" x 34", (L) Large 24" x 36", (XL) X-Large 26" x 38" and (2XL) 2X-Large 28" x 40". Specify size and color (p. 13) when placing your order. To order thyroid collar, see below.

Easy sizing — comfortable fit!

Personalize It!
See p. 13

Two-piece special procedure apron

Adjustable waist

Hook and loop closure with buckle.

Protect-A-Tech™ Premium Full Wrap Vest & Skirt

- Women's vest armholes are cut higher to protect breast tissue
- Men's vest has generous armholes for easier movement
- The skirt and kilt both have neoprene foam for long-term comfort without the bulk

Product Number	Description	Weight
TPF-RE	Vest, Women	Regular Lead
TPM-RE	Vest, Men	Regular Lead
SKF-RE	Skirt, Women	Regular Lead
SKM-RE	Kilt, Men	Regular Lead
TPF-LI	Vest, Women	Lightweight Lead
TPM-LI	Vest, Men	Lightweight Lead
SKF-LI	Skirt, Women	Lightweight Lead
SKM-LI	Kilt, Men	Lightweight Lead
TPF-NO	Vest, Women	No Lead
TPM-NO	Vest, Men	No Lead
SKF-NO	Skirt, Women	No Lead
SKM-NO	Kilt, Men	No Lead

	Vest: Woman Chest Size	Vest: Man Chest Size	Skirt: Woman Waist	Kilt: Man Waist
(S) Small	35"-36"	37"-38"	28"-29"	28"-29"
(M) Medium	37"-38"	39"-40"	30"-31"	30"-31"
(L) Large	39"-40"	41"-42"	34"-35"	34"-35"
(XL) X-Large	41"-42"	43"-44"	38"-39"	38"-39"
(2XL) 2 X-Large	43"-44"	45"-46"	—	—
(3XL) 3 X-Large	45"-46"	47"-48"	—	—

Unattached Thyroid Collar

Without binding

With binding

Attached Thyroid Collar

Attached collar stores in second apron pocket when not in use*.

Thyroid coverage at its most convenient and comfortable

- Optional pre-attached, binding-free thyroid collar
- An extra pocket will be added to your apron for the storage of your thyroid collar
- Optional detached thyroid collar with or without binding
- Color and weight options available

Weight	Product Number			
	Binding - Unattached	Binding - Attached	No Binding - Unattached	No Binding - Attached
Regular Lead	TTC-RST	TTC-RAT	TTC-RSN	TTC-RAN
Lightweight Lead	TTC-LST	TTC-LAT	TTC-LSN	TTC-LAN
No Lead	TTC-NST	TTC-NAT	TTC-NSN	TTC-NAN

Specify color (p. 13). *Attached collar must be ordered at the same time as apron.

PROTECTIVE APPAREL

APRONS/VESTS/SKIRTS

Protect-A-Tech™ Premium Hip Hugger Apron

- Standard hook and loop closure or front buckle closure
- Optional elasticized back adjusts the frontal protection for optimal weight distribution

Product Number	Description	Weight
HHA-RE	Standard Back	Regular Lead
HHA-E-RE	Elasticized Back	Regular Lead
HHA-LI	Standard Back	Lightweight Lead
HHA-E-LI	Elasticized Back	Lightweight Lead
HHA-NO	Standard Back	No Lead
HHA-E-NO	Elasticized Back	No Lead

Size: Width (hip to hip) x Length (measured from sternal notch) (S) Small 20" x 32", (M) Medium 22" x 34", (L) Large 24" x 36", (XL) X-Large 26" x 38", (2XL) 2 X-Large 28" x 40". Specify closure, size and color (see below) when placing your order. *To order thyroid collar, see p. 12.



HHA-E-RE

Protect-A-Tech™ Budget Apron

- Optional back strap for a comfortable fit
- Belt can be buckled or tied

Product Number	Weight
ECA-RE	Regular Lead
ECA-LI	Lightweight Lead
ECA-NO	No Lead

Size: Width (hip to hip) x Length (measured from sternal notch) (S) Small 20" x 32", (M) Medium 22" x 34", (L) Large 24" x 36", (XL) X-Large 26" x 38", (2XL) 2 X-Large 28" x 40". Specify size and color (see below) when placing your order. To order thyroid collar, see p. 12

Most Affordable!



ECA-RE

Back

Belt tied

Lead Equivalence Rating (Pb Eq.) for Protect-A-Tech™ Aprons

All aprons are manufactured to meet or exceed government safety requirements. With proper care, aprons last at least 2-3 years.

Kilovolt Potential	Regular Lead	Lightweight Lead	No Lead
80 kVp	0.51mm	0.52mm	0.55mm
100 kVp	0.52mm	0.53mm	0.53mm
120 kVp	0.49mm	0.52mm	0.48mm

Lightweight - Up to 10% lighter than regular lead.
 No Lead - Up to 18% lighter than regular lead.
 No Lead - Environmentally friendly - 100% recyclable.

Protect-A-Tech™ Aprons Personalization

Select from these colors at no extra charge. Contrasting binding color shown as stripe.


Monogram Available

Choose up to two lines; in either block or script style, in upper and lower case, or block style in all caps.



Nylon Fabric

Liquid repellent



All custom aprons are non-returnable.

Call (800) 435-4242 to place your order

PROTECTIVE APPAREL

PROTECTIVE APPAREL

APRONS/PORTABLE SHIELDS/APRON PEG BOARDS

Protect-A-Tech™ Half Apron

- Sewn-in tether that allows for vertical and lateral adjustments to any waist height
- Designed for even weight distribution during those long procedures
- Generous slit pocket has reinforced seams for added durability



Protect-A-Tech™ Semi-Guard				
Weight	Small (10" x 8-1/2")	Medium (12" x 12")	Large (24" x 18")	X-Large (24" x 24")
Regular Lead	ASG-RES	ASG-REM	ASG-REL	ASG-REX
Lightweight Lead	ASG-LIS	ASG-LIM	ASG-LIL	ASG-LIX
No Lead	ASG-NOS	ASG-NOM	ASG-NOL	ASG-NOX

Protect-A-Tech™ Demi-Guard Half Apron with Ring					
Weight	Small (8" x 10")	Medium (12" x 12")	Large (16" x 14")	X-Large (20" x 16")	Set with Rack*
Regular Lead	ADGR-RES	ADGR-REM	ADGR-REL	ADGR-REX	ADGR-RRS
Lightweight Lead	ADGR-LIS	ADGR-LIM	ADGR-LIL	ADGR-LIX	ADGR-LRS
No Lead	ADGR-NOS	ADGR-NOM	ADGR-NOL	ADGR-NOX	ADGR-NRS

Protect-A-Tech™ Demi-Guard Half Apron with Buckle					
Weight	Small (8" x 10")	Medium (12" x 12")	Large (16" x 14")	X-Large (20" x 16")	Set with Rack*
Regular Lead	ADGB-RES	ADGB-REM	ADGB-REL	ADGB-REX	ADGB-RRS
Lightweight Lead	ADGB-LIS	ADGB-LIM	ADGB-LIL	ADGB-LIX	ADGB-LRS
No Lead	ADGB-NOS	ADGB-NOM	ADGB-NOL	ADGB-NOX	ADGB-NRS

All custom aprons are non-returnable.
Specify color (p. 13) when placing your order.

*Set includes 1 each: S, M, L, XL plus steel wall rack.

