

QUALITY CONTROL

SPEE-D-WIPES™/SPEE-D-CLEAN™ PROCESSOR SHEETS

PDC Healthcare Spee-D-Wipes™ — Absorbent, lint-free and cost effective!

- Made from a soft specially blended anti-static cellulose/polyester material
- Rid screens of dirt and dust particles which can cause artifacts

Product Number	Description	Qty.
MAM-340	Lint-Free Wipes (9" x 9")	100/Box

The ACR recommends:

Use a "lint-free gauze pad, equivalent lint-free cloth, or camel's hair brush to assure that cassettes and screens are free of dust and dirt particles that may degrade image quality."

ACR Committee on Quality Assurance



Clean processor transport rollers without removing them or using any hazardous chemicals.

- Unique moisture activated cleaning formulation strips away the "scuzzies" from any size transport rollers in a single pass without removing the rack
- Eliminates the need for hazardous chemicals use


Spee-D-Clean™ Processor Sheets

- Same quality as Easy Clean™
- 9% Lower cost

Product Number	Size	Qty.
SC1417	14" x 17"	50 Sheets/Box

Easy Clean™ Processor Sheets

Product Number	Size	Qty.
EC-1417	14" x 17"	50 Sheets/Box



Run a sheet through your film processor for "scrub brush clean" rollers and artifact free images. Economical to use – less wasted film and less downtime to reduce your operating costs.

See our new website with new features!

- Access all PDC, St. John, Pharmex and TimeMed products on one website!
- My Account 360° empowers you with instant account information
- Easier online ordering saves time
- Streamlined navigation & advanced search help you find information fast
- User-friendly tools help you explore products
- Instant online support at your fingertips

