



The complete source for all your skin marking needs.


Switch to PDC Healthcare and Save over Beekley Skin Markers!

Why buy skin markers from anyone else and pay more?


See how one customer saves a total of \$30,536 – visit us online at www.pdchealthcare.com/sdm-costanalysis for a FREE cost-savings analysis.

- High quality — 100% Unconditional Ironclad Guarantee
- Analog/Digital — markers available for both technologies
- NO BURNOUT™ — perfect for MLO and dense breasts
- Reduce retakes — reduce patient anxiety
- No minimum order required — buy one box and save
- No bulky packaging — easy-to-use recyclable boxes
- 100% Lead-free & Latex-free — environmentally friendly and patient safe

Radiopaque



Radiolucent — NO BURNOUT™



U.S. Patent No. 7,602,883

NEW! Swirl Collection™ Skin Markers



U.S. Patent No. 7,602,883



A portion of our proceeds will be donated to the American Cancer Society.

Customers prefer PDC Healthcare Spee-D-Mark™ skin markers!

"The Sheila R. Veloz Breast Imaging Center has been purchasing Spee-D-Mark™ skin markers from PDC Healthcare for many years. Our facility performs primarily digital mammograms; we have been pleased to learn that all of their skin markers are designed to work with both digital and analog technologies.

PDC Healthcare has developed new products to assist us in our mission to better serve our patients. Their new radiolucent skin markers were designed to provide clear identification and prevent burnout. We will continue to purchase the Spee-D-Mark™ skin markers from PDC Healthcare for all of our skin marking needs."

Terry Bucknall
Director of Women's Imaging Services
Sheila R. Veloz Breast Imaging Center


SKIN MARKERS

MAMMOGRAPHY

NIPPLE MARKERS

Spee-D-Mark™ Skin Markers — Radiopaque

Virtually eliminates repeat examinations in mammography with clear nipple identification. Available in **Swirl Collection™**, **WILD ANIMAL PRINT™** and **Super-sticky**.



SDM-BB15 - 1.5mm BB



SDM-BB20 - 2.0mm BB



SDM-BB25 - 2.5mm BB

Qty.
100/Box




SDM-PBB15 - 1.5mm BB
SDM-PBB20 - 2.0mm BB
SDM-PBB25 - 2.5mm BB



SDM-TBB15 - 1.5mm BB
SDM-TBB20 - 2.0mm BB
SDM-TBB25 - 2.5mm BB



SDM-ZBB15 - 1.5mm BB
SDM-ZBB20 - 2.0mm BB
SDM-ZBB25 - 2.5mm BB



SDM-GBB15 - 1.5mm BB
SDM-GBB20 - 2.0mm BB
SDM-GBB25 - 2.5mm BB

Our Super-sticky Spee-D-Mark™ markers were designed to adhere to patients with more difficult skin types.



SDM-HBB15 - 1.5mm BB



SDM-HBB20 - 2.0mm BB



SDM-HBB25 - 2.5mm BB

Spee-D-Wing™ Skin Markers — Radiopaque


Provides added support for larger, more prominent nipples. Available in **WILD ANIMAL PRINT™** and **Super-sticky**.



SDM-W15 - 1.5mm BB




SDM-W20 - 2.0mm BB



SDM-W25 - 2.5mm BB


Qty.
100/Box



SDM-TW15 - 1.5mm BB
SDM-TW20 - 2.0mm BB
SDM-TW25 - 2.5mm BB



SDM-ZW15 - 1.5mm BB
SDM-ZW20 - 2.0mm BB
SDM-ZW25 - 2.5mm BB



SDM-GW15 - 1.5mm BB
SDM-GW20 - 2.0mm BB
SDM-GW25 - 2.5mm BB

Our Super-sticky Spee-D-Mark™ markers were designed to adhere to patients with more difficult skin types.



SDM-HW15 - 1.5mm BB



SDM-HW20 - 2.0mm BB



SDM-HW25 - 2.5mm BB

Spee-D-Strip™ Skin Markers — Radiopaque

Our adhesive-free center provides patient comfort.



SDM-NS15 - 1.5mm BB



SDM-NS20 - 2.0mm BB



SDM-NS25 - 2.5mm BB

Qty.
100/Box

NIPPLE MARKERS

Spee-D-Stretch™ Skin Markers — Radiopaque

Our flexible Spee-D-Stretch™ with super-stretchy material ensures precise and easy placement for larger, more prominent nipples.




SDM-SS15 - 1.5mm BB
Super-stretchy



SDM-SS25 - 2.5mm BB
Super-stretchy

Gauze backing for comfort.
Gauze is only located on the center



SDM-SS20 - 2.0mm BB
Super-stretchy


Super Stretchy
For maximum flexibility & patient comfort!

Qty.
100/Box

Spee-D-Stretch™ Repositionable Skin Markers — Radiopaque


The only repositionable skin marker for multiple views.

- Super-sticky for all skin types
- Repositionable for multiple views



SDM-RSS20 - 2.0mm BB
Super-stretchy • Repositionable

Gauze backing for comfort.
Gauze is only located on the center



SDM-RSS25 - 2.5mm BB
Super-stretchy • Repositionable

Super Stretchy
For maximum flexibility & patient comfort!

Qty.
100/Box

Digital Mammography Skin Markers — Radiopaque

Our radiopaque digital skin markers allow for clear nipple identification without distracting artifacts.

Spee-D-Mark™




SDM-BB23 - 2.3mm BB

Spee-D-Strip™



SDM-NS23 - 2.3mm BB

Spee-D-Stretch™



SDM-SS23 - 2.3mm BB
Super-stretchy

Gauze backing for comfort.
Gauze is only located on the center

Qty.
100/Box

SKIN MARKERS

FREE* "Easy-mount" Dispenser

No bulky packaging, unlike other skin marker boxes!

- Environmentally friendly recyclable packaging
- One-size-fits-all
- No assembly required
- Holds up to 5 Spee-D-Mark™ boxes
- Dimensions: 6-1/4"W x 3-3/8"H x 3-3/8"D

*With your first order of any skin marker product. Limit one per customer. Must mention code PRDSPSDM1 for free dispenser.



Space Saving Design!



Spee-D-Mark™ skin marker boxes not included.

SKIN MARKERS

MAMMOGRAPHY


MOLE MARKERS

Spee-D-Ring™ Skin Markers — Radiolucent NO BURNOUT™


Our radiolucent skin markers will not burnout. These skin markers are ideal for MLO exams as well as dense breast tissue. 1/2" inside diameter and dual density. Available in **Swirl Collection™**, **WILD ANIMAL PRINT™** and **Super-sticky**.



SDG-PM612 - 1/2" Diameter



SDG-M612 - 1/2" Diameter



Spee-D-Ring™ Dual Density Views




SDG-TM612 - 1/2" Diameter



SDG-ZM612 - 1/2" Diameter
©2009 The St. John Companies, Inc.




SDG-GM612 - 1/2" Diameter




SDG-SS612 - 1/2" Diameter
Super-stretchy

Super Stretchy
For maximum flexibility
& patient comfort!



SDG-RSS612 - 1/2" Diameter
Super-stretchy • Repositionable

Our Super-sticky Spee-D-Mark™ markers were designed to adhere to patients with more difficult skin types.



SDG-HM612 - 1/2" Diameter

Qty.
100/Box

All Mole Markers: U.S. Patent No. 7,602,883


SCAR MARKERS

Spee-D-Line™ Skin Markers — Radiolucent NO BURNOUT™

A valuable marking tool to reduce confusion from complex scar tissue. Our unique flexible marker ensures precise and easy placement. Available in **Swirl Collection™** and **Super-sticky**.



SDW-PSS612 - 1.0mm line
Super-stretchy



SDW-SS612 - 1.0mm line
Super-stretchy

Our Super-sticky Spee-D-Mark™ markers were designed to adhere to patients with more difficult skin types.



SDW-HSS612 - 1.0mm line



Spee-D-Line™
Dual Density View

Qty.
120*/Box

U.S. Patent No. 7,602,883

Spee-D-Line™ Skin Markers — Radiopaque

A valuable marking tool to reduce confusion from complex scar tissue. Our unique flexible marker ensures precise and easy placement.





SDW-NS1 - 0.2mm line

Qty.
120*/Box

PALPABLE MASS & PAIN / AREA OF CONCERN MARKERS

Spee-D-Angle™ Skin Markers – Radiolucent NO BURNOUT™

For palpable nodules, signaling the radiologist of the irregularity without obscuring area of concern on film.



SDM-T612
U.S. Patent No. 7,602,883


Qty.
100/Box



A portion of our proceeds will be donated to the American Cancer Society.

Spee-D-Arrow™ Skin Markers – Radiolucent NO BURNOUT™

Indicates point of entry, foreign bodies and pain without obscuring area of concern.



SDM-A612
U.S. Patent No. 7,602,883

Qty.
100/Box

Call (800) 435-4242 for samples

See how much you can save with PDC Healthcare Spee-D-Mark™ skin markers. It's simple!

Log onto www.pdchealthcare.com/sdm-costanalysis and provide us with your current skin marker usage so we can put together an annual cost-savings analysis, **GUARANTEED** to save you money. Or call us at (800) 435-4242.

ACTUAL CUSTOMER SAVINGS — \$30,536 / Year!

Item	BEFORE switching to PDC Healthcare, this customer was paying too much!						AFTER switching to PDC Healthcare, this customer SAVES \$30,536 / Year!							Annual Savings
	Previous Vendor's Part No.	No. of Boxes / Year	Markers / Box	Markers / Year	Price / Box	Marker Price	PDC Healthcare Equivalent Part No.	No. of Boxes / Year	Markers / Box	Markers / Year	Price / Box	Marker Price	Savings / Marker	
1	O-Spots® 652	240	58	13,920	\$77.00	\$1.33	SDG-M612	139	100	13,920	\$30.56	\$.31	\$1.02	\$14,226.05
2	V-Spots® 603	240	130	31,200	\$77.00	\$.59	SDM-BB25	312	100	31,200	\$30.56	\$.31	\$.29	\$8,945.28
3	A-Spots® 151	96	120	11,520	\$77.00	\$.64	SDM-T612	115	100	11,520	\$30.56	\$.31	\$.34	\$3,871.49
4	S-Spots® 650	96	129"	12,359"	\$77.00	\$.60	SDW-SS612	103	120"	12,359"	\$37.86	\$.32	\$.28	\$3,492.72
Total Savings													\$30,535.54	

Cost-savings analysis dated March 2008.

For more Mammography Products, see pages 84–91

Spee-D-Cool™ Gel Pads



Pages 84–85

Mammo Wipes



Page 86

Spee-D-Grids™



Page 87

RADIOLOGY & ONCOLOGY MARKERS

Spee-D-Mark™ CT Scan Skin Markers


- Easily Identifiable – appears bright white on the image to properly communicate an area of interest or concern to the radiologist
- Artifact-free – helps generate artifact-free images, resulting in reduced retakes
- Reduce Retakes – reduce patient radiation dose and patient anxiety
- 100% Latex-free and Hypoallergenic – safe for patients

Perfect for small areas like head and neck scans.




SDM-CT1 - 2.3mm BB

Qty.
50/Box



Perfect for large areas like chest and abdomen scans.



SDM-CT2 - 4.0mm BB

Qty.
50/Box



Spee-D-Mark™ MRI Skin Markers

- Easily Identifiable – clearly and accurately communicates an area of interest or concern to the radiologist
- Artifact-free – helps generate artifact-free images, resulting in reduced retakes
- Does not indent the skin, preventing a compromised image
- 100% Latex-free and Hypoallergenic – safe for patients

Perfect for small areas like extremities and breasts.




SDM-MR2
3/4" Capsule

Qty.
40/Box



Perfect for large areas like abdomen and spine.



SDM-MR3
1-1/4" Capsule

Qty.
40/Box



Indicator™ Radiopaque Markers for Radiology and Radiation Therapy Used in CT, Diagnostic Radiology, Angiography, Mammography and Fluoroscopy

Indicator™ Crosses

- Perfect for identifying masses, moles, birthmarks or points of interest
- Mark portal fields or set up reference points in Radiation Therapy



Product Number	Description	Qty.
SDM-RT02	10mm	130/Box
SDM-RT03	20mm	

Indicator™ Lines

- Perfect for marking scars
- Perforated every 2" for easy use



Product Number	Description	Qty.
SDW-CT05	0.5mm	138"/Box
SDW-CT10	1mm	132"/Box
SDW-CT15	1.5mm	125"/Box

Spee-D-Mark™ Skin Markers — Radiopaque

Ideal for radiation site identification. The large BB was designed to withstand radiation therapy.



SDM-BB30 - 3.0mm BB



SDM-BB40 - 4.0mm BB



SDM-BB50 - 5.0mm BB


Qty.
50/Box

Spee-D-Arrow™ Skin Markers — Radiopaque

The perfect tool when a directional marker is required. Features an arrow to indicate point of entry, foreign bodies and pain.



Our "easy-lift" tabs provide an adhesive-free access point for comfortable removal of the marker after the procedure.



SDM-A1

Qty.
100/Box

Spee-D-Left™ & Spee-D-Right™ Skin Markers — Radiopaque

Great for portable procedures, trauma cases and isolation cases.




SDM-L1



SDM-R1

Qty.
100/Box

For imaging labels, see page 35-39.



Multi-modality markers provide a clear and accurate reference on any type of scan.

Multi-modality markers are important for Image Guided Surgery because they provide reference points on the patient's scanned images without tissue indent or artifacts.



Nuclear Medicine /PET Marker

MM3003

- Visible on Nuclear Medicine, PET, CT and MRI
- 15mm outer diameter
- 3.5mm thick
- 5mm central axis hole with sealant plastic cover that forms a well that is suitable for the injection of a short-life radioisotope using a conventional hypodermic needle

Qty.
50/Jar



CT/MRI Marker

MM3007

- Appears as a bright object on CT and MRI scans and can be seen on all MRI sequences
- Perfect for Mammography and Diagnostic Radiology procedures
- 15mm outer diameter
- 3.5mm thick

Qty.
50/Jar



Radiation Therapy /Radiology Marker

MM3005

- Appears as a bright object on CT and MRI scans and can be seen on all MRI sequences
- Perfect for Mammography and Diagnostic Radiology procedures
- 15mm outer diameter
- 3.5mm thick
- 5mm central axis hole
- 2mm inner center hole suitable for passage of a conventional hypodermic needle

Qty.
50/Jar