

pdc

ScanBand® Soft

Direct Thermal
Patient ID Wristband

The Softest, Longest-Lasting
Wristband Available



ScanBand® Soft
Softest & Most Durable
ADMIT: 01/01/26
DOB: 06/30/80
MR# 123456



ScanBand® Soft

Improve Patient Safety & Satisfaction



- Print **lasts 30% longer** than competitive bands to prevent ID errors
- Softest patient ID band in the market **prevents skin cuts** and is comfortable for even the most sensitive patients, from newborn to geriatric
- Increased wristband durability to last a patient's full stay, **up to 16 days**



INCREDIBLE SOFTNESS



COMPETITIVE PRICING



CRISP PRINTABILITY



SECURE, NON-TRANSFERABLE ADHESIVE & SNAP CLOSURES



WRINKLE, WATER, & CHEMICAL RESISTANT



LAB TESTED FOR PERFORMANCE



ANTIMICROBIAL COATING



COMPATIBLE WITH ALL LEADING THERMAL PRINTERS*

Adult Thermal Wristbands



PART #	CLOSURE TYPE	SIZE	CORE	OUTER DIAMETER	BLACK BAR LOCATION	WRISTBANDS/BOX
7511C**-PDO	Adhesive	Adult 1" x 11 1/2"	1"	6.0"	1.25" from Tail	250/roll, 6 rolls
7541C**-PDO	Adhesive	Adult 1 1/8" x 11 1/2"	1"	6.0"	1.25" from Tail	250/roll, 6 rolls
7546**-PDK	Adhesive	Adult 1 1/8" x 11 1/2"	1"	4.5"	Middle	150/roll, 2 rolls
7524**-PDL	Adhesive	Adult Narrow 3/4" x 11 1/2"	1 1/2"	5.5"	Tail	200/roll, 2 rolls
7910**-PDL	Clasp	Adult 1 1/8" x 10 7/8"	1"	5.0"	Middle	200/roll, 2 rolls

COLORS			
WHITE 11	PINK 12	BLUE 13	GREEN 15

More colors available upon request

NEW Clasp Closure



Infant & Pediatric Thermal Wristbands



PART #	CLOSURE TYPE	SIZE	CORE	OUTER DIAMETER	BLACK BAR LOCATION	WRISTBANDS/BOX
7571**-PDM	Adhesive	Pediatric 1" x 7"	1"	4.5"	Tail	250/roll, 2 rolls
7527**-PDL	Adhesive	Infant Narrow 3/4" x 7"	1"	4.1"	Tail	200/roll, 2 rolls
7951**-PDL	Clasp	Infant 3/4" x 7"	1"	4.5"	Middle	200/roll, 2 rolls

COLORS		
WHITE 11	PINK 12	LIGHT BLUE 18

More colors available upon request

*Go to pdchealthcare.com/ScanBandSoft for a full list of compatible thermal printers.

Call **800.435.4242** to learn more, request samples or place an order.