

IT'S TIME FOR A CHANGE

Introducing the next generation IV change reminder that provides a visual cue and accountability at the push of a button.

TimeAlert. IV



Change is critical

- To prevent the risk of blood stream infections, the CDC established guidelines for IV tubing changes
- Human errors may occur that affect change timelines, such as illegibility, wrong date/time, or missing label

INFECTION RISK

96%

of hospital infection control professionals agree that changing IV tubing on time helps reduce the risk of infection.*

COST OF CLABSI

\$290,000

Average Annual Hospital CLABSI Costs

Estimate based upon 96-hour protocol, US hospital average infection rate 1.14/1000 central line days, average line days of 5,649, and average cost of \$45,000 per BSI.**

COMPLIANCE

86%

of nursing leaders agree that visual reminders are effective in prompting nurses to take action.***

*TechValidate. Survey of 140 infection control professionals. <https://www.techvalidate.com/tvid/4F3-D60-EDC>

**Pronovost P, Needham D, Berenholtz S, et al. An intervention to decrease catheter-related bloodstream infections in the ICU. N Engl J Med 2006;355(26):2725.

***TechValidate. Survey of 112 hospital nurse managers, directors, and CNOs. <https://www.techvalidate.com/tvid/FOF-E16-507>

TimeAlert® IV Provides a Visual Cue & Reminder When It's Time to Change


TIME-INDICATING REMINDER LABELS FOR IV TUBING CHANGES

With TimeAlert® IV, any caregiver will know at-a-glance how much time has elapsed.

- Smart indicator delivers an easy-to-read visual cue
- Reduces human error potential
- Seamlessly integrates into existing protocols
- Visual reminders assist nurse compliance
- Push-button activation plus color-changing technology clearly indicates timing
- Single reminder label can replace multiple labels
- Strong adhesive securely attaches to IV tubing

When it's this critical, it should be this simple.

TimeAlert® IV Time-Indicating Reminder Label



PRODUCT #	DURATION	LABEL SIZE	QUANTITY
TIMEALERTIV-7296	72/96 hours	4" W x 1 1/4" H	5 labels/sheet 20 sheets/box



PRODUCT #	DURATION	LABEL SIZE	QUANTITY
TIMEALERTIV-7D	7 days	4" W x 1 1/4" H	5 labels/sheet 20 sheets/box

Label is intended for use only as a secondary reminder. Healthcare professionals are to follow their organization's protocol for changing IV tubing and other use cases that utilize this technology.

PARTNER WITH THE IDENTIFICATION EXPERTS

With more than six decades of experience, PDC is a trusted leader in identification and patient safety solutions. Over 90% of hospitals rely upon PDC products today, which help reduce medical errors, minimize risk, and improve workflow efficiency.



Get free samples and information about a trial
Call 800 435 4242 or visit
pdhealthcare.com/timealert

PDC 25124 Springfield Court, Suite 200, Valencia, CA 91355
Tel 800 435 4242 | Fax 800 321 4409 | info@pdhealthcare.com | pdhealthcare.com

