

DuraProtect® LASER PATIENT ID WRISTBANDS

Superior Protection of
Patient Information




Innovative Materials Help Protect Patient Identification

Wristbands are used on every patient in your hospital.

Ensure patient ID is not at risk with paper-free wristbands.

The DuraProtect® line of laser wristbands surpasses wristbands of yesterday with technology and innovation. Engineered with advanced synthetic materials and intelligent design, DuraProtect® delivers optimized performance, reliability, and ease-of-use.



Synthetic material combined with a protective cover withstands:



Fading



Moisture



Rubbing




Alcohol



Hand Sanitizer

Superior Protection, Simplified Process

Unique print-and-peel design simplifies assembly.



High Performance, Low Cost

FEATURE	BENEFIT
Paper-free construction	Highly durable
Laminate shield	Protects for reliable scanning
Exclusive peel-off design	Easy to assemble, easy to use
Soft synthetic material	Patient comfort and durability
Drop-in replacement formats	Simple implementation
Lower cost	Supports hospital cost savings goals




Product Number	Sheet Size	Label Size	Sheets Per Box	Labels Per Sheet	Wristbands Per Sheet
9630-11-PDR	8 1/2" x 11"	2.5" L x 0.94" H	1,000	20	1
9631-11-PDR	8 1/2" x 11"	2.5" L x 0.97" H	1,000	20	1
9632-11-PDR	8 1/2" x 11"	2.5" L x 1" H	1,000	20	1
9633-11-PDM	8 1/2" x 11"	N/A	500	N/A	2 Adult, 2 Infant, 2 Infant Tags, 1 Crib Tag
9633F-11-PDM	8 1/2" x 11"	N/A	500	N/A	2 Adult, 2 Infant, 2 Infant Tags


Drop-in products for all LaserBand® 2 / LB2 part numbers available.

LaserBand® is a registered trademark of Zebra Technologies Corp., which is in no way associated or affiliated with PDC, Precision Dynamics Corporation or Brady Corporation. Reference to the LaserBand® trademark, name, or product line does not imply any co-operation or endorsement.


Adult 9630-11-PDR




Adult 9631-11-PDR




Adult 9632-11-PDR



Mother/Father/Baby 9633-11-PDM



Mother/Father/Baby 9633F-11-PDM



Get free samples and a quote today.
Visit pdchealthcare.com/duraprotect



PDC 25124 Springfield Ct., Ste 200, Valencia, CA 91355
Tel 800 435 4242 | Fax 800 321 4409
info@pdchealthcare.com | pdchealthcare.com

