




# Goodbye shield. Hello simplicity.

**It's about time.**

Faster application saves  
both time and money.



**DuraSoft**<sup>®</sup>  
LASER WRISTBANDS

# 50% Faster Application.

DuraSoft® laser wristbands may save the average health system more than 275 nurse hours per year.\*

Protect patient data and barcodes without the cumbersome laminate overlay with DuraSoft® laser wristbands.

- Quick and easy application
- Excellent barcode and print durability
- Patented design with advanced soft wristband material
- Economical laser wristband solution
- Latex-free, Phthalate-free, and moisture-resistant
- Drop-in replacements for all laser bands (40+ styles available)



\*Based on 100,000 wristbands applied and 10 seconds of time saved per wristband.

# EASY AS 1-2-3



## 1 PRINT

Easy drop-in replacement works with your existing format, so there's little to no change needed.



## 2 REMOVE

No laminate shield means a quick and easy, one-step application.



## 3 APPLY

Fast and easy application saves nursing and admission staff time.

DuraSoft® advanced synthetic material withstands:



Fading



Moisture



Alcohol



Hand Sanitizer

# Scans the first time, every time.




# DuraSoft® ADULT

PRODUCT NUMBER	SHEET SIZE	WRISTBAND SIZE	LABELS SIZE	SHEETS/ BOX	LABELS/ SHEET	WRISTBANDS/ SHEET
9500W-11-PDR	8.5" x 11"	9.7" L x 1" H	2.5" L x 1" H	1,000	20	1
9502W-11-PDR	8.5" x 11"	9.9" L x 1.1" H	2.5" L x 1" H		20	
9505W-11-PDR	8.5" x 11"	9.7" L x 1" H	2.5" L x .9" H		20	
9506W-11-PDR	8.5" x 11"	9.7" L x 1" H	2.5" L x 1" H		20	
9507W-11-PDR	3" x 11"	9.5" L x 1" H	N/A		N/A	

## 9502W-11-PDR\*


Drop-in replacement for Duraprint® 9325-11-PDR



\*AVAILABLE WITH ALERTS 9501-11-PDR


## 9506W-11-PDR

Drop-in replacement for LaserBand® (LB2-ADULT-L3)




## 9507W-11-PDR

Drop-in replacement for LaserBand® (PLS-101A)



## 9500W-11-PDR


Drop-in replacement for LaserBand® (PLS-102)



ORANGE  
17

## 9505W-11-PDR

Drop-in replacement for LaserBand® (LB2-ADULT-L2)




# DuraSoft® for BLOOD APPLICATIONS

Print barcode patient ID wristbands with matching specimen and document labels on a single sheet to optimize convenience and efficiency.

PRODUCT NUMBER	SHEET SIZE	WRISTBAND SIZE	LABELS SIZE	SHEETS/BOX	LABELS/SHEET	WRISTBANDS/SHEET
9519W-11-PDR	8.5" x 11"	9.7" L x 1" H	2" L x .8" H	1,000	35	1

## 9519W-11-PDR

Includes 10 poly labels for blood collection tubes and 25 paper labels for documentation



Ask about our Securline® Barcode Blood Bands. They feature serialized barcode IDs to help match the right patient to the right blood for blood transfusions.






# DuraSoft® TENDERCARE®

PRODUCT NUMBER	SHEET SIZE	WRISTBAND SIZE	LABELS SIZE	SHEETS/BOX	LABELS/SHEET	WRISTBANDS/SHEET
9516W-11-PDM	8.5" x 11"	Ped: 9.1" L x 1.1" H Infant: 7.2" L x 1" H Tag: 2.3" L x 1" H Extender: 7.3" L x .4" H	2.5" L x 1" H	1,000	16	2
9517W2-11-PDR	8.5" x 11"	Ped: 8.4" L x 1.1" H Infant: 6.3" L x 1" H Tag: 2.3" L x 1" H Extender: 7.3" L x .4" H	2.5" L x 1" H	1,000	12	1
9531W-11-PDM	8.5" x 11"	Tag: 2.3" L x 1" H Label 2.5" L x 1" H	2.5" L x 1" H	500	20	N/A
9540W-11-SPDM	8.5" x 11"	Adult: 9.7" L x 1" H Ped: 6.7" L x 1" H Tag: 2.3" L x .9" H	N/A	500	N/A	4

## Infant/Pediatric

### 9516W-11-PDR


- Drop-in replacement for LaserBand® (LB2-PED-L2)
- Infant and pediatric wristbands, plus two ID tags and an extender band on a single sheet to maximize workflow efficiency



## Infant/Pediatric

### 9517W2-11-PDR


- Drop-in replacement for LaserBand® (LB2-PED-L3)
- Infant and pediatric wristbands, plus two ID tags and an extender band on a single sheet to maximize workflow efficiency



## Infant/Pediatric

### 9531W-11-PDM

- Drop-in replacement for LaserBand® (LB2-BABY-TLC)




LaserBand® is a registered trademark of Zebra Technologies Corp., which is in no way associated or affiliated with PDC, Precision Dynamics Corporation or Brady Corporation. Reference to the LaserBand® trademark, name, or product line does not imply any co-operation or endorsement.

## Mother/Father/Baby

### 9540W-11-SPDM (serialized)

- Drop-in replacement for LaserBand® (LB2-L3-FAMILY-NRS)
- Helps ensure every newborn is properly identified and matched to their parent(s)



## Precision® NEONATAL FOAM BANDS

The softest ID bands for your most fragile patients.

These wrist and ankle bands were designed to positively identify the individual newborn, and eliminate the unsafe practice of taping patient information to the incubator. With extreme comfort in mind, these bands are made of a soft foam material that is gentle on sensitive, delicate skin.

PRODUCT NUMBER	STYLE	LENGTH OF FOAM BOARD	INFORMATION AREA	QTY/BOX
523-11-PDD	Newborn, Blue	6" L	2" L x 3/4" H	12 Bands
524-12-PDD	Newborn, Pink	6" L		
533-11-PDD	Preemie, Blue	5" L		
534-12-PDD	Preemie, Pink	5" L		



## PARTNER WITH THE IDENTIFICATION EXPERTS

With more than six decades of experience, PDC is a trusted leader in identification and patient safety solutions. Over 90% of hospitals rely upon PDC products today, which help reduce medical errors, minimize risk, and improve workflow efficiency.



Get free samples and a quote today  
Call **800.435.4242** or visit [pdchealthcare.com/durasoft](https://pdchealthcare.com/durasoft)

PDC 25124 Springfield Court, Suite 200, Valencia, CA 91355  
Tel 800 435 4242 | Fax 800 321 4409 | [info@pdchealthcare.com](mailto:info@pdchealthcare.com) | [pdchealthcare.com](https://pdchealthcare.com)

HC-222-0423

