

Some things fade fast

A middle school crush

Purple hair dye

A cheap spray tan



ScanBand® DR doesn't.

From admission to discharge, ScanBand® DR thermal wristband is designed to resist fading.



If it doesn't scan, nothing else matters.

The wristband is a patient's lifeline.

86%

Up to an 86% reduction
in medication errors when
barcode scanning is used.^{1,2,3}

A scannable wristband is critical to patient safety.

11%

11% of BCMA overrides
are due to unreadable
wristbands.⁴

Poor quality scans with inferior wristbands can put your patients at risk and lead to:



Patient
ID Errors



Decreased
Workflow
Efficiency



Staff
Workarounds
and Frustration



Reduced
Patient Care
Time and
Satisfaction



Re-Printing
and Re-Banding

References:

1. Malcolm B, Carlson RA, Tucker CL, et al. Veterans Affairs: eliminating medication errors through point-of-care devices. Paper presented at: Healthcare Information and Management Systems Society Conference; 2001 Apr; Dallas (TX). Vol. 2 (Session 73): 218-26. 2. Puckett F. Medication-management component of a

point-of-care information system. *Am J Health Syst Pharm* 1995 Jun 15;52(12):1305-9. 3. Johnson CL, Carlson RA, Tucker CL, et al. Using BCMA software to improve patient safety in Veterans Administration Medical Centers. *J Healthc Inf Manage* 2002 Winter;16(1):46-51. 4. Ross Koppel, PhD, Tosha Wetterneck, MD, MS, Joel Leon Telles, PhD,

Ben-Tzion Karsh, PhD; Workarounds to Barcode Medication Administration Systems: Their Occurrences, Causes, and Threats to Patient Safety. *J Am Med Inform Assoc* 2008; 15 (4): 408-423. doi: 10.1197/jamia.M2616

ScanBand® DR

Thermal Patient ID Wristbands

Scan with confidence—from admission to discharge.

ScanBand® DR is a thermal wristband that's specially designed to deliver optimal print performance and consistent durability.



- » Produces long-lasting, crisp, and readable barcodes and text
- » Helps reduce re-banding and workarounds
- » Made with smooth and comfortable material
- » Compatible with leading direct thermal printers
- » Antimicrobial coating

ScanBand® DR Resists:



Fading



Moisture



Rubbing



Alcohol



Hand Sanitizer



The competition fades. Choose ScanBand® DR for patient safety that lasts.

ScanBand® DR Adult



Product Number	Size of Wristband	Core Size	Wound	Outer Diameter	Wristbands/Box	Colors <i>Substitute ** for numbered color</i>
7458-**-PDL	1-1/8" x 11-1/2"	1.5"	Out	4.9"	200/Roll, 2 Rolls	<input type="checkbox"/> White -11 <input type="checkbox"/> Blue -13 <input type="checkbox"/> Yellow -14 <input type="checkbox"/> Green -15 <input type="checkbox"/> Red -16 <input type="checkbox"/> Orange -17
7461-**-PDO	1-1/8" x 11-1/2"	1.0"	Out	4.8"	200/Roll, 6 Rolls	<input type="checkbox"/> White -11 <input type="checkbox"/> Blue -13
7442-**-PDM	1-1/8" x 11-1/2"	1.5"	In	5.5"	250/Roll, 2 Rolls	<input type="checkbox"/> White -11 <input type="checkbox"/> Blue -13

ScanBand® DR Pediatric



Product Number	Size of Wristband	Core Size	Wound	Outer Diameter	Wristbands/Box	Colors <i>Substitute ** for numbered color</i>
7471-**-PDO	1" x 7"	1.0"	Out	4.6"	300/Roll, 6 Rolls	<input type="checkbox"/> White -11 <input type="checkbox"/> Blue -13

ScanBand® DR Infant



Product Number	Size of Wristband	Core Size	Wound	Outer Diameter	Wristbands/Box	Colors <i>Substitute ** for numbered color</i>
7481-**-PDO	1" x 6"	1.0"	Out	4.6"	350/Roll, 6 Rolls	<input type="checkbox"/> White -11 <input type="checkbox"/> Blue -13

Partner with the identification experts.

Over 90% of hospitals nationwide depend upon PDC identification solutions. No one understands the critical role of positive identification better than PDC. That's why we've been trusted by the leading hospitals for over 65 years.

Learn more. Get free samples.
Visit pdchealthcare.com/durable



PDC 25124 Springfield Ct., Ste 200, Valencia, CA 91355
Tel 800 435 4242 | Fax 800 321 4409
info@pdchealthcare.com | pdchealthcare.com



HC228-0623